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ENGLISH

Product Presentation

The product is packed in double blister packs, manufactured in nontoxic PVC (Polyvinyl Chloride), allowing contact with food and hospital products. The two blister packs are sealed with surgical grade paper, without dyes, free of holes, through the heat sealing process. The sealed double blister packs are also placed into a cardboard packaging covered with PVC Film. The implants are sterilized in ETO with 5 years validity.

TRACEABILITY: In the packaging of the implantable component there are 6 labels containing the necessary data for implant traceability.

All the models of Biomecanica Tibial Component Accessories contain the following information recorded with laser (when possible): Biomecanica Logo - Manufacturing Batch Number – Raw Material Symbol used in the manufacture of the product – Code - Size - CE 0434

Product Description:

Biomecanica Tibial Component Accessories are non-absorbable and non-active implants, used in revision surgeries.

Models of Biomecanica Tibial Component Accessories

BKS III Revision Tibial Stem

BKS III Revision Tibial Stem is intended to provide greater stability for Biomecanica BKS III Tibial Component used in revision surgery. It is manufactured from cobalt-chromium-molybdenum casting alloy, according to ASTM F75, and must be connected to Biomecanica BKSIII Revision Tibial Component, through a conical socket. Then, it must be secured to the Tibial Component with the aid of the assembly screw (made of Ti6Al4V, according to ASTM F136). The assembly

screw must be tightened with a 3.5 mm hexagonal wrench. The component surface that is in contact with the tibial canal is sprayed with alumina oxide. The upper surface of the component, that is machined, was designed to be connected to Biomecanica Revision Tibial Component in revision surgery. The BKSIII Revision Tibial Stems can be centered or decentered, and are sold in several sizes. This variation of dimensions enables the surgeon to choose the appropriate model, according to the patient’s bone characteristics and surgical technique used, optimizing the canal filling and placement of the tibial component. The BKS III Decentered Revision Tibial Stem has a different format in cases of irregular anatomy of the tibial canal.

List of the models of BKSIII Revision Tibial Stems

2328-10-060 BKS III Revision Tibial Stem - Ø 10 X 60 mm; 2328-10-070 BKS III Revision Tibial Stem - Ø 10 X 70 mm; 2328-10-080 BKS III Revision Tibial Stem - Ø 10 X 80 mm; 2328-10-100 BKS III Revision Tibial Stem - Ø 10 X 100 mm; 2328-12-060 BKS III Revision Tibial Stem - Ø 12 X 60 mm; 2328-12-070 BKS III Revision Tibial Stem - Ø 12 X 70 mm; 2328-12-080 BKS III Revision Tibial Stem - Ø 12 X 80 mm; 2328-12-100 BKS III Revision Tibial Stem -

retractor (24.5x252 mm); 5360-00-000 - Hohmann retractor (24.5x252 mm); 5367-00-000 - screw extractor; 5228-00-000 - Osteotoma 20mm; 5307-00-000 - Quick coupling; 5372-00-000 - Screw pointer; 5333-00-000 - Tibial resection guide; 5211-00-000 - Osteotome 22mm; 6003-02-032 - Drill bit with quick coupling Ø3,2x150 mm; 5381-00-000 - Patellar guide 3 pegs; 5383-00-000 - Test femoral extractor; 5382-00-000 - Paquimeter; 5343-01-000 - Patellar drill - Ø7,0 mm; 5384-00-000 - Patellar clamp; 5385-00-000 - Extramedullary tibial cutting guide; 5386-00-000 - Plateau inserter; 5313-12-000 - Rod for femoral alignment guide 120mm; 5313-15-000 - Rod for femoral alignment guide 150mm; 5313-20-000 - Rod for femoral alignment guide 200mm; 5324-10-000 - Intercondylar guide P/M; 5353-28-000 - Patellar Trial Ø28 mm (3pegs); 5353-30-000 - Patellar Trial Ø30 mm (3 pegs); 5353-32-000 - Patellar Trial Ø32 mm (3pegs); 5353-34-000 - Patellar Trial Ø34 mm (3 pegs); 5353-36-000 - Patellar Trial Ø36 mm (3pegs); 5211-16-000 - Osteotome 16mm; 5211-25-000 - Osteotome 25mm;

Indications, Contraindications, Adverse Effects and Warnings Indications

Biomecanica Tibial Component Accessories are intended to provide greater stability to Biomecanica Tibial Component in cases the surgeon considers the bone stock to be inappropriate or insufficient for the femoro-tibial anatomic alignment and fixation with bone cement, in Revision Total Knee Replacement.

Contraindications

The contraindications below do not limit the use Biomecanica Tibial Component Accessories; however they must be examined by the surgeon: Quadriceps Weakness; infections or systemic diseases near joints; Knee reartrodese after arthroplasty; Previous patellectomy; Neuropathic arthritis; Discontinuity or severe dysfunction of the extensor mechanism; Recurvatum deformity due to muscle weakness and presence of a knee arthrodese; Morbid obesity, fever or leukocytosis; Patient's inability to receive anesthesia; Leg atherosclerotic disease; Skin problems such as psoriasis; Neuropathic arthropathy; Recurrent urinary tract infections; Prevailing cases of osteomyelitis around the knee; Intense neuromuscular, muscular or vascular deficiencies in the affected area; Bone structures seriously damaged creating severe risk for the application of the prosthesis; Bone tumors near the area where the component will be implanted; Muscular atrophy; Bone immaturity; Knee extensor mechanism failure; Structural failure of knee soft tissues; Knee anklyosis; tibia, femur and patella insufficient bone structure; Patients who are unwilling or unable to follow post surgery cares and instructions.

Adverse Effects

Joint line migration; Walking standard change; Patellar tendon rupture; Patellar height decrease or increase; Hematomas; Wound dehiscence; Progressive and asymptomatic bone resorption (osteolysis) around the component; Loosening, migration or fracture of the medical product; Superficial or deep infection, vascular disorders, including thrombosis and pulmonary embolism; Limb shortening because of bone resorption; Allergic reaction to foreign body that may result in histological reactions involving various types of macrophages and fibroblasts; Pain, discomfort or abnormal sensations depending on the presence or migration of the medical product, and also according to the surgical procedure; nerve or neurological damage due to surgical trauma (including paralysis and soft tissue injuries); inability to perform normal daily physical activities; periprosthetic fractures.

Warnings

Biomecanica Tibial Component Accessories help stabilize Biomecanica Tibial Component, that is intended to replace the damaged joint and not normal structures of the human body. Surgical implants must never be reused, and the explanted implants must never be implanted again. Small defects of internal pressure patterns may occur, even if the implant appears to be intact, causing rupture and failure of the implant;

The correct selection of Biomecanica Tibial Component Accessories is extremely important and must be in accordance to the patient's bone characteristics examined by the surgeon. Infections often generated are linked to comorbidity; Aseptic failure of total knee arthroplasty can be caused by components loosening, wear of the ancillary component made of polyethylene, laxity, periprosthetic fracture, arthrofibrosis and patellofemoral complications; The results presented in revision surgeries are lower when compared with the results presented in primary surgeries; Patients who smoke are most likely to go through a longer operative period and have higher costs. Knee replacement has a finite expected survival that is adversely affected by activity level.

The patient should be informed of the need of periodic medical monitoring to observe possible changes in the condition of the implant and adjacent bone. Only by monitoring, possible component loosening or the occurrence of osteolysis can be detected; Not performing revision surgery when the component releases or osteolysis occurs may result in progressive loss of periprosthetic bone stock. The patient must fully understand and be warned that the product does not replace and does not have the same performance of the natural bone; The patient should be instructed to inform that s/he has a metal implant in case s/he needs to go through magnetic resonance imaging scans, because the metal implant impairs the diagnostic examination and it may cause displacement of the implant causing him/her trouble. The patient must be advised of the need to attend regular examinations during postoperative follow-up throughout the period in which the product remains in place.

Allowable Combination with Other Medical Products

According to NBR ISO 21534 - Non-active surgical implants - Joint replacement implants - Particular requirements, the cast cobalt-chromium-molybdenum alloy and 6-aluminum 4-vanadium titanium alloy, raw materials of Biomecanica Tibial Component Accessories, are acceptable for the manufacturing of implants and can be used with implants made of ultra high molecular weight polyethylene (UHMWPE), the raw material of Biomecanica Posterior-Stabilized Articular Insert.

A list of acceptable combinations with Biomecanica Tibial Component Accessories was formulated according to the standard "NBR ISO 21534 - Non-active surgical implants - Joint replacement implants - Particular requirements. Regardless of this list, the implants must be made by the same manufacturer because they were designed for such combinations observing surface finishing, surface treatment and other factors required on projects that may interfere with the combinations. Therefore, metallic implants from different manufacturers must not be used together because of chemical, physical, biological and functional incompatibility.

Allowable Contact Alloys with Biomecanica Tibial Component Accessories (when there is not articulation)

Titanium-based alloy (NBR ISO 5832-3, NBR ISO 5832-11, ASTM F136); Cobalt-based Alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBR ISO 5832-8, NBR ISO 5832-12, ASTM F75) - Stainless steel (NBR ISO 5832-9)

Special Storage Conditions, Conservation and Product handling:

- Open the packaging aseptically
- Transport and store it in a cool, dry place at a maximum temperature of 35°, and relative humidity around 30%
- Do not store it either directly on the floor (minimum height = 20cm) and or in very high places, close to the bulbs, which could cause dryness or damage on the packaging label.
- Do not store in places where contaminants, such as, cleaning materials, insecticides, pesticides, etc, are stored.

Do not use if package is damaged

Keep out of the sun

Fragile, handle with care

Keep dry

Product Care Before Use

Before using the product, the responsible person designated by the hospital should check the physical condition of the packaging and product. It is important to verify if the external cardboard packaging is wrapped in transparent PVC film and identified with the product information contained on the package. After this verification the transparent PVC film can be removed and the cardboard package can be opened. The product packed in a double blister should be carefully removed from the external packaging in a surgical environment. The blisters should be opened by pulling the medical grade paper. With product in hand, it is necessary to check if it has any grooves, scratches or stains. If any adverse effect that affects the user's security occurs, such as: problems with the performance of the product, damage to the implantable component, serious problems or deaths related to these components, the surgeon in charge

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must report this adverse event to the Health Authority and Biomecanica by e-mail sac@biomecanica.com.br or by phone 055 14 2104-7906

In cases of doubt, the surgeon or health professional can make the communication of the adverse event through the Brazilian Health Notification System at ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way Biomecanica Tibial Component Accessories explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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