# FIXADOR EXTERNO

# EXTERNAL FIXATOR

**Nome Técnico:** Fixadores Externos

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## ENGLISH

**Technical Name :** External Fixator

**Commercial Name:** Bio Tub External Fixator Line – BM – Non-sterile

## Product Description

The external fixators have two basic functions: a) Fuse the fractured parts of diaphyseal or epiphyseal bone tissue. That will promote dynamic compression or distraction or static stabilization between the fragments. b) Promote the elongation of the bone tissue in case of shortening of it. The external fixator is a medical product which is not implantable and it is designed to help the bone fixation and intended to: minimizing, stabilizing and fixing several types of fractures in the lower and upper limbs as well as in the human skeleton according to the indication for use. It also has the purpose of transferring the load applied to the fractured bone to the device which keeps the bone free of load and able to regenerate after the fracture. The external fixators present a structure designed exclusively for tubes and connectors which are connected by wires and pins. The external fixators are available in different designs and with proper methods of application besides general principles that direct their applications. The wires and pins are ancillary components for the external fixators, therefore they are validated in another process of registration and sold separately.

## Manufacturing Materials

The external fixators are made of materials that are well known in the medical field and that meet ASTM Standards for the steels and NBR Standards for the aluminum. Two types of raw-materials are used – ASTM F899 for 304 typed stainless steel; ASTM F899 for 420 typed stainless steel and NBR 6834 for 6065-T5 aluminum alloy. The table below brings the codes, the components of the set and the raw-material of each component of the fixator.

## Table - List of external fixators

Code	Components
50066-00-000	Bar-Pin connector set
50067-00-000	Bar-bar connector set
50068-00-000	Squeeze key
50069-27-000	Drill Ø2,7mm x 220mm
50070-10-000	Tube Ø9,5mm x 100mm
50070-11-000	Tube Ø9,5mm x 110mm
50070-12-000	Tube Ø9,5mm x 120mm
50070-13-000	Tube Ø9,5mm x 130mm
50070-14-000	Tube Ø9,5mm x 140mm
50070-15-000	Tube Ø9,5mm x 150mm
50070-16-000	Tube Ø9,5mm x 160mm
50070-17-000	Tube Ø9,5mm x 170mm
50070-18-000	Tube Ø9,5mm x 180mm
50070-19-000	Tube Ø9,5mm x 190mm
50070-20-000	Tube Ø9,5mm x 200mm
50070-21-000	Tube Ø9,5mm x 210mm
50070-22-000	Tube Ø9,5mm x 220mm
50070-23-000	Tube Ø9,5mm x 230mm
50070-24-000	Tube Ø9,5mm x 240mm
50070-25-000	Tube Ø9,5mm x 250mm
50070-26-000	Tube Ø9,5mm x 260mm
50070-27-000	Tube Ø9,5mm x 270mm
50070-28-000	Tube Ø9,5mm x 280mm
50070-29-000	Tube Ø9,5mm x 290mm
50070-30-000	Tube Ø9,5mm x 300mm
50070-31-000	Tube Ø9,5mm x 310mm
50070-32-000	Tube Ø9,5mm x 320mm
50070-33-000	Tube Ø9,5mm x 330mm
50070-34-000	Tube Ø9,5mm x 340mm
50070-40-000	Tube Ø9,5mm x 400mm
50071-00-000	Schanz pin guide
50071-01-000	Drill guide Ø3,3mm
50069-35-000	Drill Ø3,5mm x 220mm
50072-00-000	Wrench 10/13

## Symbology of the Traceability tag (Label) according to NBRISO 15223 and

	No Sterile		XXXXX Code		Validity
	XXXXX Lot n°		Manufacturing date		Manufacturing data
	REP	European representative			

## TECHNICAL INFORMATION

**Indications, Precautions, Restrictions/Warnings, Adverse effects and Contraindications:** The external fixator is indicated for bones of lower and upper limbs in corrective surgeries or in the treatment of traumas because the fixators keep the reduction of the fractured bone and stabilizes the bone or adjacent structures in order to promote or correct the bone fusion. The surgeon must be qualified and familiar with osteosynthesis and its limitations including the pre and postoperative, surgical techniques to be adopted, precautions and potential risks. Do not use the product if it is damaged. The pre-operative and surgical procedures as well as the knowledge about the surgical technique and proper selection and implantation of the implant are important considerations for a successful use of the provisional devices for external fixation. Do not use the product in case the packaging is damaged. Some potential adverse effects that may occur in the use of the external fixators are: postoperative superficial or deep infection; pain;

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Registro ANVISA: 80128580113

B046B - REV01 - 09012017

discomfort or abnormal sensations due to the presence or migration of the medical device or due to the surgical procedure and migration; limb edema; psuedarthrosis; movement limitations; failure of fracture healing. Due to the incompatibility of the materials, design and dimension the use of components from different manufacturers is forbidden. The use of different manufacturers may cause problems to the product and bring risks. The company disclaims all responsibility of any adverse event that may occur. Do not use implants from different manufacturers. Therefore it is recommended that the products have the same provenance.

**Sterilization of the external fixators:** The external fixators are distributed not sterile. Before using the external fixators they should be sterilized. It is recommended the steam sterilization in Autoclave at hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices)

**AUTOCLAVE:** It is a device for sterilization through saturated steam under pressure.

**INSTRUCTIONS:** Before starting the device check: if the switch is on, if the water valve is open and the flush valve is closed.

1 – open the door of the equipment; 2 – properly settle the material to be sterilized; 3 – close the door of the equipment; 4 – select the desired cycle according to the material to be sterilized; 5 – Turn on the switch; 6 – the cycle will run automatically in the sequence; 7 – when the lamp "cycle end" turns on partially open the door for about ten (10) minutes to cool down the material. It is recommended to follow the parameters of physical sterilization in autoclave (saturated steam) showed below.

#### Table – Parameters of physical sterilization in autoclave

Ciclo	Temperature	Time of exposure
Conventional (1 atm of pressure)	121°C (250°F)	30min.
Conventional (1 atm of pressure)	132°C (270°F)	15min.
Gravity	132°C (270°F)	45min.
High vacuum	132°C (270°F)	7min.

**Note:** Register the time when the heat of the sterilization chamber reaches the desired temperature.

**For further information check the instruction manual that comes with each autoclave.**

**Other method of sterilization:** Besides the autoclave the sterilization process may be in Ethylene Oxide (E.T.O) – parameters and procedures established by the validation protocol and ISO 11135– Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

**NOTE:** The methods of sterilization, the equipments, the controls and the instructions on the applied sterilization under responsibility of the hospital. Storage conditions of the sterilized articles: - as the environment: it must be cleaned; ventilated and dried and it should be restricted to the team of the sector, - as the articles: after the sterilization process do not settle them on cold surfaces (Stone or stainless steel), use baskets or perforated recipients until they cool down; - wrapping (raw-cotton fabric, nonwoven fabric, surgical paper, crepe paper, film paper, tyvec or perforated metallic boxes) must remain intact and little handled in order to avoid the packages to tear or lose the seal; - store it in cabinets with shelves; - identify the shelves in order to make it easier for the removal of the material; the material should be stored according to the expiration date of the sterilization in order to make it easier for the distribution of the material and to prevent having expired material in the stock; store the sterile materials separately from the non-sterile ones in order to reduce the level of external contaminants.

**Storage and transport of the product:** The storage location should be clean, dry and illuminated in a way that the proper conditions of storage and transport are preserved as well as the physical and chemical integrity. The product must be stored and transported in dry and fresh environment at room temperature (Max. 35°C) and relative humidity around 30% to 70% and away from direct sunlight. The effects of vibration, shocks, corrosion, defective settlement during the handling and transport, improper piling in the storage must be avoided. The product must be carefully handled in order to avoid damages that may impair the quality of the material and the safety of the patient. The products can not be stored in direct contact with the floor (minimum height = 20cm). They can not remain in very high places close to lamps because it may respect the packaging or damage the label. They should not be stored with contaminant substances such as cleaning materials, insecticide, pesticide, etc. The transporter must be advised of the content and delivery deadline. Carefully transport the products in order to avoid fall and attrition in order not to damage the device.

**Identification:** The external fixators come with labels that present the following information: Data of the manufacturer; Authorized representative in the European community; Technician in charge; nº of the register with ANVISA; Technical name of the product; Commercial name of the product; Commercial model of the product; Raw material; Quantity of items per package; Lot number; Code of the product; Manufacturing date; Expiration date; Address of the manufacturer; The external fixators come with the following engraved information for traceability: Biomecânica logo; Number of the manufacturing lot and code.

**Customer Service:** In case it is necessary to make a complaint about the use of the external fixators related to some adverse effect that may affect the safety of the user the responsible surgeon must inform the Health Authorities and Biomecânica through the email - e-mail [saca@biomecânica.com.br](mailto:saca@biomecânica.com.br) – or by the phone - 0xx14 2104 7900. In case of doubts the responsible surgeon or the health Professional may inform about the adverse event through NOTIVISA on the webpage of ANVISA: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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MDSS GmbH