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# CEMENTATION KIT

Português

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English

## Product Description

The Cimentation Kit accessories are used in surgery to improve the quality of the cement and to facilitate the cementation technique and improve the bone fixation.

## Construction Materials

Table 1 summarizes the code, set components and raw materials of each component of the Cimentation Kit.

Cimentation Kit Model

Code	Description	Set Components	Raw material	Illustration - Applicators
1931-00-000	Vacuum Cement Kit	01 Mixer Tube / Vacuum Cementation Injector	Polypropylene	
		01 Injection Tube Piston	Polypropylene	
		01 Funnel	Polypropylene	
		01 Filtro	Polypropylene	

## Symbology of the Traceability Label (plate) as defined by NBRISO15223 and EN980

LOT XXXXX	Lot No.	XXXX-XX	Validity	EC REP	European Representative	REF XXXX-XX	Code
	Manufacturer Data		XXXX-XX	Date of Manufacture	STERILE R	Gamma Radiation	

**Indications, Precautions, Restrictions / Warnings, Contraindications and Adverse Effects:** The Cimentation Kit is used to assist the surgeon during the bone cementation in for hip arthroplasty surgery. Devices that assist the surgeon in bone cementation, where manual filling would not be possible. They should be handled in appropriate environments and with extra care. Only trained professionals should handle these devices. Surgical techniques vary according to the decision made by the surgeon, leaving him the final choice of method as well as the criteria for evaluating the results of surgery. They should be applied and adapted in accordance with the requirements and appropriate surgical techniques. Do not use the product if it is damaged. If the product is damaged, in poor condition or with a damaged package it should not be used and discarded immediately. The surgeon must be familiar with and have sufficient knowledge of the use of component, surgical technique used, precautions and potential risks. Preoperative and surgical procedures, including knowledge of surgical technique are important

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# CEMENTATION KIT

# KIT CEMENTACIÓN

considerations for a successful use of the Cementation Kit.

**Sterilization:** This product is supplied in sterilized form. Sterilization is done by Gamma Radiation. If the product has a sterilization date that is expired or a damaged packaging, it must be returned to the supplier or directly to Biomecânica.

**Storage and transportation of the product:** The site should be clean, dry and equipped with adequate lighting in order to maintain the ideal conditions for storage and transportation, as well as its physical and chemical integrity. The product should be stored and transported in a cool dry place with room temperature (max. 35 °C), relative humidity around 30% to 70% and away from direct sunlight. Cannot be stored directly on the floor (minimum height = 20cm). Cannot stay in very high places, near lamps, which can cause dryness or damage to the packaging or label. Should not be stored in locations in which are stored contaminants such as, for example, cleaning materials, insecticides, pesticides, etc.. Carry them with care to avoid falling and friction in order not to damage the packaging and the device. The effects of vibration, shock, corrosion, defective seating during handling and transportation, improper stacking during storage should be avoided. The product should be handled carefully to avoid damage that might affect the quality of the material and patient safety. The carrier should be informed about the content and delivery. The opening of the sealed blister packaging can be done manually, in an appropriate location (operating room), pulling the Tyvek. This opening should only be performed during surgery to prevent product contamination and by qualified personnel. During handling scratches and nicks must be avoided or anything that may damage or mark the components. Carry them with care to avoid falling and friction in order not to damage the packaging and component, thereby ensuring sterility of the material.

**Identification:** The Cementation Kit comes with labels which contain the information: Manufacturer Data; Authorized Representative of the European Community; Technical Manager; ANVISA Registration No.; Technical Product Name; Commercial Product Name; Commercial Model Product, Raw Material, Amount of items in the package; Batch Number; Code product, date of manufacture, date of expiry, address of the manufacturer.

**Customer Service:** If there is need for any complaints regarding the use of the Cementation Kit **related to any adverse effect that affects the safety of the user or surgeon should send feedback of this adverse event to the competent health care body and to Biomecânica via the email [tecnovigilancia@biomecânica.com.br](mailto:tecnovigilancia@biomecânica.com.br) or by phone 0xx14 2104 7900.** In case of doubt the surgeon or health care professional can send feedback of the event through the Notifications in Health Surveillance System on the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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