HASTES/PREGOS PARA SÍNTESE ÓSSEA

Sistema de Qualidade Certificado
Quality System Certified
ISO 9001 ISO 13485

STEM AND NAILS FOR

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ENGLISH

PRODUCT DESCRIPTION

Device in a tubular way, with several diameters and lengths of intramedullary use, in the treatment of fractures, he has to be fastened with screw aid as much proximal as distal. The responsible surgeon, in agreement with the adopted surgical technique, should establish the dimension.
Placement performed through instrumental aid in the surgical act.

FORM OF PRESENTATION

Non-sterile plastic packing.
Trade name, validity and the production lot are suitable in the external label of identification.
It is advisable to position the stems in tray boxes of stainless steel. Boxes projected for sterilization in autoclave.

INDICATIONS

The device for the osteosynthesis is use as auxiliary element for bone fixation in the treatment of diaphyseal fractures, medium third and for most of fractures transverse and oblique.

CONTRAINDICATIONS

All visible infection is an absolute contraindication. All effort should be made to govern the possibility of pré-operative asepis in a patient with one or more of the following abnormalities: (1) fevers and/or signs of local inflammation; (2) fast destruction of the articulation or bone absorption in roentgenograms; and (3) elevation of the rate of sedimentation inexplicable through other diseases, elevation of the count of WBC, or demarcation change in WBC with differential count.
Additionally, in distant focus of infection, such as genitourinary, lung, skin (chronic lesions or ulceration), and other places, it is a contraindication relative, due hematogenesis spill that can happen to the place of the implant. Like this, the infection focuses should be treated before, during and after the implantation.
The use of this implant it is a contraindication where there are musculature loss or compromising neuromuscular assumed in the affected member, what would make the procedure unjustifiable.

WARNINGS

Do not to reuse the product. Orthopedic implants should be use in a single patient. Although it can seem they are not damaged, previous tensions can create imperfections that can reduce the success of the implant.
The responsible surgeon, in agreement with the surgical technique adopted by the same, should establish the dimension.
Due to presentation form and its characteristics, this product does not present collateral effect; however, it is suggested to monitor the patient's re-establishment through the X Ray, according to International Norms of Safety.
Do not use magnetic resonance, due to interferences.

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PRECAUTIONS

As they are temporary, implant that will serve in the aid for the treatment of the bone fractures should accompanied of cares: in transport, manipulation, and sterilization. Do not use the product if it is damaged.

The surgeon should be aware of the possibility of the patient's physical development, could have revision need of the it implant.

ADVERSE EFFECTS

Mechanical loosening can be the result of defective fixation or infection hides. Sensibility reactions to metal in patients were rarely informed.

Implantation of strange material in woven results in histological reactions that involve several macrophages sizes and fibroblasts. The chemical importance of this effect is uncertain, as well as similar changes can happen as a precursor, or during the healing process.

STERILIZATION

Product sterilized by vapor in Autoclave.

RECOMMENDATIONS

Qualified people should only manipulate this product.

Components should not be used if your implant is from different manufacturers. Therefore, we recommended that the products have the same origin.

All implants are for single use. Therefore, when retired of the patient, it should be clean and discarded.

WBC: White cells

We recommend filling in the certificate of product quality, and even the return of Biomechanica.

PRODUCT DISCARD

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way stem and nails explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

CONTAMINATION RISKS

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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