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## PORTUGUÊS

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The mixture of different biomedical materials not compatible and/or similar may raise the risk of corrosion. It must be remembered that in spite of having biocompatible material, the product, when deployed, is subject to constant environmental changes promoting the concentration of salts, acids and alkalis and they can cause corrosion. The clinical results and durability of the implant depend on the existence of a precise surgical technique.
Manufacture date, expiration date and batch number of the product: SEE LABEL.

#### CONTRAINDICATIONS

The use of this implant is contraindicated in the following situations:

Sensitivity to metal; fever or indication of local inflammation, osteoporosis, obesity (excessive spinal load on the system produced by obese patients may compromise the setting of the spinal implant), especially patients over 102 kilograms; inadequate tissue coverage in the region to be operated; pregnancy; bone affected and injured by disease and/or infections, alcoholism or drug addiction; unexplained raising of sedimentation rate by other illnesses, increase in WBC, or suspected difference in WBC count, visible infection, bone immaturity; early ambulation and other situations that may preclude the success of the surgery of spinal implant.
WBC: White Cells

#### ADVERSE EFFECTS

a) Metal sensitivity
Reactions in patients have rarely been informed. Implantation of strange material in tissues results in histological reactions involving various sizes of macrophages and fibroblasts. The importance of this chemical effect is unclear, as well as similar changes may occur as a precursor, or during the healing process;
b) Delay or no consolidation of bone that can lead to breakage of the implant;
c)– Mechanical loosening, twisting, dismemberment, breach of the components - may be the result of defective fixing or occult infection;
d) - superficial and/or profound infection;
e)– Shortening of the segment and the spinal column due to bone resorption;
f)– Pain, discomfort or abnormal sensations due to metal implant
nonspecific placement or due to migration of implantable mechanical components, and thrombosis;
g)– Pseudo-arthritis;
h)– Fracture;
j)– Paralysis;
j)– Hematomas and bleeding in blood vessels;
k)– Pulmonary embolism;
m) – Neural or neurological problems caused by surgical trauma such as: gastrointestinal disorders;
n) – Early or late failure of the implantable metal components;
o) Chronic pressure sores;
p)– Superficial or muscular sensitivity in patients with inadequate tissue coverage in the operated area;
q)– Difficulty in performing physical activities;
r)– Allergy;
s)– Illness and Death: in all surgical procedures there is an incidence of disease and death. The patient must be warned (by the surgeon) about this possible incidence before undergoing the surgery.

#### PATIENT INFORMATION

#### THE INFORMATION ABOUT THE FOLLOWING ITEMS MUST BE SUPPLIED TO PATIENTS:

Indications – Contraindications – Information Use Policies – Possible Adverse Effects – Warnings and Precautions.
The patient must be informed about the importance of post surgery monitoring. The lack of monitoring slows the detection of post-surgical problems such as loosening of components or occurrence of osteolysis. The non-realization of revision surgery when the loosening of components or osteolysis occurs may result in progressive loss of periprosthetic bone stock.

The patient must be informed about the type of implant used, its brand, code, name, batch, the surgeon's name and surgery date.

Performance foreseen in General Requirements of ANVISA Regulation, which provides the essential safety and efficacy requirements for products, described in the resolution - RDC 56, April 6th 2001.

The non-expansive Cages - BM are manufactured in Titanium Alloy, according to ASTM F136 proven by laboratory testing reports and quality certification of the product. These materials are acceptable for the manufacture of implants as defined in Annex A of the NBRISO21534 standard. The manufacture and packaging of non-expansive Cages – BM are based on specific procedures.

The non-expansive Cages – BM are distributed on a non-sterile form and packed in transparent and nontoxic PVC plastic envelope, sealed in machines assuring its integrity. The non-expansive Cages – BM are distributed with sufficient information focusing on indications and sufficient information for the product correct handling as well as the safety of the ones involved with its handling.

#### ALLOWABLE COMBINATIONS WITH OTHER MATERIALS

Metallic alloy implants of different chemical compositions are from different manufacturers are not recommended for reasons of chemical, physical, biological and functional incompatibility.

#### ANCILLARY COMPONENTS

The following related ancillary components must be purchased separately, because they are not parts of this product.
- Sine Fixation System – Pedicular - ANVISA Registration Number 80128580081 (neither an object of this register, nor part of this product)
- Implant-System for Posterior Spine Setting - BM Plate - ANVISA Registration Number 80128580083 (neither an object of this register, nor part of this product)

Bone Graft (neither an object of this register nor part of this product). The implants may be filled with any bone graft. Such choice is the surgeon's responsibility. This bone graft begins to grow through the existing holes in the walls of the cages forming a solid bone (fusion) which shall hold the vertebrae together. This process is known as intervertebral fusion.

#### ACCESSORIES

The non-expansive Cages – BM do not have any accessory with the aim of integrating the medical product.

#### IMPLANTED PRODUCT EVALUATIONS

After the implantation during the surgery, the concerning professional must execute radiological control to verify the correct positioning of the product. The professional must be in charge of leading clinical and radiological evaluations after the surgery, since it is his/her own responsibility to do so. It has to be done according to his/her own decision toward frequency in order to avoid the implantation condition and bone consolidation evolution. In case the product is found out of the correct positioning, or indicate any non-compliance, it is the surgeon responsibility to take the most appropriate corrective action.

#### USEFUL INFORMATION FOR PREVENTING RISKS DERIVING FROM IMPLANTATION

To reduce the risks of implantation contraindications, instructions for use and all information contained in the "Instructions for Use" of the product must be strictly followed.

#### PRODUCT DISCARD

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way non-expansive Cages explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implants is also under responsibility of the hospital. Biomechanca recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

#### ADDITIONAL INFORMATION ABOUT THE PROCEDURE BEFORE USING THE PRODUCT, STERILIZATION METHOD RISK OF CONTAMINATION

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the non-expansive Cages – BM metal components get in contact with tissue and body fluids. Explanted products must be treated as highly contaminants.

#### DECONTAMINATION, CLEANING AND STERILIZATION

The product comes previously cleaned in its original packaging, and ready for sterilization. It must be handled with care to avoid contamination. If the product is subjected to conditions that undermine their state of cleanliness, the implant must be cleaned before sterilization. Below we suggest some procedures for cleaning:
The cleaning can be developed through manual or mechanical methods.

#### EQUIPMENT FOR CLEANING:

-ultrasonic Washer-
-sterilizer Washer -
-disinfecter Washer -

Manual Cleaning: the procedure is performed manually, where the dirt is removed through physical action with aid of detergent, water and artifacts such as a sponge and a brush. The manual cleaning is strongly recommended because it causes less damage to surgical implants and instruments. In manual cleaning, water must be used at ambient temperature, nylon brush should be used instead of steel or abrasive one because it can damage the material protective layer. Decontamination should be done by using a bactericide and broad spectrum antifungal solution, avoiding aggressive cleaning products, brushes made of metal so that the quality of metallic implants are not damaged. In case of using cleaning machines, the implants and surgical instruments must be placed in drawers so that they do not bump each other avoiding damage to the material.

The ones in charge of cleaning the implants and surgical instruments must be careful about the type of cleaning product used because the concentrations specified by the manufacturer. The implants and surgical instruments must be thoroughly rinsed until all residues are removed and they must be dried immediately after cleaning.

#### TYPES OF CLEANERS

Enzymatic cleaners: they are basically composed by enzymes, surfactants and solubilising. A balanced combination of these elements means that product can remove the organic matter of the material in a short period of time.

Enzymatic solutions: have excellent action for cleaning, but do not have bactericidal and bacteriostatic activity.

Enzymes: are substances produced by living cells and chemical reactions that govern the process. Once produced by cells, an enzyme can be isolated and will maintain its catalytic properties, if certain conditions are maintained in its manufacture. Enzymes are classified in three major functional groups depending on the type of substrate that they will affect: protease, lipase and amylase that act on protein substrates, fats and carbohydrates, which tend to solubilize themselves and come off the items. Currently, cleaning of articles of complex configuration is recommended to ensure cleanliness.

#### DISINFECTION:

This process that destroys micro-organisms, pathogenic or not, of the items, except for bacterial spores, by physical or chemical means.

#### LEVELS OF DISINFECTION:

High level: destroys all microorganisms except the high number of spores
-2% Glutaraldeide- 20-30 minutes.

Indication: preferably in hospitals.

Middle level: eliminates vegetative bacteria, most viruses, fungi and mycobacteria

-1% sodium hypochlorite - 30 minutes.

Indication: For UBS, day care, rest homes, nursing homes.

Low level: eliminates most bacteria, fungi and some viruses, but does not eliminat mycobacteria

-Sodium hypochlorite 0,025%

Indication: nutrition.

#### STERILIZATION

This product supplied is not sterile. Before using, it should be sterilized. We recommend the sterilization steam autoclave in the hospital (ISO 1134: 1994 Sterilization of health care products - Requirements for validation and routine control - Industrial moist heat sterilization).

#### AUTOCLAVE

It is an apparatus (as for sterilizing) using superheated steam under high pressure

#### OPERATING INSTRUCTIONS

To start the equipment, check:

if the circuit breaker is connected;
if the waterworks control device is open;
if the flush control device is closed.

1 - open the equipment door;
2 - accommodate the material to be sterilized properly;
3 - close the equipment door;
4 - select the desired cycle in accordance with the material to be sterilized;
5 - connect the electric keys;
6 - the cycle will run automatically in sequence;
7 – when the light goes on it indicates "end of cycle", partially open the door for ten (10) minutes, in order to cool the material.

# CAGE NO EXPANSIVO

It is recommended that they apply the following parameters of physical sterilization in an autoclave (saturated steam):

| Cycle                                | Temperature   | Time of Exposure |
|--------------------------------------|---------------|------------------|
| <b>Conventional (1 atm pression)</b> | 121°C (250°F) | 30 minutes       |
| <b>Conventional (1 atm pression)</b> | 132ºC (270ºF) | 15 minutes       |
| <b>Gravity</b>                       | 132ºC (270ºF) | 45 minutes       |
| <b>Vacuuous</b>                      | 132ºC (270ºF) | 7 minutes        |

Note: The time should be marked when the heat of the sterilization chamber reaches the desired temperature.

FOR BEST CLARIFICATION READ THE INSTRUCTION MANUAL OF THE AUTOCLAVE.
Another method of sterilization: the following sterilization method can be used instead of the autoclave as defined by the hospital:
- Sterilization by ethylene oxide (ETO) - parameters and procedures established in the protocol for validation and EN550: 1994 Sterilization of Medical Devices - Validation and Routine Control of Ethylene Oxide Sterilization.

#### NOTE

It is the responsibility of the hospital, the method of sterilization, equipment, controls, and instructions for sterilization used.

#### CARE WITH STERILIZED PRODUCTS

Conditions for storage of sterile items related to the environment: it must be clean, airy and dry; it should be restricted to staff; related to the product: after the sterilization process do not place it on cold surface (stone or stainless steel), use baskets or containers with small openings until it cools down; Mante (unbleached cotton fabric, nonwoven, surgical grade paper, creped paper, film paper, lycoc or perforated metal boxes) should remain intact and not be handled to prevent packages from tearing or releasing the seal; they should be stored in clean cabinets with shelves; the shelves must be identified in order to facilitate the removal of the material; the material should be stored according to the expiration date of the sterilization to facilitate distribution and not to keep products that are not valid in stock; they should be stored separately from non-sterile products to reduce the level of external contaminants.

#### REMLING OF MATERIAL TO BE ANALYZED BY THE MANUFACTURER

The implants to be analyzed by the manufacturer must be cleaned in the hospital using a bactericide and broad spectrum antifungal solution. Then, they should be disinfected or sterilized in a steam auto clave or using ethylene oxide. They should be sent to Biomecânica in a packaging identified with the cleaning and sterilization method, and data of the product.

#### AFTER SALE SUPPORT (CLIENT COMPLAINT)

If there is a need to make a complaint of non expansive cages - BM related to any adverse effect that affects the safety of the user, such as if the product is not working, implantable metal component damage, death or serious problems related to these components, the responsible surgeon must communicate this adverse event for the national health agency and Biomecânica through e-mail or phone sac@biomecânica.com.br +55 14 2104 7900. In cases of doubt, the responsible surgeon or the health professional can make the communication of adverse events through the notification system in the Health Surveillance of ANVISA: http://www.anvisa.gov.br/hotste/notiwa

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