COMPONENTE PATELAR

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ENGLISH

Product Presentation The product is packed in a double blister packaging, manufactured in nontoxic PVC (Polyvinyl Chloride). The main features are: good impact resistance, excellent dimensional stability, non-toxic product that allows contact with food and medical products.

The two blister packagings are sealed with medical grade paper without dyes and free of holes. This paper is sealed in the blister through the heat sealing process.

The sealed double blister packaging is also placed into an external cardboard packaging. The implantes are sterilized in ETO with 5 years validity.

Traceability

On the packaging of the implantable component there are 6 traceability labels containing the data that permits its traceability. One of these labels must be affixed to the patient’s records and another must be given to the patient together with his/her report; the third one must be attached to the fiscal document or invoice, and the others must be given to suppliers (so they can keep it in the Historical Record of Distribution), to the surgeon, and to the Health Insurance (if applicable). In the patient’s record, the following information is essential: name of the implant used; surgery date; product code; batch number; name of patient who received the implant; and surgeon’s name.

None of the BKS III (Biomecânica Knee System III) - Patellar Component is engraved due to their anatomy. Still, the traceability of the BKS III - Patellar Component is guaranteed by the batch number available on the product traceability labels.

Product Description

The BKS III - Patellar Component is a non-absorbable and non-active implant, used in primary or revision surgeries. It has a special format that permits articulation with Biomecânica BKS III - Femoral Component with or without ligament preservation, and can be used for the right or left knee.

The BKS III - Patellar Component is used in total knee arthroplasty, therefore, it should be used in conjunction with other implants. The BKS III - Patellar Component is sold in different dimensions with 1 or 3 pegs for fixation. This variation in size and peg-types enables the surgeon to choose the appropriate model according to the patient’s bone characteristics.

List of BKS III - Patellar Components

3 PEGS:

2360-26-000 - BKS III - Patellar Component 026 mm; 2360-28-000 - BKS III - Patellar Component 028 mm; 2360-30-000 - BKS III - Patellar Component 030 mm; 2360-32-000 - BKS III - Patellar Component 032 mm; 2360-34-000 - BKS III - Patellar Component 034 mm; 2360-36-000 - BKS III - Patellar Component 036 mm.

1 PEG:

2346-26-000 - BKS III - Patellar Component 026mm; 2346-28-000 - BKS III - Patellar Component 028mm; 2346-30-000 - BKS III - Patellar Component 030mm; 2346-32-000 - BKS III - Patellar Component 032mm; 2346-34-000 - BKS III - Patellar Component 034mm; 2346-36-000 - BKS III - Patellar Component 036mm

The BKS III - Patellar component is fixed to the patella with the use of bone cement and is in articular contact with the BKS III - Femoral Component with or without Ligament Preservation, which is fixed to the distal femur.

The only difference between models with three pegs to the one peg for fixation is that the patellar component with three pegs offers greater

fixation since it presents a larger contact area between the implant, the patella and the bone cement. However, both models work under pressure of the patellar ligament.

3 pegs 
1 peg

artificial joint that has the same functionality of the natural knee. The total knee arthroplasty aims the elimination of painful symptoms, correction of deformities, allows functional range of motion and stabilizes the knee.

In the case of elderly patients with advanced osteoarthritis, the best treatment is total knee replacement. In younger patients, the total knee replacement is indicated only in cases of tricompartmental illness (medial and lateral tibiofemoral, and patellofemoral), with the need of a future revision surgery. In cases of rheumatoid arthritis, total knee replacement is recommended for all ages because of the complications presented by the disease.

The BKS III - Patellar component is indicated in the procedures for knee arthroplasty due to:

- Osteoarthritis:** Known as the “wear-and-tear” kind of arthritis, it is a chronic condition characterized by the breakdown of the joint’s cartilage, causing the bones to rub against each other, causing stiffness, pain and loss of movement in the joint.
- Rheumatoid arthritis:** is an autoimmune disease that causes chronic inflammation of the joints. It can also cause inflammation of the tissue around the joints, as well as in other organs in the body.
- Comminuted fracture:** a fracture in which the bone involved in the fracture is broken into several pieces, making it difficult to achieve a good coaptation.

Contraindications

The contraindications below do not limit the use of BKS III - Patellar Component, but must be examined by the surgeon:

- Quadriceps weakness, infections or systemic diseases near joints; previous patellectomy; neuropathic arthritis; In younger patients, a total knee replacement is indicated only in case of tricompartmental illness (medial and lateral tibiofemoral, and patellofemoral), and requires a future revision surgery, Morbid obesity, fever or leukocytosis; intense neuromuscular, muscular or vascular deficiencies in the affected region;

Severe risk of bone structures making it infeasible the application of the prosthesis; inappropriate ligament system state for the type of component chosen; bone tumors near the site where the component will be placed; muscular atrophy; bone immaturity; failure of the knee extensor mechanism; structural failure of the knee soft tissues; Ankylosis of the knee; insufficient bone structure in the tibia, femur and patella; allergy to metal medical devices; Patients who are unwilling or unable to follow the care and instructions in the postoperative period.

Adverse Effects

Vascular damage during knee replacement is extremely rare and potentially destructive, Asymptomatic and progressive local bone resorption (osteolysis) around the component may occur; Loosening, migration or fracture of the medical product; Superficial or deep infection, vascular disorders, including thrombosis and pulmonary embolism; Shortening of limb because of bone resorption; material sensitivity or allergic reaction to foreign body that may result in histological reactions involving various types of macrophages and fibroblasts; Pain, discomfort or abnormal sensations due to the presence or migration of the medical product and surgical procedure; nerve or neurological damage due to surgical trauma (including soft tissue paralysis and injuries); inability to perform normal daily physical activities; hemorrhage of blood vessels and / or bruising;

Warnings

Biomechanica patellar component should be implanted with the help of Biomechanica instruments. Components of implants from different manufacturers should not be used together; therefore it is recommended that the products have the same provenance.

The BKS III - Patellar Component aims to replace the damaged joint and not normal structures of the human body. The BKS III - Patellar Component will be sold already sterilized with ethylene oxide, and therefore it is of single-use, and its reprocessing is forbidden; Under the packaging aseptically, the models containing three pegs offer greater fixation since it presents a larger contact area between the implant, the patella and the bone cement. The patient’s weight, activity level and attention to the instructions regarding the limitations of these conditions, take effect on load and number of cycles to which the medical device is exposed. The surgeon must be thoroughly acquainted with both medical and surgical aspects of the medical device. The correct selection of BKS III - Patellar Component model is extremely important; the BKS III - Patellar Component can not support activity levels and / or loads equal to those incurred in normal and healthy bones; Replacement of patellofemoral joint, in cases of osteoarthritis, is not necessary.

The replacement of the patellofemoral joint can lead to possible complications, such as the creation of two more interfaces (bone-cement and cement-prosthesis), changes in bone vascularization by the decreasing patellar resistance and patellar fracture caused by osteonecrosis. However, the replacement of the patellofemoral joint in rheumatoid arthritis patients should be performed because of the disease evolution. Patellar luxation may occur even where there was no joint replacement due to rotational bad positioning of the tibia, loosening of the reconstruction of medial capsule to the patella, comorbidity and severe osteoporosis. Infections, often generated, are linked to comorbidity.

The patient should be informed of:

There is a need of a periodic medical monitoring to observe possible changes in the state of the implant and adjacent bone. Only the monitoring can detect possible components loosening or the occurrence of osteolysis; Failure to perform revision surgery related to osteolysis components loosening, may result in progressive loss of prosthethic bone stock. The patient should completely understand and alert that the product does not replace or have the same performance of the normal bone and, therefore, can break, deform or loosen due to excessive efforts or activities, over the permitted ones and if applying excessive load, etc. Information regarding the topics: indications, contraindications, side effects and warnings; The patient must be advised of the need to attend regular follow-up examinations after surgery during the period in which the product remains deployed.

Allowable Combination with Other Medical Products

A list of acceptable combinations with BKS III - Patellar Component in accordance with standard NBR ISO 21534 - Implants for Surgery non-active surgical implants – Joint Replacement Implants - Particular requirements was formulated. However, independently of this combination, the implants must be of the same brand because they were designed for such combinations observing the surface finishing, surface treatment and other factors required on projects that may interfere with the combinations. Therefore, metal implants from different manufacturers are not recommended because of chemical, physical, biological and functional incompatibility.

Allowable Contact Legases

Implant - BKS III - Patellar Component:

Articulation (joint)

- Stainless Steel Conformend (NBR ISO 5832-1 or similar (ASTM F138));
- Wrought High Nitrogen Stainless Steel (ISO 5832-9);
- Cobalt-Chromium-Molybdenum Casting Alloy (ISO 5832-4) or similar (ASTM F75);

COMPONENTE PATELAR

- Wrought cobalt-chromium-tungsten-nickel alloy (ISO 5832-5);
- Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-iron alloy (ISO 5832-7);
- Wrought cobalt-nickel-chromium-molybdenum-tungsten-iron (1 S O 5832-8);
- Wrought titanium 6-aluminum 4-vanadium (ISO 58 323) or similar (ASTM F136);
- Wrought titanium 6-aluminum 7-niobium alloy (ISO 5832-11);
- Ceramic Materials Based on Alumina (ISO 6474);
- Ceramic Based on Zirconia (ISO 13 356); Wrought cobalt-nickel-chromium-molybdenum alloy (ISO 5832-6);
- Cobalt-chromium-molybdenum casting alloy (ISO 5832-4)

** For implant articulating surfaces of joint replacement, they were considered acceptable in specific applications; however, adequate attention must be given to the project, surface finishing and surface treatment.*

In accordance with standard NBR ISO 21534 - Non-active surgical implants – Joint replacement implants - Particular requirements, the Ultra High Molecular Weight Polyethylene (UHMWPE) is acceptable for the manufacture of implants and can be combined with Cobalt-Chromium-Molybdenum Casting Alloy, raw material of the BKS III - Femoral Component with or without Preservation.

Open Packaging Aseptically

Transport and store it in a cool dry place, with room temperature (max 35 ° C) and relative humidity around 30%
Do not store directly on the floor (minimum height = 20cm) and in very high places, close to the bulbs, which could cause dryness or damage to the packaging label.
Do not store in places where contaminants, such as, cleaning materials, insecticides, pesticides, are stored etc.

	Do not use if package is damaged
	Keep out of the sun
	Keep dry
	Fragile, handle with care

Product Care Before Use

Before using the product, the responsible person designated by the hospital should check the physical condition of the packaging and product. It is important to verify if the external cardboard packaging is wrapped in transparent PVC film and identified with the product information contained in the package. After this verification the transparent PVC film can be removed and the cardboard package can be opened. The product packed in a double blister should be carefully removed from the external packaging in a surgical environment. The blisters should be opened by pulling the medical grade paper. With product in hand, it is necessary to check if it has any grooves, scratches or stains.

If it occurs any adverse effect that affects the user’s security, such as: product is not working, implantable components are damaged, serious problems or death related to these components, the surgeon should report this adverse event to the national health authority and Biomecanica by e-mail: sac@biomecanica.com.br – or phone: (14) 2104-7900 In cases of doubt, the surgeon or health professional can make the communication through the adverse event reporting system on the website of the Health Surveillance ANVISA: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn’t look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the BKS III - Patellar Component explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implants is also under responsibility of the hospital. Biomecanica recommends the implants to be mechanically destroyed with the use of a hammer or press and immediately identified with the sayings “improper for use”. According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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