implantation of metallic component for implantation or due to migration of mechanical components for implantation thrombosis, pseudothrombosis, fracture, paralysis, hematoma and hemorrhage in the blood vessels. pulmonary Embolism, neural or neurologic problems caused by surgical trauma such as gastrointestinal disturbance, premature or late failure of the metallic components for implantation, eschar, surface sensitivity or muscle tenderness in patients with improper tissue coverage of the operated area, difficulties in performing physical activities, allergic reactions and death. In all surgical procedures there is a probability of sickness and death. The surgeon must warn the patient about those incidences before the surgery.

Information to be provided to the Patients

Sayings of this instruction in the articles:
Indications – Contraindication – information of use – Possible Adverse effects – Precautions and Warnings.
 The patient must be informed about the importance of a postoperative follow up. The absence of a follow up prevents the detection of postoperative problems such as the loosens of components or osteolysis and that may result in a progressive loss of the periprosthetic bone stock.
 The patient must be informed about the type of implant to be used, its brand, code, name and lot, the surgeon's name and date of the surgery.

The performance is foreseen in the General Requirements of ANVISA regulations that determine the Essential requirements of efficient safety as well as the safety applicable to the products as described in Resolution - RDC 56, from 06/04/2001. The Cylindric Cages – BM are manufactured in Titanium Alloy according to ASTM F136. They are verified by the laboratory reports of the tests and Certificate of Quality of the Product (Annex I). Those materials are acceptable for the implant manufacturing as established in Annex A of NBR021534 standard.
 The manufacturing and packaging of the Cylindric Cages – BM are based on the specific internal procedures. Cylindric Cages – BM are provided non-sterile and packaged in transparent and nontoxic PVC plastic envelopes and sealed in sealing machine in order to maintain the integrity of the products.
 The Cylindric Cages – BM come with the necessary information for the correct handling of the product as well as the safety of people involved with it.

Admissible combinations with other materials

Metallic alloy implants of distinct chemical composition and from different manufacturers are not recommended due to chemical, physical, biological and functional incompatibility.

Ancillary Components

The related Ancillary Components below must be bought separately because they are not integrated with this product.
 -System for Spine Fixation – Pedicular – Register with ANVISA nº 80128580081.
 - Implant System for Posterior Spine Fixation – Plaque BM – Register with ANVISA nº 80128580083.
 - Instrumental Set for Placement of the Cages – BM – Register with ANVISA nº 80128580080.
 The Bone Graft: The implants may be filled with any bone graft; the choice of the graft is responsibility of the surgeon. The bone graft starts to grow through the existent orifices on the walls of the cages and turns into a solid bone (fusion) which becomes to hold the vertebrae together. This process is known as Intervertebral Fusion.

Accessories

The Cylindric Cages – BM don't have any accessories with the purpose of composing the medical product.

Evaluation of the implanted Product

After implantation, in the intraoperative, the Professional in charge must do the radiological control in order to verify the correct positioning of the product. The Professional in charge must do the clinical and radiologic evaluations after the surgical procedure, in the frequency stipulated by him, to verify the conditions of the implant and the evolution of the bone consolidation. In case the product is out of its correct positioning or presents some unconformity the surgeon must take appropriate corrective actions.

Useful information to Avoid Risks due to Implantation

In order to reduce the risks from implantation one must strictly follow: Contraindications, instruction of use and all information in "instruction of use" of the product.

Intrinsic Risks of Implantation

Cylindric Cages – BM are made of well known biomedical use materials based on the requirements of ASTM F136 – 02 – rule - Standard Specification For Wrought Titanium–aluminum–vanadium Eli (Extra Low Interstitial) Alloy For Surgical Implant Applications (Uns R56401).

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the Cylindric Cages – BM explanted from the patients should be properly discarded by the hospital. The complete characterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "improper for use". According to Resolution nº 2605 from 06/11/2006, it is forbidden to repress implantable devices of any nature classified as single use.

Risks of Contamination

There are risks of biological contamination and transmission of viral diseases such as HIV and Hepatitis because the metallic Cylindric Cages – BM get in contact with the tissue and body fluids. The explanted products must be treated as highly contaminant.

Decontamination, Cleaning and Sterilization

To guarantee the reduction of the risks of infection, the implants (which were not used in surgery but were placed on the operating table or have lost its sterilization) and Surgical instruments must be cleaned after the surgery and before the sterilization.

The cleaning may be performed through manual or mechanic methods.

Equipment for cleaning:

- Ultrasonic Cleaner
- Cleaner - Sterilizer
- Washer - Desinfector

Manual cleaning: It is a manual procedure, where the dirt is removed by physical action, detergent, water, sponge and brush. The manual cleaning is most recommended because it is less aggressive to the surgical implants and instruments. For the manual cleaning, water at room temperature and nylon brush should be used. Wire and abrasive brush should never be used because they may harm the protective layer of the material.

For decontamination, a bactericide antifungal solution of broad-spectrum must be used and aggressive cleaning agents and metal brushes must be avoided in order not to preserve the quality of the metallic implant.

In case of using washing machines, the surgical implants and instruments must be arranged in drawers in a way that they don't touch each other avoiding damage to the material.

People in charge of the cleaning of the surgical implants and instruments should be aware of the type of cleaning product used while attempting to use the concentrations indicated by the manufacturer. The surgical implants and instruments must be carefully rinsed until all residues are washed away. They also must be immediately dried after the cleaning.

Types of Cleaners

Enzymatic Cleaners: The enzymatic cleaners are basically composed of surfactant and solubilizing enzymes. The balanced combination of those elements enables the removal of the organic matter from the material in a short period of time. Enzymatic solutions: The enzymatic solutions present an excellent cleaning action but don't have any bactericidal and bacteriostatic activities. Enzymes: They are substances produced by living cells that guide the chemical reactions of the process. Once the enzymes are made by cells they may be isolated that are still going to keep their catalytic properties since certain conditions are preserved in their development process. The enzymes are classified in three larger functional groups according to the type of Aconitic Acid that they will affect: Protease, Lipase and Amylase that act in the Aconitic Acid such as fat and carbohydrate, which tend to solubilize itself and loose from the articles. Nowadays it is recommended the cleanup of complex configuration article to guarantee the cleaning.

Disinfection:

It is a process that destroys the microorganisms, pathogenic or not, except the ones of bacterial spores, of the articles through physical or chemical means. Disinfection Levels: High Level: It destroys the microorganisms except the ones with a large number of spores => Glutaraldehyde 2% - 20 – 30 minutos. Indication: More likely for Hospital area. Medium level: it eliminates vegetative bacteria. The viruses in their majority are fungus and microbacterias =>1% of sodium Hypochlorite – 30 minutos. Indication: For UBS, Day care center, Asylums and nursing home. Low level: it eliminates the majority of the bacteria, some viruses and fungi; however it doesn't eliminate microbacterias => 0,025% of Sodium Hypochlorite. Indication: Nutrition.

Sterilization

This product is provided not sterile. Before its use it should be sterilized. We recommend the Steam Sterilization in Autoclave at hospital (ISO 11134: 1994

Sterilization of health care products – Requirements for validation and routine control –Industrial moist heat sterilization).

AUTOCLAVE

It is an equipment for the sterilization through saturated steam under pressure.

OPERATING INSTRUCTIONS

For the starting of the equipment, check:
 If the circuit breaker is on; If the water valve is open; If the flush valve is closed.
 1-Open the door of the equipment; 2-Place the material to be sterilized properly; 3-Close the door of the equipment; 4- Select the desirable cycle according to the material to be sterilized; 5-Turn on the Switch; 6- The cycle Will run automatically in the sequence; 7- When turning on the lamp "final of the cycle" open the door partially for about ten (10) minutes for the cooling of the material. The next physical sterilization parameters in autoclave are recommended to be applied (saturated steam);

Cycle	Temperature	Exposure Time
Conventional (1 atm of pressure)	121°C (250°F)	30 minutos
Conventional (1 atm of pressure)	132°C (270°F)	15 minutos
Gravity	132°C (270°F)	45 minutos
High Vacuum	132°C (270°F)	7 minutos

Physical Sterilization Parameters in Autoclave

Note: The moment must be timed when the heat of the sterilization chamber reaches the desired temperature.

FOR CLARIFICATION CONSULT THE INSTRUCTIONS FOR USE THAT FOLLOWS EACH AUTOCLAVE.

Other sterilization method: Besides the Autoclave, the following sterilization method may be used as determined by the hospital:
 -Sterilization by Ethylene Oxide (E.T.O.) – Parameters and Procedures established in the validation protocol EN550:1994 Sterilization of Medical Devices – Validation and Routine Control of Ethylene Oxide Sterilization.

NOTE

The sterilization method, the equipments, the controls and the sterilization instructions are under responsibility of the hospital.
 Care with the Sterilized articles
 Storage Conditions of the Sterilized articles.
 •About the environment: It should be cleaned, ventilated and dried; it should be restricted to the use of the department.
 •About the article: After the sterilization process do not place it on cold surface (Stone or stainless steel), use baskets or perforated containers until it gets colder.
 •Wrapper (cotton fabric, nonwoven fabric, medical grade paper, crepe paper, role film, tyvec or perforated metal boxes) must remain intact and be little handled to prevent it to torn or lose the Seal.
 •Store it in closed cabinets;
 •The shelves must be identified to facilitate the removal of the material;
 •The material must be stored according to the expiring date of its sterilization in order to facilitate the distribution and prevent the storage of expired material; Store the non-sterile materials separately to reduce the level of external contaminants.

The Sending of Material to the Manufacture for Analysis

In case the materials are sent to the manufacturer for analysis, they should be sanitized at Hospital by a bactericidal and antifungal solution of broad spectrum. Then they should be disinfected or steam sterilized in Autoclave or in Ethylene Oxide. They should be taken to Biomechanica in intact packages and the packages should also be identified with the cleaning method, sterilization and data of the product.

After Sales (Customer's complaint)

In case it is necessary to make a complaint about the Cylindric Cages – BM related to some adverse effect that would affect the user's safety, as if the product is not working or the implanted metallic component is damaged, or in case of serious problem or death related to these components, the surgeon in charge must inform the competent Health authorities and Biomechanica about this adverse effect by the e-mail sac@biomechanica.com.br or by phone 0xx14 2104 7900. In case of questions the surgeon or the health care professional in charge can inform about the adverse effect through the Sistema de Notificações em Vigilância Sanitária on the website of ANVISA: <http://www.anvisa.gov.br/hotsite/notvisa/index.htm>

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ENGLISH

Product information and description

The Fusion Cylinders are implants used in several countries due to its recognized application and characteristics in Vertebral Surgery. The Fusion Cylinders are small, cylindrical and hollow instruments and usually made of titanium. These Fusion Cylinders are generally used for substitution of the vertebral bodies and discs. The fusion Cylinders must have inner space to be filled with bone graft between both vertebrae. The implants can be filled with any bone graft; it is the surgeon's responsibility to choose the graft to be used. This bone graft starts growing through the holes of the walls of the cages and forming a solid bone (fusion) which from that moment on keeps the vertebrae united. This process is known as intervertebral Fusion.

The Cylindric Cage - BM is a medical product for implantation. They are devices for substitution of the vertebral bodies and indicated for use in cervical or thoracic spine. It is also used for substitution of a damaged or unstable vertebral body due to trauma or tumor. The Cylindric Cages are indicated to be used as an accessory fixation system, like the cervical plaques. The accessory fixation may be posterior or anterior and it is intended to guarantee the stability of the spine and maintain the compression in the Cylindric Cage. The interior of the Cylindric Cages are filled with bone in order to provide cervical fusion and increase the stability of the system which is designed to give support to the anterior spine even in the absence of fusion for an extended period.

The Cylindric Cages are used with the Spine Fixation System – Pedicular – and the Implant System for Fixation of Posterior Spine – BM Plaques. The Cylindric Cages – BM are available in the following diameters - Ø10, Ø12, Ø14, Ø16, Ø19, Ø22 and Ø25mm.

The implants are finished by mechanical polishing and then go through a pickling process which consists of chemical etching and polished surface of the metal that removes the residual impurities, such as oils, grease and other manufacturing materials, from the manufacturing process. Then, this product goes through an ionization process, in other words, a surface treatment through which it is intended to strengthen the superficial layer of the products made of titanium.

The Cylindric Cages - BM are packed in plastic packages properly identified. The graphic information presented below is only illustrative and refers to the Cylindric Cages - BM.

The Cylindric Cages - BM are manufactured in Titanium alloy based on the specifications of the ASTM F136 standard. The Cylindric Cages - BM are distributed in units according to the table below:

8120-10-010	Cylindric Cage Ø10 by 10mm long, titanium alloy
8120-10-015	Cylindric cageØ10 by 15mm long, titanium alloy
8120-10-020	Cylindric CageØ10 by 20mm long, titanium alloy
8120-10-025	Cylindric CageØ10 by 25mm long, titanium alloy
8120-10-030	