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## ENGLISH

**Technical Name:** Container in General (Boxes, trays, vats, etc)

**Commercial Name:** Boxes and Trays for Transport and Sterilization of the Instruments, Trays Instrumental - Polymeric

**Symbology of the Traceability Tag (Label) according to NBRISO15223 and EN980**

LOT XXXXX Lot REF XXXXX Code EC REP European Representative

Validity Not sterile Manufacturing Date Data of the Manufacturer

## TECHNICAL INFORMATION

**Description of the Product:** The boxes and trays are designed for neatly packing the instrument in order to help the surgeon in the surgical instrumentation as well as in the sterilization of the instruments. The trays are used for settlement of those instruments through a silicon support which allows the instruments to be fixed in a way that even if the trays are moved they remain in place and are prevented from falling and damage. Each tray has a specific layout – a configuration to couple the instruments of determined implant- and the design of each instrument that shows where each of them should be settled in order to guarantee an easy handling in the surgery.

The boxes come in different dimensions that vary according to the number and size of the instruments to be settled for each type of implant. They are rectangular with holes to help in the sterilization process in their inner side. They are used to settle those instruments in a way that even if the trays are moved they remain in place and are prevented from falling and damage.

The boxes and trays are distributed non-sterile and need to be sterilized with the instruments coupled before their use

**Manufacturing Material:** Most of the trays and instruments are manufactured in stainless steel 304 as specified by NBR 13911:1997 which characterizes the chemical composition, mechanical properties and the final behavior of stainless steel. However there are also instrument trays manufactured in Polypropylene – TECAPRO MT brand – which is a low-density resin that offers good equilibrium of thermal, electric and chemical properties and good resistance besides being lighter.

The instrument boxes are manufactured in stainless steel 304 as specified by NBR 13911:1997 which characterizes the chemical and mechanical properties and the final behavior of stainless steel.

**Indications, Precautions, Restrictions/Warnings, Adverse effects and**

**Contraindications:** The boxes and trays are indicated as an auxiliary for the organization of the instruments during surgery besides being used to store, sterilize and transport the instruments in the Operating room. Do not use the product in case it is damaged. The boxes and trays must be carefully handled: during the process of cleaning and sterilization as well as during their transport. Biomecânica Boxes and Trays can only be used to help with the assembly of Biomecânica instruments because they have been manufactured and designed to be used in conjunction. They must be only used for their intended purpose. The boxes and trays must be sterilized before their use. It is necessary to wash and sanitize the boxes and trays before re-sterilizing them. This process must be followed every time they are used in a new surgery. Re-using the trays and boxes before submitting them to the sanitation and re-sterilization process is extremely forbidden. Due to the incompatibility of the materials, design and dimensions it is forbidden the use of components from different manufacturers. The use of components from different manufacturers may impair the quality of the product and bring risks to the patient. In this case the company is disclaimer of any adverse effect that may occur. In case the product is damaged, in bad conditions or if the packaging is damaged, it can not be used and should be immediately discarded.

**Cleaning:** The indispensable cares with the boxes and trays of the surgical instruments start with an effective cleaning. Besides the cares with the cleaning process this instruction mentions some enemies of the boxes and trays such as: blood, surgical residue, stains and discoloration. Water and humidity also harmfully affect the boxes and trays because they dry or are absorbed in the trays and boxes causing undesirable stains. Other enemies present in the washing process are: soap, decolorant, disinfectant and other solutions that are not recommended. For a proper maintenance of the boxes and trays it is important to apply the recommended methods of cleaning and understand the causes of the undesirable effects such as stains. The idea is to assure a proper care for eliminating or extinguishing those causes.

Particularly the trays manufactured in Polypropylene present good resistance to cleaning agents and disinfectants however their maintenance is not exempted from special cares and effective cleaning.

The cleaning process must precede sterilization and sanitation. The organic matter protects the microorganisms and that makes the subsequent steps ineffective. Proper mechanical actions and solutions promote a more effective cleaning. If the article is not properly cleaned the sanitation and sterilization processes will be impaired. When the initial number of contaminant bodies (bioburden) is very high the probability of eliminating such agents in the end of the cleaning process is lower. Check the general conditions of the boxes and trays: if they present porosity, fissure, corrosion, damaged edges, risks or other deterioration. In case of problem suspicion or evidence an identification tag should be pasted on the box and/or tray and they must be separated to be evaluated or discarded. The

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visual inspection is fundamental to verify if the product is damaged, marked or improper to be re-used and guarantee a perfect operation of the product.

**Sterilization of the Boxes and Trays:** The boxes and trays are distributed non-sterile. They must be sterilized before use.

**Recommendations the steam sterilization in Autoclave at the Hospital (ISO 17665-1:** 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices)

**AUTOCCLAVE:** It is a device for sterilization through saturated steam under pressure.

**OPERATING INSTRUCTIONS:** To turning on the device, verify: if the switch is on; if the water valve is open and if the flush valve is closed.

1 – open the door of the equipment; 2 – properly settle the material to be sterilized; 3 – close the door of the equipment; 4 – select the desired cycle according to the material to be sterilized; 5 – Turn on the switch; 6 – the cycle will run automatically in the sequence; 7 – when the light “End of the cycle” turns on, partly open the door for about 10 minutes for cooling down the material.

The following parameters of physical sterilization in autoclave (saturated steam) are recommended to be applied:

Table – Parameters of Physical sterilization in autoclaves		
Cycle	Temperature	Exposure Time
Convencional (1 atm de presión)	121°C (250°F)	30 minutos
Convencional (1 atm de presión)	132°C (270°F)	15 minutos
Gravity	132°C (270°F)	45 minutos
High Vacuum	132°C (270°F)	7 minutos

Note: The time when the heat of the sterilization chamber reaches the desired temperature must be recorded.

**For a better clarification read the instruction manual that comes with the autoclave.**

**Sterilization in ethylene oxide (E.T.O) is another method of sterilization besides the autoclave** - parameters and procedures established in the validation protocol and ISO 11135-1: Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

**NOTE:** It is under responsibility of the hospital the method of sterilization, the equipments, the controls and instructions for sterilization. Care with the sterilized articles.

Storage conditions of the sterilized articles: - as the environment: it must be clean; dry; restricted to people of the sector; - as the article: after the sterilization process do not lay the product on a cold surface (Stone or stainless steel), use baskets or perforated container up to the product cools down; - the wrapping (raw cotton, nonwoven fabric, surgical paper, crepe paper, film paper, tyvec or perforated metallic boxes) should remain intact and be little handled in order to prevent the packs to torn or lose the seal; - it must be stored in closed cabinets with shelves; - the shelves must be identified with the expiration date of the sterilization in order to facilitate the distribution and avoid expired material in the stock; - store the non-sterile products separately from the sterile ones in order to minimize the risks of external contaminants.

**Note:** The polypropylene trays – TECAPRO MT brand – may be sterilized as recommended above since the maximum temperature of 134°C is not exceeded. They can not be submitted to sterilization in incubators and/or other methods of sterilization that requires the temperature to be higher than 140°C because they may be damaged or deformed.

**Storage:** The non-sterile boxes and trays must be stored in a way that the packaging is not damaged. The area of storage should not present impurities, dust or any type of pollutant and the products should not be exposed to extremes temperature and humidity. The flow in this area must be as little as possible. The boxes and trays must be stored in ventilated and dry places, away from the weather conditions, at 35°C of maximum temperature and minimum relative humidity around 30% to 70%.

**Identification:** The boxes and trays come with labels that contain the following information: s Name of the manufacturer; Commercial name; Technical Name; Number of Lot; Number of the register with ANVISA; Product code; Quantity; Description of the product in the packaging (model); Manufacturing date; Expiration date; Description of the raw-material used in the manufacturing of the product (composition); Address of the manufacturer; Name of the technician in charge; the sayings “Non-sterile product”; Vide instructions before use; Warnings/Precautions/Special cares/Symbology; Vide instructions for Use. The boxes and trays come with the following laser engraved information in order to allow the traceability: Biomecánica Logo; Number of the manufacturing lot and code.

**Customer Service:** In case it is necessary to make a complaint related to the use of the boxes and trays and this complaint is associated with some adverse effect that affect the safety of the user, the surgeon in charge must inform the health authorities and Biomecánica about this adverse effect through the email sac@biomecánica.com.br or by the phone 0xx14 2104 7900. In case of doubts the surgeon in charge or the health Professional may notify the adverse event through NOTIVISA on the website of ANVISA: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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