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5315-00-000	Anterior femoral cutting guide	Stainless steel 304
5327-00-000	Tibial cutting guide	Stainless steel 304
5312-00-000	Femoral rotational alignment guide	Stainless steel 304
5363-00-000	Tibial retractor	Stainless steel 304
5362-00-000	Patellar tendon retractor	Stainless steel 304
5317-00-000	Distal femoral cutting guide	Stainless steel 304
5314-00-000	Impactor	Stainless steel 304
5359-00-000	Universal extractor	Stainless steel 304
5361-00-000	Hohmann retractor (24.5x252 mm)	Stainless steel 304
5360-00-000	Hohmann retractor (24.5x252 mm)	Stainless steel 304
5367-00-000	Screw extractor	Stainless steel 420
5228-00-000	Osteotoma 20 mm	Stainless steel 420
5307-00-000	Quick coupling	Stainless steel 304
5372-00-000	Screw pointer	Stainless steel 304
5352-00-000	Obial ressection guide	Stainless steel 420
5211-00-000	Osteotomo 22 mm	Stainless steel 420
6003-02-032	Drill bit Ø3,2x150 mm with quick coupling	Stainless steel 420
5381-00-000	Patellar guide 3 grips	Stainless steel 304
5383-00-000	Femur test extractor	Stainless steel 304
5382-00-000	Paquimeter	Stainless steel 304
5381-01-000	Patellar drill Ø7,0 mm	Stainless steel 420
5384-00-000	Patelar Clamp	Stainless steel 304
5385-00-000	Extramsudador tibial cutting guide	Stainless steel
5386-00-000	Plateau inserter	Stainless steel 420
5313-12-000	Bar for Femoral alignment guide 120mm	Stainless steel 304
5313-15-000	Bar for Femoral alignment guide 150mm	Stainless steel 304
5313-20-000	Bar for Femoral alignment guide 200mm	Stainless steel 304
5324-10-000	Intercondylar guide P/M	Stainless steel 304
5353-28-000	Patellar trial Ø28 mm (3grips)	Polyphenylene
5353-30-000	Patellar trial Ø30 mm (3grips)	Polyphenylene
5353-32-000	Patellar trial Ø32 mm (3 grips)	Polyphenylene
5353-34-000	Patellar trial Ø34 mm (3 grips)	Polyphenylene
5353-36-000	Patellar trial Ø36 mm (3 grips)	Polyphenylene
5211-16-000	Osteotomo 16mm	Stainless steel 420
5211-25-000	Osteotomo 25mm	Stainless steel 420

Indications, contraindications, adverse effects and warnings

Indications: The accessories for femoral components – BM have the purpose of correcting the asymmetric segment defects of the bone as well as giving more stability to the femoral component without preservation – BM designed for revision surgeries. They are used in revision surgeries for restitution or replacement of the femorotibial and patellofemoral joints (knee) in conjunction with their ancillary components. That procedure is known as “total knee arthroplasty” and tries to create an artificial articulation with the same functionality of the natural knee. The total knee arthroplasty has the purpose of eliminating the pain symptoms, correcting the defects, allowing a functional range of motion and stabilizing the knee.

Contraindications

The contraindications below do not limit the use of the accessories for femoral components – BM, however they should be analyzed by the surgeon: Serious disability of the quadriceps; infection or systemic disease or close to the joints; Arthrosis of the knee after the arthroplasty, previous patellectomy, neuropathic arthritis; serious discontinuity or dysfunction of the extensor mechanism; Recurvatum deformity due to muscle weakness and knee arthrosis; obesity, fever or leukocytosis; the patient is unable to receive anesthesia; rheumatic disease in the leg; skin problems such as psoriasis; urinary tract infection; osteomyelitis history close to the knee; intense neuromuscular, muscular or vascular deficiency on the affected region; serious injury of the bone structure which generates risk and becomes inviable for the application of prostheses; Bone tumor close to where the component was implanted; muscle atrophy; Bone immaturity; insufficiency of the extensor mechanism of the knee; failure of the structural soft tissues of the knee; Knee ankylosis; insufficient bone structure in the tibia, femur and patella; Patients is unable or unwilling to follow the cares and instructions for the postoperative.

Adverse Effects

Migration of the articular line; change in the human gait pattern; Decrease or increase of the patellar height; hematomas; wound dehiscence; progressive assymptomatic bone reabsorption (osteolysis) located towards the component; looseness, migration or fracture of the medical device; surface or deep infection; vascular disorders including thrombosis and pulmonary embolism. Limb shortening due to bone reabsorption; allergic reaction to foreign bodies that may lead to histological reactions involving several types of macrophages and fibroblasts; pain, discomfort or abnormal sensations due to the presence or migration of the medical device and due to the surgical procedure; brain or neurological damage due to surgical trauma (that includes paralysis and injury of the soft tissues); incapacity to perform regular daily activities; per prosthetic fractures.

Warnings

The accessories for femoral components – BM help in the stabilization of the femoral component without preservation – BM which has the purpose of replacing the damaged articulation and not the healthy structures of the body. The surgical implants should never be reused and the explanted ones should never be implanted again. Small defects and inner pressure patterns should be present even if the implant looks intact and that may cause the failure or breakage of the implant. The proper selection of the models of the accessories for femoral components – BM is extremely important and it should be done accordingly to the bone characteristics of the patient, the analysis of the surgeon and to the technique adopted by the surgeon. The generated infections are usually linked to comorbidity. The aseptic failure of the total knee arthroplasty may be caused by the component looseness, wear of the ancillary component manufactured in polyethylene, ligament looseness, periosseotic fracture, arthrofibrosis and patellofemoral complications. The results of the inferior surgeries are inferior compared to the result of the primary surgeries. Smokers have a greater chance of having longer operative times and higher costs. The replacement of the knee has a finite useful life which may be adversely affected by the activity level.

The patient should be informed as:

Need a periodic medical follow-up for diagnosis of potential changes in the conditions of the implant and adjacent bones. Only a follow-up is already able to detect potential looseness of the component or osteolysis. The failure to perform a revision surgery in case of component looseness or osteolysis may result in progressive bone stock. The patient should be present even if the implant looks intact and that may cause the failure or breakage of the implant. The proper selection of the models of the accessories for femoral components – BM is extremely important and it should be done accordingly to the bone characteristics of the patient, the analysis of the surgeon and to the technique adopted by the surgeon. The generated infections are usually linked to comorbidity. The aseptic failure of the total knee arthroplasty may be caused by the component looseness, wear of the ancillary component manufactured in polyethylene, ligament looseness, periosseotic fracture, arthrofibrosis and patellofemoral complications. The results of the inferior surgeries are inferior compared to the result of the primary surgeries. Smokers have a greater chance of having longer operative times and higher costs. The replacement of the knee has a finite useful life which may be adversely affected by the activity level.

Operations in conjunction with other medical devices: According to “NBR ISO 21534 standard– Non-active surgical implants– joint replacement implants – Particular requirements, the wrought cobalt-chromium-molybdenum alloy which is the raw-material of the accessories for femoral components – BM is acceptable for manufacturing of the implants. It has been elaborated a list of the acceptable combinations with the femoral components – BM according to “NBR ISO 21534 – Non-active surgical implants– joint replacement implants– Particular requirements. Independently of that list it is also necessary that the combinations are between same brand implants because they have been designed for such combinations based on the surface finishing and treatment and other factors required in the project that may interfere in the combinations. Therefore the use of metallic implants from different manufacturer is not recommended due to chemical, physical, biological and functional incompatibility.

Admissible contact alloy (when there is no articulation)

- Titanium alloy (NBR ISO 5832-3, NBR ISO 5832-11, ASTM F136);
- Cobalt alloy (NBR ISO 5832-4, NBR ISO 5832-5, NBR ISO 5832-6, NBR ISO 5832-7, NBR ISO 5832-8, NBR ISO 5832-12, ASTM F75)
- Stainless steel (NBR ISO 5832-9)

Special conditions for storage, conservation and/or handling of the product (described in the external protective packaging of the product).

Open the Packaging Aseptically

Temperature Limits

Transport and store it in dry and cool places and at room temperature (maximum 35°C) and Relative humidity (around 30% and 70%.

Do not store it in direct contact to the floor (minimum height = 20cm) and neither in very high places and close to lamps because that may lead to the dryness and damage to the label.

Do not store it in places where contaminants substances are stored such as cleaning materials, insecticide, pesticide, etc.

Attention! Instructions, warnings and precautions: Ver Instrucción de Uso

Do not use it if the packaging is damaged

Keep it dry

Fragile, carefully handle it

Keep it away from the Sun

Vide Instruction For Use

with all information related to the material present in the package. After the checking the responsible person should remove the PVC transparent film that covers the outer package and open the cardboard packaging. The product packed in double blister should be carefully removed from the outer package in a clean environment. The first blister should be open by pulling the surgical paper from the edge of the blister. The second blister should also be open by pulling the surgical paper from the edge of the blister. With the product in hands the surgeon should check if it does not present any notches, scratches or marks.

Load supportable by the system associated with the weight of the patient

Some tests showed that the Tibial components – BM for revision may be used with no restrictions in patients weighing up to 106 Kg and in cases of patients weighing over 80Kg activities such as getting up the stairs is not recommended.

Post-sales (Client’s complaint): in case it is necessary to make a complaint of the accessories for femoral components – BM related to some adverse event that may affect the safety of the user such as if the product is not working, if the metallic component is damaged or in cases of serious problems or death related to the components the responsible surgeon should inform the competent health authorities and Biomecânica about that event through the email: tecnovigilancia@biomecânica.com.br or by the phone 0xx14 2104 7900. In case of doubts the responsible surgeon or health Professional may inform about the adverse event through **NOTIVISA** on the website of ANVISA: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn’t look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the accessories for femoral components – BM explained from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecânica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings “Improper for use”. According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Contamination Risks

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

Patient’s card: A card that comes with the Modular heads and Femoral stems. It should be handed to the patient together with the identification tags of the product for arthroplasty. With the card and tags in hands the patient may obtain information related to the quality of the implant. Further information is available on www.biomecânica.com.br/paciente

Send the material to the manufacturer for analysis: In case of sending the material to the manufacturer for analysis that material should be sanitized at hospital with the help of a broad spectrum bactericide and antifungal solution. Then it should be disinfected steam sterilized in autoclave or ethylene oxide. It should be sent to Biomecânica in intact packaging and identified with the method of cleaning and sterilization as well as with the data of the patient.

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