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STRAIGHT PLATES SMALL AND LARGE FRAGMENTS - FREE-BLOCK

An orthopedic implant can only be implanted in one single patient, for one single time, through it doesn't risk damage, previous stress may create imperfections that may impair the success of the implant. The improper selection of the implant may cause uncommon stress on the implant and may result in a subsequent fracture of the implant. Due to its presentation and characteristics the product does not present any side effect. However, it is suggested to monitor the patient's re-establishment through x-ray method, in order to avoid magnetic resonance because the product in stainless steel causes interference in the exam and impairs the diagnosis. It may also cause the detachment of the implant and that may bring problems to the patient. The use of rigid fixation systems of straight plates for large and small fragments – BM are designed to be an auxiliary element for the bone fixation and stabilization in the treatment of fracture of upper and lower limbs and not to replace the healthy structures of the human skeleton. The surgeon should be aware of the chance of a patient's physical development which may ask for a revision of the implant. The surgeon should be familiar with and have enough knowledge of osteosynthesis and its limitations including the pre and post-operative techniques, precautions and risks. Attention should be paid to the bone preparation for the insertion of the rigid fixation system of straight plates for large and small fragments – BM because the stability and the success of the fixation may be compromised. During the handling of the implant, scratches, grooves or any other event that may damage or mark the implant should be avoided because such defects, such as stress nucleations areas may be sites of nucleations and cracks impairing the resistance to corrosion, and that may result in the fracture of the implant. In certain situations it is needed to mold the implant in order to perfectly adapt it to the bone anatomy it is treated. In case the modeling of the implant is needed the bone anatomy itself is the limit for that modeling. In order to guarantee a proper implantation it should be used the instruments manufactured by Biomecnica being they were designed and manufactured specifically to be used in conjunction with those implants. The use of instruments from different manufacturers may compromise the surgery besides being in discordance with the register of the product with ANVISA.

Do not use the product in case it is damaged: NON-STERILE PRODUCT – STERILIZE BEFORE USE; SINGLE USE PRODUCT – THE REUSE OF THIS PRODUCT IS EXTREMELY PROHIBITED;

Special storage conditions: Keep it in ventilated and dry places away from the light and weather conditions. Do not use the product in case the packaging is damaged. Note: components of implants from different manufacturers should not be used; therefore it is recommended that the products have the same provenance. Manufacturing date, expiration date and product lot: VIDEO LABEL.

Adverse effects: the mechanical looseness may be the result of a defective fixation, instable reconstruction of the fractured bone or occult infection. Metal sensitive reactions have been rarely informed. The implantation of a foreign material in the tissues results in histological reactions that involve several sizes of macrophages and fibroblasts. The chemical importance of that effect is uncertain as similar changes may occur as a precursor or during the healing process.

Information to be given to the patient: it is available in the items: indications-contraindications-information for use-potential adverse effects-precautions and warnings. The patient should be informed as the importance of a postoperative follow up. The failure of a postoperative follow up prevents the detection of postoperative problems such as the component loosens or osteolysis. The lack of a revision surgery in case of implant loosens or osteolysis may result in a progressive loss of the periprosthetic bone stock. The patient should be informed as the type of implant to be implanted, the brand, code, name, lot, surgeon's name and date of the surgery.

Evaluations of the implanted device: after implantation, in the intraoperative, the person in charge should perform a radiological control in order to check the correct positioning of the device. The Professional in charge should do the clinical and radiological evaluations after the surgical procedure at a frequency stipulated by him in order to check the conditions of the implant and the evolution of the bone healing. In case the product is out of the correct positioning or presents some non-conformity it is the surgeon's responsibility to take the most appropriate corrective action.

Useful information to avoid the risks resulted from implantation: In order to minimize the risks due to implantation it should be strictly informed: contraindications, instruction for use and all information available in the "instruction for use" of the product.

Intrinsic risks of implantation: the rigid fixation system of straight plates for large and small fragments – BM are manufactured with well known biomedical use materials which follow the requirements of ASTM F138 Standard -Standard Specification For Wrought 18 Chromium-14 Nickel-2.5 MolybdenumStainless Steel Bar And Wire For Surgical Implants (UNS S31673) for the manufacturing of the plates and screws and ASTM F136 standard -Standard Specification For Wrought Titanium-6Aluminum-4Vanadium EL Extra Low Interstitial Alloy For Surgical Implant Applications (UNS R56401) for the manufacturing of the rings.

Risks of contamination: There are risks of biological contamination and transmission of viral diseases such as HIV and Hepatitis because the metallic components of the rigid fixation system of straight plates for large and small fragments – BM get in contact with the bone tissues and fluids. The explanted devices should be treated as highly contaminants.

Decontamination, cleaning and sterilization: The product is distributed in its original packaging and previously clean ready to be used. It should be carefully handled in order to avoid contamination. In case of contact with the product, it is recommended to compromise its cleaning state it should be clean before its sterilization. Some cleaning procedures are suggested: The cleaning may be through manual or mechanical methods.

Equipments for Cleaning:
-Ultrasound cleaner
-Washer – sterilizer
-Washer – disinfectant

Manual cleaning: it is a manual procedure where the dirt is removed through physical action with the help of detergent, water and artifacts such as sponge and brush. The manual cleaning of the implants and instruments is indicated in the manual of the instruments less. For the manual cleaning it should be used water at room temperature, nylon sponge and never abrasive or steel ones because they may harm the protective layer of the material.

For the decontamination a bactericide and antifungal solution of broad spectrum should be used in order to aggressively cleaning agents and metal brushes should be avoided so that the quality of the metallic implants is not impaired.

In case of cleaning machines the surgical implants and instruments should be settled in drawers in a way that they don't touch each other in order to avoid damage to the material.

People responsible for the implant and surgical instrument cleaning should be aware of the correct use of the product used while respecting the concentrations indicated by the manufacturer. The implants and surgical instruments should be immediately dried after the cleaning and carefully rinsed up to all residues are washed off.

Types of cleaners
Enzymatic cleaners: they are basically composed of enzymes, surfactants and solubilizers. The enzymatic combination of those elements allows the removal of the organic matter from the material in a short period of time.

Enzymatic solutions: they present excellent cleaning action however they do not present any bactericide or bacteriostatic activity.
Enzymes: they are substances produced by live cells that govern the chemical reactions of the process. Once they are produced by the cells, the enzymes may be isolated by the way they will stop their catalytic properties if certain conditions are kept during their manufacturer. The enzymes are classified in three major functional groups depending on the substrate type that they will affect: proteases, lipases and amylases that act on protein substrates, fats and carbohydrates which tend to solubilize and detach from the articles. Currently it is recommended the cleaning of complex articles for the cleaning guarantee.

Disinfection: It is a process that destroys the microorganisms – pathogenic or not- from the articles except for the bactericide spores through physical and chemical means.

Disinfection levels:
High level: it destroys all microorganisms except for the high number of spores -> Glutaraldehyde 2%-20-30 minutes.
Indication: Rafter hospital.

Medium level: It eliminates vegetative bacteria, most of the viruses, fungi and microbacterium ->Sodium Hypochlorite 1%-30minutes
Indication: for UBS, Day cares almshouses and nursing homes.
Low level: it eliminates most of the bacteria, some viruses and fungi but it does not eliminates microbacterium -> Sodium Hypochlorite 0,025%
Indications: Nutrition.

Sterilization: This product is provided non-sterile. It should be sterilized before use. It is recommended the steam sterilization in Autoclave at ISO 11134:1994 Sterilization of health care products – Requirements for validation and routine control – industrial moist heat sterilization).

AUTOCULVE: It is a device for sterilization through saturated steam under pressure.

sequence: 7-when turning on the lamp "final of the cycle", partially open the door for about ten (10) minutes for cooling the material.
The physical sterilization parameters in autoclave (saturated steam) specified in the table below are recommended:

Cycle	Temperature	Exposition time
Conventional (1 atm of pressure)	121°C (250°F)	30 minutos
Conventional (1 atm of pressure)	132°C (270°F)	15 minutos
Gravity	132°C (270°F)	45 minutos
High Vacuum	132°C (270°F)	07 minutos

Note: The time when the heat of the sterilization chamber reaches the desired temperature should be recorded.

FOR FURTHER INFORMATION CHECK THE INSTRUCTION MANUAL THAT COMES WITH EACH AUTOCLAVE.

Other sterilization method that may be applied besides the autoclave is the following as defined by the hospital:
- Ethylene Oxide sterilization (E.T.O.) – parameters and procedures established in the validation protocol and in ENISO-1994 Sterilization of Medical Devices – Validation and Routine Control of Ethylene Oxide Sterilization.

NOTE
It is under responsibility of the hospital the sterilization method, the equipments, the controls and the sterilization instructions used.
Cares with and storage conditions of the sterilized articles. As the environment: it should be clean, ventilated and dry, it should be restricted to the staff of the sector. As the article: after the sterilization process do not place the article on cold surfaces (Stone or stainless steel), use perforated baskets, avoid hanging articles; coil down; the wrapping (raw-cotton fabric, nonwoven fabric, surgical paper, crepe paper, film paper, tyvek or perforated metallic boxes)should remain intact and be little handled in order to avoid the packaging to tear or to lose the seal; it should be stored in close cabinets on shelves, the shelves should be identified in order to facilitate the removal of the material; the material should be stored according to the expiration date of the sterilization in order to facilitate the distribution and to prevent having expired material in the inventory; do not store sterile materials with non-sterile materials in order to reduce the level of external contaminants.

Product discard: No surgical implant should be re-used. Any metallic implant should be discarded after its use. Even though the implant doesn't look damaged it may present little defects that may be due to contamination or its leakage due to fatigue. This way the Rigid Fixation System for Straight Plates for Large and Small fragments – BM explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecnica recommends the use of a mechanical destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Traceability: The implantable device comes with to traceability tags with the data of the used implant. That tag is equal to the one of the outer labeling of the product. The information available in that traceability tag is related to the item "Presentation" of the technical report. One of the tags should be pasted on the patient's report and another one is for the control of the hospital. The hospital should register that information in its files and forward it to the patient. We inform that the patient's report must have all pertinent information to the traceability of the implanted product registered. By pasting that tag on the patient's report the traceability of the used product is allowed. Among the major information we highlight as indispensable the name of the implant manufacturer, the name of the implanted device, the code of the device, the lot number and the number of the product register at ANVISA. That information is described in the traceability tags that come with the product and in the outer labeling. The following information should be also available as a displacement such as the date of the phone box 12,104.7906, in case the patient had the device implanted, the name of the surgeon, the patient's weight, the age of the patient, and other information required in the patient's report should also be filled out.

Marking: The rigid fixation system of straight plates for large and small fragments – BM presents the following laser engraving information to allow the traceability and patient follow up in the postoperative according to the space available in each product: Straight plates with a locking system: Biomecnica logo; Lot number; acronym of the material (stainless steel) used for the manufacturing of the material (S1); Code; Dimension and CE mark. The washer is laser engraved with the following: Biomecnica logo; lot number; acronym of the dimension and the manufacturing of the product (S1). Dimension and CE mark. The locking rings and sleeves is laser engraved with the following information due to the little space: Biomecnica logo; lot number and acronym of the material (Titanium) used for the manufacturing of the product (T3). Due to the little space, the screws are laser engraved with the following information: Biomecnica logo; lot number; acronym of the material (S1) and dimension.

Post-sale (Client complaint): In case it is necessary to make a complaint about some adverse effect that affects the safety of the user such as if the product is not working, or if the product is damaged, or in case of serious problems or death related to those components, the responsible surgeon should inform that adverse event to the competent health authorities. Biomecnica through the email hotline@biomecnica.com.br or by the phone box 12,104.7906, in case the responsible surgeon or health professional may communicate that adverse event to NOTIVISA through ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Send of the material to the manufacturer for analysis: In case of sending the implant to the manufacturer for analysis it should be sanitized at hospital with the help of a bactericide and antifungal solution of broad spectrum. Then the product should be disinfected or steam sterilized in autoclave or in ethylene oxide. It should be sent to Biomecnica in intact packages and identified with the cleaning method, sterilization and date of the product.
Description of the efficiency and safety of the medical product in accordance with the regulations of ANVISA that rules the essential requisites of efficiency and safety of the medical product: The risks associated with the rigid fixation system of straight plates for large and small fragments – BM are in accordance with the essential requisites of the medical products according to the regulation of the ANVISA for the products classified in rule number class 3. The risk management of the product and process are in accordance with NBRISO14971 as foreseen in the procedure - PRC002 - Risk analysis of the medical devices.

Locking system of the ring and stability of the system
Locking system of the ring: The locking of the ring in the plate happens through interference. The ring used in the hole of the plate is internally conical and has a lateral slit to allow the expansion and locking of that ring in the plate when it receives the FB screw. The ring has a conical thread in its interior that fits the head of the FB screw. The ring is externally convex and therefore it fits perfectly into the concave hole of the plate and that allows adjustment movement according to the direction of the fracture. Besides, the area of the hole of the plates is large enough to allow the locking of the screw in the ring in any direction. That movement is only possible before it is finally tightened in the system. After the tightening the locking is no longer possible due to the locking feature of the system. In order guarantee the moldability of the ring generated during the fixation of the screw, the ring is manufactured in titanium alloy ASTM F136. The titanium presents high moldability due to its small size. The use of screws are the ones specifically used by its extreme chemical passivity and consequently excellent biocompatibility and also because it has physical properties suitable for a good biomechanical behavior in a long term.
Stability of the system: The verification of locking of that system was given through a pull-out test using an assembled plate and rings in the polyurethane block. The polyurethane block used for the pull-out test presents characteristics similar to the bone which simulates situations as similar as possible to the screw fixed to the bone. The used block was the 40° one following the recommendations of ASTM F1839 standard. The plates were fixed in three parts of the polyurethane block. Only the small fragment plate had a 4 point fixation due to its small size. The use of screws are the ones specifically used for the tested plates. The issue of the project was to know what the pull-out load the assembled system would be able to support before disassembling. The screw pull-out as a common failure usually caused by a defective fixation of the screw to the bone. After the tests were finished it was verified that the assembled system showed an excellent fixation. The plates were removed without any signs of the screw pull-out. At the end of the test it was possible to observe the quality of the fixations- screw-bone and plate-screw because the system remained intact with no pull-out beginning. Through the pull-out test it is possible to verify that the dimensions stipulated in the design of the holes of the plates and the locking ring are in accordance with the intended ones because after the test the system remained assembled. Therefore the pull-out test showed that the locking ring fulfilled its role of locking the screw in the plate in such a way that after 01 end of the test, with the deformation of the plate, the system remained assembled. In other Word, an absolute stabilization was achieved through the locking of the threaded screw into the locking ring of the plate.

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