CIMENTO ORTOPÉDICO VERTEBROPLASTIA

REF 2010

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REF 2010

PRODUCT DESCRIPTION

Orthopedic Cement for Spine Surgery is an acrylic composition, auto-polymerized, obtained by a powder of Polymethyl-Methacrylate and Sulfate of Barium with a polymerization initiator, Peroxide of Benzoyl; and by a liquid component constituted by a monomer Methyl methacrylate with NN Dimethyl-p-toluidine activator. When the monomer and Polymethyl-Methacrylate are in contact (mixed conveniently) to NN Dimethyl-p-toluidine activates Peroxide of Benzoyl which begins the polymerization. The process is progressive, and the product liquid in the initial stage is going turning pasty until reach the solid state. The reaction is exothermic and with maximum temperatures that can vary depending on the volume and dissipation of heat, could reach 110°C

PRESENTATION FORM

01 liquid component (20ml):

Colorless Monomer, inflammable filled in glass flask.

Composition:

Methyl Methacrylate	98,215% w/w
NN Dimethyl-p-toluidine*	0,816% w/w
Alcool Ethylic	0,945% w/w
Acid Ascorbic	0,022% w/w
Hydroquinone**	0,002% w/w

** Hydroquinone is added to prevent premature polymerization which can happen under certain conditions, such as: high temperatures and exposition to the light (see pharmaceutical precautions).

* N, N-Dimethyl-p-toluidine is added to promote the cold cure of the cement

OBS: Supplied sterile in double blister.

01 powder component (29g):

Composition:

Polymethyl Methacrylate	68,96% w/w
Barium Sulfate***	29,32% w/w

ORTHOPEDIC CEMENT SPINE SURGERY

REF 2010

Peroxide of Benzoyl

1,72%w/w

*** Barium Sulfate supplies the radiopacity characteristic to the orthopedic cement.

OBS: Supplied sterile in double blister.

INDICATIONS

They are suitable to Spine Surgery, patients with osteolytic lesion of vertebral body resultant of metastases, myeloma, and vertebral hemangioma and for osteoporosis lesion, in agreement with the following criteria:

- Hemangioma of aggressive clinical behavior and/or radiologic.
- Severe Pain and incapacitating, resistant to the analgesic use.
- Progression of deformity and increase of the pain for osteolytic lesions of vertebral body.
- Vertebral body cannot totally be destroyed, should possess at least 1/3 of its original height.

PREPARATION AND APPLICATION


Before the mixture

Store the product in atmosphere with temperature close to 23°C, during 3 to 4 hours. (In pharmaceutical precautions, are described important details on the storage process).

SPECIAL CARES

The times of catalysis can suffer significant alterations when components are in temperature different from 23°C.

APPROXIMATE MEDIUM TIMES AND POLYMERIZATION PHASES IN 23°C ORTHOPEDIC CEMENT FOR SPINE SURGERY

CONTRAINDICATIONS

The Procedure is contraindicated:

- In the presence of infection
- In the disturbances of non compensated coagulation.
- Destruction of the subsequent wall of the vertebral body should be evaluated, because the potential of extravasation of the injection material for the space of the rack channel, it can produce compression of structures nervous.
- Where there is hypersensitivity to any one of the components of the formula.

INTERACTION WITH OTHER MEDICINES

Unknown.

WARNINGS

In spite of the low level of complications of the percutaneous vertebroplastia, it should be attentive to the following potential risks:

- Bleeding in the place of the puncture usually solved with simple compression of the place.
- Temporary worsens of the pain and fever as reaction of the temperature generated by the polymerization of the biomaterial.
- Infection.
- Lesion of roots or marrow for inadequate extravasation of contrast, with potential risk of paralysis and /or radiculopathy.
- Fistula liquor
- Passage of the biomaterial to the veined systems with lung clot.

ADVERSE EFFECTS

Listed below are the most serious adverse effects, which can happen with the use of the bone cement. The surgeons should be conscious of those reactions and prepared to treat them when they happen.

SERIOUS:

- Infarct of the myocardium, 2) Accident brain-vascular, 3) heart failure, 4) sudden death, and 5) lung Embolism.

MORE FREQUENT:

- Transitory fall of the sanguine pressure, 2) Tromboflebite, 3) Hemorrhage and hematoma, 4) Superficial Infection or deep lesion, 5) Irregularity in the short-term heart conduction.

PRECAUTIONS:

- The use of the Orthopedic Cement requests collaboration and consult between surgeon and anesthetist.
- The anesthetist should be informed, during the operation, the placement of cement.
- The monomer can cause hypotension; and this fall blood pressure can produce heart arrhythmias or a myocardium ischemia .
- The hypotensive effects of the methyl methacrylate can be potentialized if patient suffers of hypoxemia.
- The preparation of the cavity of the bone marrow, results in medullar contents entering in the sanguine current. Enlarging a hole of the medullar cavity can had similar effects on the blood pressure as well the introduction of the Orthopedic Cement
- The implantation of a strange body in the tissues increases the normal risk of infection after surgical interventions.
- The routine use of appropriate prophylactic antibiotic can be desirable in the doctor's prescription.

USER PRECAUTIONS

The components liquid and powder of the Orthopedic Cement is to be used in **unique dose**. Always add completely the liquid component to the powder component while manipulating the mixture of the Orthopedic Cement.

The liquid component is volatile. Ventilate the surgical center adequately. The continuous inhalation of the vapor of the liquid can have a soporific effect or also cause irritation of the breathing area and eyes therefore should be avoided.

Liquid component is a powerful solvent and can cause contact dermatitis in sensitive people. The possibility of hypersensitive reactions reduces, using a second pair of gloves and following the preparation instructions strictly.

If the liquid component has direct contact with the eyes, wash them with abundant water.

Concentrate vapors of the liquid component can have adverse reaction with soft contact lenses.

PHARMACEUTICAL PRECAUTIONS

Store the product in a dry place, airy and protected of the light, where the temperature is below 25°C. Otherwise, it can have polymerization of the liquid component spontaneously.

Special cases of polymerization of the liquid component: (1) lack of stabilizer; (2) in constant exposition to the solar light, direct or indirectly; and (3) in contact or exposition to the heat (consider additional cares, since the product is inflammable).

STERILIZATION

Final Packing - ETO - Validity 3 years.

Powder component - Gama Ray

Liquid component - Filtration

This product should **not be re-sterilized**. If the surgical packing is damaged, he should be discarded.

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