# ACETÁBULOS

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</

(expired).

**Disclaimer:** Biomecnica, as a manufacturer of this device, does not practice medicine or recommend this or any other surgical technique for use in a particular patient. The surgeon performing any implant procedure is responsible for defining and using the proper techniques for implanting the medical device in each patient. Biomecnica is not responsible for the selection of the appropriate surgical technique to be used in each patient. The surgical technique depends on the technique used by the surgeon, who is also responsible for choosing the method, type and size of products to be placed.

**Ancillary components:** The ancillary components listed below must be purchased separately:

- Biomecnica Metal Femoral Head
- Cemented Modular Stem
- Acrillex Orthopedic Cement

**Accessories:** Biomecnica Polymeric Cemented Monocomponent Acetabular Cups do not have any accessories with the purpose of integrating the medical product.

**Allowable Combination with Other Medical Products**
The table below "allowable combinations with other materials" was based on the standard "NBR ISO 21534 - Non-active surgical implants - Joint replacement implants". Particular requirements of the alloys listed in the table "allowable alloys", the implants must also be made by the same manufacturer because they were designed for such combinations observing surface finishing, surface treatment and other factors required on projects that may interfere with the combinations. Therefore, metallic implants from different manufacturers must not be used together because of chemical, physical, biological and functional incompatibility.

**Allowable Contact Alloys**

Implant	Allowable contact alloys* articulation
Monocomponent Acetabular Cup	Wrought Stainless steel (ASTM F138, Cemented Polietileno (NBRISO5832-2)) <p>Wrought high nitrogen stainless steel (NBRISO5832-9) Cobalt-chromium-molybdenum casting alloy (ASTM F75, NBRISO 5832-4) Wrought cobalt-chromium-tungsten-nickel alloy (NBRISO 5832-5) Forgeable and cold-formed cobalt-chromium-nickel-molybdenum-tungsten-iron alloy (NBRISO5832-7) Wrought cobalt-chromium-nickel-chromium-molybdenum-tungsten-iron (NBRISO 5832-8) Wrought titanium 6-aluminum 4-vanadium (ASTM F136, NBRISO5832-3) Wrought titanium 6-aluminum 7-niobium alloy (NBRISO 5832-11) Ceramic Materials Based on Alumina (NBRISO 6474) Ceramics Based on Zirconia (NBRISO 13356)</p>

\* For articulating surfaces, they were considered acceptable in specific applications; however, degrading surface must be given to the project, surface finishing and surface treatment.

**Warnings, Restrictions and Precautions**
**Precautions: PRODUCT FOR SINGLE-USE. DESTROY AFTER REMOVAL. DO NOT REUSE THE PRODUCT. PROHIBITED TO REPROCESS.**

An orthopedic implant must not be used again. Even if they appear undamaged, they may have small defects and stress patterns which may lead to fatigue failure. Check the sterilization validity and if the product packaging is intact. Open the packaging aseptically, after certifying that the dimension is the one chosen. Sterilization validity: 5 years (indicated in internal and external product packaging). The surgeon should be aware of the physical development of the patient, a revision of the implant may be needed.

**Contraindications:** Total hip arthroplasty is a major procedure and it is associated with a significant number of complications and a mortality rate of 1% to 2%. Thus, when indicating a total hip replacement, patients should be carefully evaluated especially those with systemic diseases. Cardiac, pulmonary, hepatic, genitourinary or metabolic diseases, hypertension and malignancy should be treated and controlled before total hip arthroplasty. Other situations are also contraindicated in case of total hip replacements, they are: Osteomyelitis in the upper femur, pyogenic infection of the hip joint, or evident infections are absolute contraindications. Every effort should be made to guarantee the asepis preoperatively in patients who have one or more of the following abnormalities: fever, local inflammation, rapid joint destruction, apparent bone absorption shown in roentgenogram; unexplained high sedimentation rate through other diseases, increase of the WBC rate, or suspected striking changes in the count. However, an infection recurrence has a higher probability to occur in revisions than in primary surgeres. Under no circumstances it is recommended that the insertion is done where there has been extensive infection involving the bone; or when involved soft tissue cannot be removed; or where the organism is gram negative, antibiotic resistant Staphylococcus aureus; or where highly virulent or mixed organisms are involved. Additionally, in distant site infection, such as genitourinary, lung, skin or in other organs, it is relatively contraindicated due to transmission of the infection to the hip joint site. The focus of infection must be treated before, during and after the implantation. The use of this implant is contraindicated when there is loss of musculature or neuromuscular disease in the affected limb, which would make the procedure unjustifiable.

**Contraindications:** Total hip arthroplasty is a major procedure and it is associated with a significant number of complications and a mortality rate of 1% to 2%. Thus, when indicating a total hip replacement, patients should be carefully evaluated especially those with systemic diseases. Cardiac, pulmonary, hepatic, genitourinary or metabolic diseases, hypertension and malignancy should be treated and controlled before total hip arthroplasty. Other situations are also contraindicated in case of total hip replacements, they are: Osteomyelitis in the upper femur, pyogenic infection of the hip joint, or evident infections are absolute contraindications. Every effort should be made to guarantee the asepis preoperatively in patients who have one or more of the following abnormalities: fever, local inflammation, rapid joint destruction, apparent bone absorption shown in roentgenogram; unexplained high sedimentation rate through other diseases, increase of the WBC rate, or suspected striking changes in the count. However, an infection recurrence has a higher probability to occur in revisions than in primary surgeres. Under no circumstances it is recommended that the insertion is done where there has been extensive infection involving the bone; or when involved soft tissue cannot be removed; or where the organism is gram negative, antibiotic resistant Staphylococcus aureus; or where highly virulent or mixed organisms are involved. Additionally, in distant site infection, such as genitourinary, lung, skin or in other organs, it is relatively contraindicated due to transmission of the infection to the hip joint site. The focus of infection must be treated before, during and after the implantation. The use of this implant is contraindicated when there is loss of musculature or neuromuscular disease in the affected limb, which would make the procedure unjustifiable.

**Warning:** improper implant selection, insertion, placement, or excessive muscular activity may cause unusual strain on the implant and may result in subsequent fracture of the prosthesis. Where there is acetabular failure or loss of acetabular bone stock it is advised that the use of grafts or other procedures to provide support to the acetabular component. Acceptable techniques in postoperative care must be followed. The patient should be advised about the limitations of partial joint reconstruction and his/her atual condition. Excessive physical activity and trauma that affect the partially replaced joint implicate in premature failure of the reconstruction caused by loosening and fracture in the use of implants. Excessive demands on the implant may occur because of excess of activities or for not controlling body weight. It is suggested to monitor the patient's recovery through X-ray. The acetabular components have a metallic ring that facilitates the visualization of the acetabular component on the X-ray. Do not use magnetic resonance imaging, because the ring manufactured in stainless steel causes interference in the examination impairing the diagnosis, and may also cause displacement of the implant causing trouble for the patient.

**Testing (Load Restrictions / Mechanical Strength):** It was observed roughness and roundness requirements aiming at eliminating the lack of congruency between Biomecnica Polymeric Cemented Monocomponent Acetabular Cup and Metal Femoral Head.

**Adverse Effects:** Early release of components of partial hip replacement may occur. Premature mechanical loosening may be the result of impaired fixation or hidden infection; late loosening resulting from biological complications or mechanical problems (localized high strain). Peripheral neuropathy and heterotopic bone formation were informed in Total Hip Arthroplasty. Subclinical nerve damage occurs more frequently, possibly as the result of surgical trauma. Reactions caused by metal sensitivity in patients followed by Total Articular Replacement were rarely reported. Implantation of foreign material in tissues may result in histological reactions involving various sizes of macrophages and fibroblasts. Similar changes may occur as a precursor, or during the healing process in the tissues. Dislocations and subluxation were informed, as the result of improper positioning of implant components. Muscle and fibrous tissue laxity can also contribute to these conditions.

**Information to be provided to the Patient:** Talk to your orthopedic surgeon about any medical condition that may affect the surgery. Total Hip Replacements are successful in 90% of patients. When complications arise, most can be successfully treated among the complications that may arise are: pain, infection; infection may occur in the wound or it can be deep (around the prosthesis). It may appear while patients are in hospital or at home. It may appear even years later. Small surgical wound infections are usually treated with antibiotics. Large or deep infections may require further surgery (or deep cleansing) or even the removal of the prosthesis. Any infection in your body (bladder, throat, teeth, ears, etc.) may take germs through the blood to your prosthesis and cause infection.

# ACETABULOS

**Thrombosis:** is the formation of a blood clot inside a blood vessel, obstructing the flow of blood through the circulatory system. Blood clots are the result of several factors, including its decreased mobility, and cause blood to flow more slowly in the veins, which may facilitate the development of thrombosis. A common sign of blood clots is pain and edema in the thigh or calf. If this occurs, your orthopedist will do tests and exams to evaluate the veins of your leg and may require monitoring and evaluation of a vascular surgeon. Several measures should be used to reduce the possibility of blood clots - Medication to thin the blood (anticoagulants) - Elastic Stockings - Exercise to improve circulation in the legs - plastic boots that inflate and compress the soles of the feet and calf, increasing venous return. Even with the use of these preventive measures, blood clots may still occur. If edema (swelling), redness or pain in the calf after leaving hospital is noticed, you should contact your orthopedic surgeon.

**Loosening:** Prosthesis loosening may occur after surgery. This can cause pain if loosening is significant. A revision surgery (implant exchange) may be required. New materials and new methods of attachment should minimize this problem.

**Luxation:** eventually, after a hip prosthesis, the head of the prosthesis may come out (dislocate) of the acetabulum. Most of the time, the dislocated femoral head may be pulled back in place without operation. To prevent dislocation, muscles need be strong (do exercises recommended by your orthopedic surgeon) and the leg must not be flexed above the hip (over 90 degrees) in the first month.

**Wear:** some type of wear may be detected in any type of prosthesis. Excessive wear may contribute to loosening and may require further surgery (Revision Surgery).

**Prosthesis Breakage:** Once implanted, the biomaterials used in prostheses necessarily come into contact with bodily fluids. These fluids, while seemingly harmless can, over time, significantly deteriorate the majority of materials of considerable "chemical inertness." Additionally, most orthopedic implants are subjected to static and /or cyclic mechanical strain, which may cause the failure of the material, that is, the breakage of the prosthesis.

With the materials currently used in implants, fracture of the implant is rare. However, if this occurs, further surgery to replace the fractured prosthesis will be required.

**Nerve damage:** nerves around the prosthesis can be damaged during surgery, although this is very rare. This injury is more likely to occur when the surgery involves correction of major hip deformities or it is necessary to lighten a leg that is too short due to wear. Over the time, these nerve injuries often improve and may fully recover. Eventually, a surgical exploration of the involved nerve may be required.

**Implanted Product Evaluations:** After implantation, in the intraoperative, the professional in charge must perform a radiological control to verify the correct positioning of the product. It is this professional's responsibility to make the clinical and radiological evaluation after the surgical procedure in the frequency stipulated by him/her to check the condition of the implant and the bone development. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecnica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

**Useful information to avoid risks derived from implantation:** To reduce the implantation risks the instructions contained in the Contraindications, Instructions for use and "Product instructions for Use" must be strictly followed.

**Implantation Intrinsic Risk:** Biomecnica Polymeric Cemented Monocomponent Acetabular Cups are made from materials of recognized biomedical use: NBR ISO5834-2 - Implants for surgery - Ultra-high-molecular-weight polyethylene - Part 2: Moulded forms, ASTM F138 - ASTM F138 - Standard Specification for Wrought 18 Chromium-14 Nickel-2.5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (Uns S31673), and NBR ISO 5833 - Implants for surgery - Acrylic resin cements. Contamination Risks: There are risks of biological contamination and the transmission of viral diseases such as HIV and hepatitis because the metallic components of Biomecnica Polymeric Cemented Monocomponent Acetabular Cup come in contact with body tissues and fluids. Explanted products should be treated as highly polluting.

**Sterilization:** This product is sterilized with Ethylene Oxide Gas ETO. Its sterilization expires in 5 years. If the validity of sterilization is overdue or packaging is damaged, the product should be returned to the distributor or directly to Biomecnica.

**Product discard:** No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way Biomecnica Polymeric Cemented Monocomponent Acetabular Cup explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecnica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

**TRACEABILITY:** On the packaging of the implantable component there are 6 labels containing the necessary data for implant traceability.

**Patient Card:** The patient card that accompanies Biomecnica Polymeric Cemented Monocomponent Acetabular Cup should be given to the patient, together with identification labels of products used in arthroplasty. With the card and labels in hand, the patient may obtain information regarding the quality of implant used. More information can be found at www.biomecnica.com.br / patient

**MARKING:** Biomecnica Polymeric Cemented Monocomponent Acetabular Cup contain the following information recorded with laser to trace and monitor the patient after surgery: Batch Number and Size

**Post-Sales (Customer Complaint):** If there is a need to make a complaint about Biomecnica Polymeric Cemented Monocomponent Acetabular Cup related to any adverse effect that affects the security of the patient, such as problems with the performance of the product, damage to the implantable metal component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the national health authority and Biomecnica by e-mail sac@biomecnica.com.br or by phone 55 14 2104-7906. In case of doubts, the surgeon in charge or health professional can make the communication of the adverse event through the Public Health Notification System at ANVISA website: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

**Sending the product for the manufacturer's Analysis:** when the implant is sent for the manufacturer's analysis, it should be cleaned at the hospital using a bactericidal and broad-spectrum antifungal solution. Then, it should be disinfected or sterilized in a steam autoclave or with ethylene oxide. It should be sent to Biomecnica in an intact packaging, identifying the method of cleaning, sterilization, and product data.

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