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

ENGLISH

The Family of Non-absorbable anchors – BM is a non-absorbable and non-active implant according to NBRISO 14630 manufactured from a polished and annealed round bar in titanium-6aluminium-4vanadium alloy (Ti6Al4V) in accordance with ASTM F136. Its chemical and mechanical compositions are within the limits established by the standards. The Family of Non-absorbable anchors – BM is composed of the models: Anchor screw in diameters of 2.0mm, 2.7mm, 3.5mm and 5.0mm. The Anchor screws present a pointed end that helps in the insertion of the screw into the bone wall. The threaded portion makes possible for the screw to be securely and easily inserted providing stability and tension to the fixation of the screw into the bone. A comparative study showed that the torsion strength of the threaded anchor screw is 382% higher than the torsion strength of the anchor screws with wings, what proves that the models with thread are superior to the models with wings when it comes to torsion strength.

Compositon:
The models of the screws that compose the Family of Non-absorbable anchors – BM are manufactured in titanium alloy in accordance with ASTM F136.

Presentation of the Medical Product:
The implants are packed in plastic envelopes with a non-sterile tag and outer labeling. They are distributed in single units. The product comes with the following information: Name of the manufacturer; Commercial name of the product; Technical name; Lot number; Registration number with the Ministry of Health; Code of the Product; Quantity; Product description; Manufacturing date; Expiration date; Not determined; Description of the raw material used in the manufacturing of the product; Address of the manufacturer (available in the plastic envelope); Name of the Responsible Technician; The sayings: "Non-sterile product", "Do not reprocess", "Single use Product", "Do not reuse", "Before use read the instructions".

Symbols of the Traceability tag (Label

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discomfort and/or abnormal sensations due to the Anchor screw and the Sutura wire and surface in deep infection.

Warnings:

The incomplete insertion of the Anchor screw will not allow an effective fixation; the degree of the postsurgical mobility will be determined by the doctor and shall be informed to the patient.

During the implantation of a device there is a risk of inserting a foreign and particulate material including glue talc, lint and cleaning material besides other contaminants of the surface. All efforts shall be made to limit the handling of the fixation devices; the doctor is responsible for choosing and implanting the device as well as for deciding on leaving or removing the implants after they fulfill their clinical functions; The implants shall never be reused and the explanted devices shall never be implanted again, because they may present little defects and inner stress that may cause their failure and breakage even if the implants look intact. If the patient is involved in some activity or job that cause stress to the implant (substantial walk, race, weight lifting or muscle stress) there is a chance of failure occurrence to the device.

The patient shall understand that the Anchor screw is not as strong as a healthy bone and, therefore, it may break in case it receives over stress. A patient that is not able to assimilate that information runs serious risks during the recovery process. It is suggested using two anchors of smaller diameter than one of greater diameter due to the anchors of smaller diameters cause less damage to the receptive bone and better distribute the loads.

Information to be provided to the patient:

It is available in the instructions in terms: Indications, contraindications, information for use, Potential adverse effects, Precautions and warnings, the patient shall be informed of the importance of a postsurgical follow-up. The lack of a follow-up prevents the detection of postsurgical problems. The failure to submit the patient to a revision surgery in case of the component loosness or osteolysis may result in a progressive loss of the periprosthetic bone stock. The patient shall be informed of the type of the implant, brand, code, name, lot, name of the surgeon and date of the surgery. The patient shall fully understand that the product does not replace or have the same performance of a healthy bone and, therefore, it may break, deform or detach due to over stress or activity, early loading, etc. The patient shall be advised of the need of performing regularly postoperative follow-up while the device remains implanted.

Important:

The surgeon shall guide the patient to the cares to be taken and the necessary physical support necessary for the time the implant remains implanted. Before a perfectly stabilized fracture, the load may be liberated. If the fracture is not stabilized, no load can be applied to the implantable system.

Evaluation of the implanted product:

After implantation in the intraoperative, the responsible Professional shall do the radiological control to check for the correct position of the product. The responsible Professional shall do all the clinical and radiological evaluations after the surgical procedure at a frequency he stipulates in order to check the conditions of the implant and evaluate the evolution of the healing. In case the product is out of its correct position or presents any non-conformity, it is the surgeon's responsibility to take the most appropriate corrective actions.

Useful information to avoid risks resulted from implantation:

In order to minimize the risks resulted from implantation, the following shall be strictly obeyed: contraindications, instructions for use and all other information available in the "Instructions for Use" of the product.

Intrinsic Risks of Implantation:

The Family of Non-absorbable anchors – BM are manufactured with well known biomedical use materials which are in accordance with the requirements of the Standards ASTM F136-Standard Specification for Wrought Titanium-6 Aluminum-4 Vanadium ELI (Extra-Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401).

Risks of Contamination:

There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis, because the metallic components of the Family of Non-absorbable anchors – BM get in contact with tissues and body fluids. The explanted products shall be treated as highly contaminating.

Decontamination, cleaning and Sterilization:

The Family of Non-absorbable anchors – BM is distributed non-sterile; therefore, the manufacturer is responsible for the sterilization and shall use highly regulated methods for it. The product is provided in its original packaging and previously cleaned. Before use, the packaging shall be inspected for breakage and damages. The product shall only be unwrapped before sterilization in order to preserve the surface finishing and configuration intact. It shall be handled the least as possible as recommended in the standard NBR ISO 8528. In order to guarantee the reduction of risks of infection, the implants (the ones that were not used but were set on the surgical table or lost sterilization) shall be cleaned before sterilization. In case of the surgical instruments, they shall be cleaned after surgery and before sterilization. The cleaning may be performed through mechanical or manual means. Cleaning equipments: -Ultrasonic cleaner; -Washer –sterilizer; Washer –desinfector. Manual cleaning: it is the procedure that removes the dirtiness through physical action with the help of detergent, water and artifacts such as sponge or brush. The manual cleaning is mostly recommended due to it is less aggressive to the implants and surgical instruments. For the manual cleaning, it shall be used water at room temperature, nylon brushes and never steel or abrasive ones, because they may impair the protective layer of the material. For the decontamination, it shall be used a bactericide and antifungal solution of broad spectrum and avoid aggressive cleaning agents and metal brushes in order not to impair the quality of the metallic implants. In case of cleaning machines, the implants and surgical instruments shall be set in drawers in a way they do not hit each other to avoid damage to the material. People responsible for the cleaning of the implants and surgical instruments shall be aware of the type of the cleaning product respecting the concentrations indicated by the manufacturer. The implants and surgical instruments shall be carefully rinsed up to all residue is out and immediately dried after cleaning.

Type of cleaners

Enzymatic cleaners: basically, they are composed of enzymes, surfactants and solubilizers. A balanced combination of those elements makes it possible for the product to remove all the organic matter from the material in a short term.

Enzymatic solutions: they present excellent cleaning action but do not present any bactericide or bacteriostatic activity.

Enzymes: they are substances produced by living cells that govern the chemical reactions of the process. Once produced by the cells an enzyme may be isolated and will still keep its catalytic properties if certain conditions are kept in its manufacturing. The enzymes are classified in three major functional groups depending on the type of substrate they will affect: protein, lipases and amylases that act on the protein substrates, fats and carbohydrates which tend to solubilize and detach from the articles. Currently, it is recommended the cleaning of complex configuration articles for the guarantee of the cleaning.

Desinfection:

It is a process that destroys pathogenic or non-pathogenic microorganisms from the articles except for the bacterial spores through physical and chemical means.

Level of desinfection:

High Level: it destroys the microorganisms except for the ones with a high number of spores => 2% Glutaraldehyde –20–30 minutes.

Indication: specially hospital area

Medium level: it eliminates vegetative bacteria, most of the virus, fungi and microbacteria => 1% Sodium hypochlorite –30 minutes.

Indication: for UBS, Day cares, asylum and nursing home.

Low level: it eliminates most of the bacteria, some virus and fungi but not the microbacteria => 0,025% Sodium hypochlorite.

Sterilization:

This product is provided non-sterile. It must be sterilized before use. We recommend steam sterilization in Autoclave at hospital ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1. Requirements for the development, validation and routine control of a sterilization process for

medical devices)

Autoclave:

It is an equipment for sterilization through saturated steam under pressure.

Operating instructions:

Before starting the equipment, check:

if the water valve is open;

if the flush valve is closed;

1 – open the door of the equipment;

2 – properly set the material to be sterilized;

3 – close the door of the equipment;

4 – select the desired cycle according to the material to be sterilized;

5 – turn on the general switch; 6 – the cycle will automatically run in the sequency;

7 – when the lamp "end of cycle" turns on, partially open the door for tem (10) minutes for the material to cool down.

It is recommended applying the following parameters of physical sterilization in autoclave (saturated steam):

Cycle	Temperature	Exposure Time
Conventional (1 atm of pressure)	121°C (250°F)	30 minutes
Conventional (1 atm of pressure)	132°C (270°F)	15 minutes
Severity	132°C (270°F)	45 minutes
High Vacuum	132°C (270°F)	7 minutes

Note.: The time when the heat of the sterilization chamber reaches the desired temperature shall be recorded.

FOR FURTHER CLARIFICATION CHECK THE INSTRUCTIONS MANUAL THAT COMES WITH EACH AUTOCLAVE.

Other methods of sterilization:

Besides the autoclave the following method can also be used as defined by the hospital:

- Sterilization in Ethylene oxide (E.T.O) – parameters and procedures established in the validation protocol and ISO 11135-1. Sterilization of health care products – Ethylene oxide -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices.

NOTE: It is under responsibility of the hospital, the method of sterilization, the equipments, the controls and the instructions for the sterilization used. Cares with the sterilized articles.

Storage conditions of the sterilized articles:

- As the environment: it shall be cleaned, ventilated, dried and restricted to the staff of the sector.

- As the article: after the sterilization process do not place it on a cold surface (stone or stainless steel), use perforated baskets or recipients up to they cool down;

- Wrapper (raw cotton fabric, nonwoven fabric, surgical paper, crepe paper, paper with a film, tyvec or perforated metallic boxes) shall remain intact and be little handled to avoid that the package tears or loses its seal;

- Store it in closed cabinets with shelves;

- Identified shelves to facilitate the removal of the material;

- The hospital shall register that information in its files and forward it to the patient. We inform that it is mandatory that the patient's chart brings all information to allow the traceability of the implanted product. The tag pasted on the chart allows the traceability of the used product. Among the main information we consider indispensable the name of the manufacturer, the name of the implant, the code of the implant, lot number and product registration with ANVISA. That information is listed in the traceability tags that come with the product and in the outer labeling. The following information also be considered: date of the surgery, patient's name, surgeon's name, weight and age and other information required in the patient's chart. Other information required in the patient's chart shall also be filled out.

Post-sale (Patient's complaint): in case it is necessary to make a complaint of the Family of Non-absorbable anchors – BM related to an adverse effect that affects the safety of the user such as the product is not working, damage to the implantable component, serious problems or death related to those components, the responsible surgeon shall inform the competent health body and Biomecnica of the adverse event through the e-mail sac@biomecnica.com.br or by the phone 0xx14 2104 7900. In case of questions the responsible surgeon or health professional shall inform the Health Surveillance Notification System of the adverse event on the website of ANVISA: http://www.anvisa.gov.br/hotite/notivisa/index.htm

Traceability:

Together with the implantable components come two traceability tags with the data of the implant. Those tags bring the same information as the ones in the outer labeling of the product. The information available in that traceability tag is listed in the item "Presentation" of the technical report. The traceability tag shall be pasted on the patient's chart. The chart for the hospital control, registration and to keep the information in its files and forward it to the patient. We inform that it is mandatory that the patient's chart brings all information to allow the traceability of the implanted product. The tag pasted on the chart allows the traceability of the used product. Among the main information we consider indispensable the name of the manufacturer, the name of the implant, the code of the implant, lot number and product registration with ANVISA. That information is listed in the traceability tags that come with the product and in the outer labeling. The following information also be considered: date of the surgery, patient's name, surgeon's name, weight and age and other information required in the patient's chart. Other information required in the patient's chart shall also be filled out.

Post-sale (Patient's complaint): in case it is necessary to make a complaint of the Family of Non-absorbable anchors – BM related to an adverse effect that affects the safety of the user such as the product is not working, damage to the implantable component, serious problems or death related to those components, the responsible surgeon shall inform the competent health body and Biomecnica of the adverse event through the e-mail sac@biomecnica.com.br or by the phone 0xx14 2104 7900. In case of questions the responsible surgeon or health professional shall inform the Health Surveillance Notification System of the adverse event on the website of ANVISA: http://www.anvisa.gov.br/hotite/notivisa/index.htm

Send the Material for the manufacturer's analysis: in case the implants are sent for the manufacturer's analysis, they shall be cleaned at hospital with the help of a bactericide and antifungal solution of broad spectrum. Then, it shall be disinfected with steam sterilized autoclave or ethylene oxide. It shall be sent to Biomecnica in intact packages identified with the cleaning methods, sterilization and data of the product.

Description of the effectiveness and safety of the medical product in accordance with the regulations of ANVISA that rules on the Essential requisites for Effectivity and Safety of the Medical product
The risks associated with the Family of Non-absorbable anchors – BM are in accordance with the Essential Requisites for Safety and Effectiveness of the Product in accordance with the RDC 56/01 for products classified in rule 8, risk class 3. The risk managements of the product and process are in accordance with the NBR ISO 14971 as listed in the PRCE02 procedure – Risk analysis of the Medical device. The Family of non-absorbable anchors – BM is manufactured in Titanium alloy in accordance with ASTM F136 as verified through laboratory testing reports and quality certificate of the product. Those materials are acceptable for the manufacturing of implants in accordance with annex A of the standard NBRISO21534.

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