FIXADORES EXTERNOS

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ENGLISH

Product Description

The external fixators have two basic functions: a) Fuse the fractured parts of diaphyseal or epiphyseal bone tissue. That will promote dynamic compression or distraction or static stabilization between the fragments. b) Promote the elongation of the bone tissue in case of shortening of it. The external fixator is a medical product which is not implantable and it is designed to help the bone fixation and intended to: minimizing, stabilizing and fixing several types of fractures in the lower and upper limbs as well as in the human skeleton according to the indication for use. It also has the purpose of transferring the load applied to the fractured bone to the device which keeps the bone free of load and able to regenerate after the fracture. The external fixators present a structure designed exclusively for tubes and connectors which are connected by wires and pins. The external fixators are available in different designs and with proper methods of application besides general principles that direct their applications. The wires and pins are ancillary components for the external fixators, therefore they are validated in another process of registration and sold separately.

Manufacturing Materials

The external fixators are made of materials that are well known in the medical field and that meet ASTM Standards for the steels and NBR Standards for the aluminum. Two types of raw-materials are used – ASTM F899 for 304 typed stainless steel; ASTM F899 for 420 typed stainless steel and NBR 6834 for 6065-T5 aluminum alloy. The table below brings the codes, the components of the set and the raw-material of each component of the fixator.

List of external fixator

FBM 800 Fixator - Fêmur / Tibia (50051-00-000)

Components of the Set: 02 Aluminium Bar - 07,96 x 340mm; 02 Stainless Steel Bar - Ø4,8 x 220mm; 12 Platform - 8mm; 01 Guide - 1 hole; 12 Steel Connector; 12 M6 Nut; 12 M6 Washer

Femoral Proximal Fixator (50040-00-000)

Components of the Set: 01 130° Platform; 01 80° Platform; 02 M8 x 150mm Bar; 01 Stabilizer; 01 Disposable Tensor (Guide + Screw); 01 Aluminium Guide - 2 holes; 04 M8 Washer

Iliac - Anterior Fixation (50031-00-000)

Components: 01 M6 x 270mm Threaded Bar; 06 Connector with no bar; 04 M6 x 21mm Connector with Bar; 02 Pin holder platform; 02 Lid for Humerus Platform; 02 Platform for Pubis; 27 M6 Nut; 05 M6 x 30mm Screw; 02 M6 x 100mm Threaded Bar; 02 M6 x 180mm Threaded Bar

JS FX - Femoral Diaphyseal (50024-00-000)

Components: 02 Platform - 90 mm; 02 M8 x 340 Threaded Bar; 04 M6Connector - normal P; 04 M6 Connector - Long G

JS FX - Tibial Diaphyseal (50023-00-000)

Components: 02 Platform - 60 mm; 02 M8 x 300 Threaded Bar; 04 M6 Connector - normal P; 04 M6 Connector - Long G; M8 Nut - narrow; M6 Nuts - Calote

Diaphyseal Fixator - Radio / Ulna (50018-00-000)

Components: 02 M6 x 200 mm Bar; 02 Pin holder Platform; 04 M6 x 25mm Connector; 04 M6 x 21 mm Connector; 04 M6 x 25 Threaded Bar; 04 M6 x 21 Threaded Bar; 16 M6 Nut

Articulating Fixator for Fist (Colles type) (50019-00-000)

Components: 02 Platform for Colles Fixator; 02 Head for Colles Fixator; 01 M8 x 53mm Threaded Bar; 01 Colles Ring; 01 M8 Nut; 01 Articulating Nut; 01 Block Ring; 01 Axis for Colles Fixator; 06 M6 x 8mm Headless Allen Screw

FBM 820 Fixator - Rádio / Ulna (50053-00-000)

Components: 12 Platform - 6mm; 12 Steel Connector - 6mm; 12 M6 Washers; 12 M6 Nut; 02 Stainless Steel Bar - Ø3,0 x 120mm; 02 Aluminium Bar - Ø6,3 x 120mm; 01 Aluminium Guide - Hole - Ø3,0mm

FBM 810 Fixator - Humerus (50052-00-000)

Components: 02 Aluminium Bar - 07,96 x 230mm; 02 Bar - Ø4,0 x 220mm; 12 8mm Platform; 01 Guide - 1 hole; 12 M6 Nut; 12 M6 Washer; 12 Steel Connector

Metacarpal/Proximal/Distal Phalange / Fixator (50017-00-000)

Components: 02 Platform - "T" Base; 02 Platform - "T" Lid; 01 M4 x 60mm Threaded Bar; 02 Allen M4 x 14mm Screw; 04 M4 Nut

Symbology in the traceability tag (Label) according to NBRISO15223 and En980

XXXXX	Lot nº	XXXX-XX	Code
XXXX-XX	Validity		Sterilized in ETO
XXXX-XX	Manufacturing date		European representative
	Manufacturing data		

TECHNICAL INFORMATION

Indications, Precautions, Restrictions/Warnings, Adverse effects and Contraindications: The external fixator is indicated for bones of lower and upper limbs in corrective surgeries or in the treatment of traumas because the fixators keep the reduction of the fractured bone and stabilizes the bone or adjacent structures in order to promote or correct the bone fusion. The surgeon must be qualified and familiar with osteosynthesis and its limitations including the pre and postoperative, surgical techniques to be adopted, precautions and potential risks. Do not use the product if it is damaged. The pre-operative and surgical procedures as well as

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Sistema de Qualidade Certificado
Quality System Certified
ISO 13485

FIJADORES EXTERNOS

the knowledge about the surgical technique and proper selection and implantation of the implant are important considerations for a successful use of the provisional devices for external fixation. Do not use the product in case the packaging is damaged. Some potential adverse effects that may occur in the use of the external fixators are: postoperative superficial or deep infection; pain; discomfort or abnormal sensations due to the presence or migration of the medical device or due to the surgical procedure and migration; limb edema; pseudarthrosis; movement limitations; failure of fracture healing. Due to the incompatibility of the materials, design and dimension the use of components from different manufacturer is forbidden. The use of different manufacturers may cause problems to the product and bring risks. The company disclaims all responsibility of any adverse event that may occur. Do not use implants from different manufacturers. Therefore it is recommended that the products have the same provenance.

Sterilization: this product is provided sterile. The sterilization is in Ethylene oxide (E.T.O). In case the sterilization date is expired or the packaging is damaged the product must be returned to the respective distributor or directly to Biomecánica.

Storage and transport of the product: The storage location should be clean, dry and illuminated in order to maintain the proper storage and transport conditions as well as the physical and chemical integrity. The product must be stored and transported in dry and fresh environment at room temperature (Max. 35°C) and relative humidity around 30% to 70% and away from direct sunlight. The effects of vibration, shocks, corrosion, defective settlement during the handling and transport and improper pilling during storage must be avoided. The product must be carefully handled in order to avoid damages that may impair the quality of the material and the safety of the patient. The products can not be stored in direct contact with the floor (minimum height = 20cm). They can not remain in very high places close to lamps because it may respect the packaging or damage the label. They should not be stored with contaminant substances such as cleaning materials, insecticide, pesticide, etc. The transporter must be advised of the content and delivery deadline. Carefully transport the products to avoid fall and attrition in order not to damage the packaging and device; that guarantees the sterility of the material. The blister sealed packaging may be manually open in a proper place (surgical room) by pulling the surgical paper from the peeling around the blister. The packaging should only be open at the moment of the surgery in order to avoid the contamination of the product.

Identification: The external fixators come with labels that present the following information: Data of the manufacturer; Authorized representative in the European community; Technician in charge; n° of the register with ANVISA; Technical name of the product; Commercial name of the product; Commercial model of the product; Raw-material; Quantity of items per package; Lot number; Code of the product; Manufacturing data; Expiration date; Address of the manufacturer; The external fixators come with the following engraved information for traceability: Biomecánica logo; Number of the manufacturing lot and code.

Customer Service: In case it is necessary to make a complaint about the use of the external fixators related to some adverse effect that may affect the safety of the user the responsible surgeon must inform the Health Authorities and Biomecánica through the email - **e-mail sac@biomecánica.com.br** - or by the phone - **0xx14 2104 7900**. In case of doubts the responsible surgeon or the health Professional may inform about the adverse event through **NOTIVISA on the webpage of ANVISA:** <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

Product discard: No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the external fixator explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecánica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Do not re-use the product. An orthopedic implant may be implanted in one single patient and for one single time. Although it may not look damaged, previous stress may have created imperfections which may impair the success of the implant. Due to its presentation and characteristics this product does not present any side effects. However it is suggested that the patient's reestablishment is monitored through x-ray in accordance with the International Safety Standards.

Do not submit the patient to magnetic resonance because the products in stainless steel cause interference in the exam and that impairs the diagnosis and causes the displacement of the implant besides harming the patient.

Contamination Risks: There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the implant gets in contact with body tissues and fluids. The removed implants should be treated as highly contaminants.

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MDSS GmbH