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ENGLISH

SPECIAL PLATES FOR BONE SYNTHESIS

Special Plates for Bone Synthesis are produced in stainless steel, according to the ASTM F139 standard specifications. The stainless steel used in the factoring of the Special Plates for Bone Synthesis is analyzed in a third part laboratory to the material conformity check.

T* 4,5 mm Plate

Plates in T* 4,5 shape are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being used to support metaphyseal fractures of the proximal humerus, of the tibia medial plateau and distal of shin-bone. Its shape has a double fold allowing an adaptation on the indicated bone region. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting up and down.

Semi-Tubular plate

Semi tubular plates can be easily adapted. Its rotational stability is due to the edges that "bite" the sub adjacent surface when the screws are held tight. These plates have been used in round bones such as the radio or in crests and just as in tension band plate, as its maximum resistance is in tension. The main indications to these plates are radio fractures and proximal ulna, particularly comminuted fractures of the olecranon, and fractures of the top boots of distal tibia.

****T* 3,5 mm Plate**

Plates in T* 3,5 shape are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being indicated for setting in the distal radio region. Its shape has a double fold allowing an adaptation of the distal radio. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting the sagittal plan.

Left or Right L' Plates

L' shaped plates are relatively thick plates, destined to support a thin cortex or to prevent a collapse of a spongy bone defect, being indicated for setting in the tibia proximal region. Its shape that has a double fold allows the adaptation to the tibia lateral plateau. It has an oblong hole that allows a temporary fixation of the plate to the bone, also allowing some adjusting the sagittal plan.

Composition: Special Plates for Bone Synthesis are manufactured in stainless steel, according to the ASTM F139 standard specifications.

Presentation Form: Implants are conditioned in plastic envelopes with a non-sterile tag and external label. The external label is made of a product tracking tag, in which contains the following information:

a) Producer's Name; b) Commercial Name of the Product; c) Technical Name; d) Lot Number; e) ANVISA registration number; f) Product Code; g) Quantity; h) Description of the product that is included in the package (model); i) Production Date; j) Expiration Date; Undetermined Expiration Date; k) Description of the raw material used in the production of the product (Composition); l) Producer's Address; m) Name of the Responsible Technician; n) The Saying: "Non-sterile Product"; o) The Saying: "PROHIBID TO REPROCESS"; p) "Single Use Product, do not

SPECIAL PLATE

reutilize"; q) Before using, read the usage instructions; r) Warnings/Precautions/Special Cares/Symbology: Read the usage instructions.

	Manufacturer Date		Do not use if package is damaged		Expiration date
	"Read Instructions"		Medical Product for unigue use. Do not re-use.		Keep dry
	Keep away from sunlight		Fragile. Handle it carrefully		

Description of the packaging material

The non-sterile product is identified with a green tag of 20mm diameter with the phrase "Non-Sterile"

Attention: Relation of instruments for use in Special Plates for Bone Synthesis. For use of the Special Plates for Bone Synthesis it is necessary the use of the specific instruments mentioned below. The charts below describe the relation of the instruments that compose the trays in the cases 4,5 and 3,5.

Relation of the instruments that compose the tray of the instrument for placing the Plates and Screws to the Bone Synthesis 4,5
Tray of instruments for placing the Plates and Screws of the Bone Synthesis 4,5
Hohmann Retractor Large (6019-0G-000), Drilling Guide (6011-00-000), Drill - Stainless Steel Ø3.2 X 150mm (6003-32-150), Drill - Stainless Steel Ø4.5 X 150mm (6003-45-150), Hexagonal Wrench (6015-00-000), Compressor Wrench (6007-00-000), Clamp Compressor (6005-00-000), Countersink (6016-00-000), Handle Double Guide Of Perforation N/ox/std (6013-00-000), Drill Guide 4,5mm - STD (6013-01-000), Drill Guide Ø4,5 - N (6013-02-000), Drill Guide Ø4,5 - EX (6013-03-000), Tap Protector Guide 4,5 (6014-00-000), Tap For Cortical Screw Ø 4,5mm (6002-00-000), Tap For Cancellous Screw 6,5mm (6001-00-000), Cortical Measure Ø4,5 (6004-00-000), Self Centering Tweezer 4.5 (6043-00-000), Reduction Tweezer 4.5 (6042-00-000),bending Iron Plate 4.5 (Pair) (6017-45-000)

Relation of the instruments that compose the tray of the instrument for placing the Plates and Screws to the Bone Synthesis 3,5

Tray of instruments for placing the Plates and Screws of the Bone Synthesis 3,5

Hohmann Retractor Medium (pair) 3,5 (6019-0M-000), Hohmann Retractor Small (pair) 3,5 (6019-0P-000), Drill - Stainless Steel Ø2.0 X 150mm (6003-20-150), Hexagonal Wrench 3,5mm (6036-00-000), Small Countersink 3,5mm (6033-00-000), Guide Of Perforation 4,5 N/EX/STD (6013-00-000), Drill Guide STD/NE/EX Ø3,5 (6003), Drill Guide 2,5 (6039-00-000), Drill Guide 3,5mm (6029-00-000), Tap For Cortical Screw Ø3,5mm (6021-00-000), Cortical Measure 3,5 (6024-00-000), Bending Iron Plate (pair) (6017-00-000), Tweezer For Screw (6006-00-000)

Care For Handling and Transportation of the Medical Product:

-Do not use the product if it is damaged.
-The implants must be handled with care avoiding damage that can spoil the quality of the material and the safety of the patient.
-Preserve it in a cool, dry place, away from light and away from the action of weather.
- While transporting it, shocks and inadequate stacks must be avoided.

FOR BONE SYNTHESIS

standard Conventional cortical screw Ø 4,5mm, total spiral manufactured in austenitic stainless steel according to ASTM F138 standard Spongy screw Ø 6,5mm, total or partial spiral, manufactured in austenitic stainless steel according to ASTM F138 standard

Accessories: Special Plates for Bone Synthesis do not have any accessory with the purpose of integrating the medical product.

Allowable Combinations with other Materials: The chart "Allowable Alloy of Contact" was based on the "NBRI5021534 – Implants for Non-active Surgery – Implants for Articulation Substitution – Particular Requirements" standard, but independent of the alloy mentioned on the "Of Contact" item, it is necessary that the combinations are from implants of the same brand, because they were projected for such combinations observing the superficial finish and superficial treatment and other demanded factors in projects that can interfere in such combinations. Because of that, metallic implants from different manufacturers are not recommended for chemical, physical, biological and functional incompatibility reasons.

Chart: Allowable Alloy of contact	
Implant	*Of Contact (when there will be no articulation)
Plate (without ring) (Stainless Steel ASTM F138 - 18 Cromo-14 Nickel - 2,5 Molybdenum)	Stainless steel adapted 18 Cromo-14 Nickel -2,5 Molybdenum (NBRI50 5832-1 and ASTM F138) <p>Stainless steel adapted of high nitrogen (NBRI505832-9) Adapted league of titanium 6-aluminio 4-vanadium (NBRI50 5832-3 and ASTM F136) Adapted league of titanium 5-aluminio 1,2,5-iron (NBRI505832-10) Adapted league of titanium 6-aluminio 7-niobium (NBRI505832-11)</p>

*For application where a metal or alloy is in contact with others and is not intended joint, provided that adequate attention is paid to the design, surface finish, surface treatment and metallurgical conditions.

The Special Plates for Bone Synthesis from Biomecânica brand must be used only with Biomecânica instruments.

Precautions, restrictions e warnings

PRECAUTIONS: SINGLE USE. DESTROY IT AFTER EXPLANTED. DO NOT REUSE THE PRODUCT. REPROCESS IS FORBIDDEN. An orthopedic implant can only be used in one patient, one single time. Although it may appear undamaged, previous stress may create imperfections that may reduce the implantation success. Improper implant selection may cause uncommon stress on the implant and can result in subsequent fractures of the implant.

Due to its shape of presentation and to its characteristics, this product does not show side effects, however, it is suggested to monitor the patient's recovery through x-ray according to the International Safety Rules. Do not use magnetic resonance imaging, as the product is made of stainless steel causes interference in the examination harming the diagnosis.

The use of Special Plates for Bone Synthesis is conceived to be an assist element to fixation and bone stabilization on the superior and inferior members fracture treatment and not to substitute normal human skeleton structures.

The surgeon must be sure of the possibility of the patient's physical development, being able to have a necessity of reviewing the implant.

The surgeon must be familiarized and have sufficient knowledge of osteosynthesis and its limitations, including the pre and post-operative, adopted surgery technique, precautions and potential risks. The surgeon must have attention to the bone preparation for the introduction of the special plates, because the stability and fixation success may be compromised.

During handling of the implant, there must be no scratches, nicks, or anything that may damage or mark the implant, because these defects are stress concentrators and may be sites of crack nucleation and decrease the corrosion resistance and may result in fractures on the implant. In certain situations it is necessary to shape the implant to fit it perfectly to the bone anatomy to be implanted. Should it prove necessary to modeling the implant to the bone anatomy is the limit of modeling. To guarantee an adequate implantation, should be used in instruments manufactured by Biomecânica, because they are designed and manufactured specifically for the use of these implants. Using instruments from different manufacturers may compromise the surgery as well as not being in harmony with the registration of the product in Agencia Nacional de Vigilância Sanitária; do not use the product if it is damaged;

NON-STERILE PRODUCT – STERILIZE BEFORE USE
SINGLE USE PRODUCT – IT IS HIGHLY FORBIDDEN TO REUSE THIS PRODUCT;

Special storage conditions: Store in a cool, dry place, away from light and away from the action of the weather.

Do not use the product in case the package is damaged;

Observation: Implants' components from different producers shall not be used; so we recommend that products have the same origin.
Manufacture Date, expiration date and product lot: SEE LABEL.

Contraindication: The use of this implant is contraindicated under the following situations:

-Sensitivity to metal;
-Fever or sign of local inflammation;
-osteoporosis;
-Obesity;
-Inadequate tissue cover in the region to be operated;
-Pregnancy;
-Bone compromising from illness and/or infections;
-Elevated sedimentation rate unexplained by other diseases, raising the rate of WBC or suspected striking in the differential count
-Other infection.

Adverse Effects: mechanical slackening can be a result of defective fixation, unstable reconstruction of the fractured bone or hidden infection. Reactions to metal sensitivity in patients were rarely informed. Implantation of strange material in tissues results in histological reactions involving macrophages and fibroblasts of several sizes. The chemical importance of this effect is uncertain, but similar changes can occur as a precursor, or during the healing process.

Information to be given the patient: The patient must be informed about the importance of the post-surgical attendance. The lack of attendance obstructs the detection of post-surgical problems such as components looseness or osteolysis occurrence. Not performing the review surgery when there is component looseness or osteolysis can result in progressive loss of periprosthetic bone stock. The patient must be informed as to the type of implant being used, brand, code, name, lot, surgeon's name and date of surgery.

Evaluation of the implanted Product: after the implantation, intra-operative professional must perform a radiological control to verify the correct positioning of the product. The responsible professional must perform, and it is his responsibility, clinical and radiological evaluations done after the surgical procedure at a frequency stipulated by him to check the state of the implant and its bone consolidation evolution. In case the product is found out of the correct positioning, or presents any non-agreement; it is the surgeon's responsibility to take the most appropriate corrective measure.

Useful Information to Avoid Recurrent Implantation Risks: to decrease the risks due to implantation this should be strictly followed: contraindications, instructions for use and all the information contained in the product's "Use Instructions".

Implantation Intrinsic Risks: The Special Plates for Bone Synthesis of materials are recognized biomedical use, which follow the requirements of the ASTM F139 Standard, Standard Specification for Wrought 18Chromium-14Nickel-2.5 Molybdenum Stainless Steel Sheet and Strip for Surgical Implants.

Contamination Risks: there are risks of biological contamination and viral diseases conveyance like hepatitis and HIV, because the metallic components of the Special Plates for Bone Synthesis get in contact with body tissues and fluids. Explanted products must be treated as highly contaminating.

Sterilization: this product supplied non-sterile. It must be sterilized before use. We recommend the steam sterilization in Autoclave at a hospital (ISO 11134: 1994 Sterilization of health care products – Requirements for validation and routine control – Industrial moist heat sterilization).

Product discard: No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the Special Plates for Bone Synthesis explanted

from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecânica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Tracking: together with the implantable component's package are 2 tracking tags containing the data of the implant used. This tag is the same tag contained in the external label of the product. The information that is in this tracking tag is related in the "Presentation Form" item of this technical report. This tag must be attached to the patient's handbook and the other for hospital control. The hospital must register this information in its files and give this information to the patient. We inform that the patient's handbook must have this information register which allows the implanted product to be tracked. Attaching this tag at the handbook allows that tracking of the implanted product used. Among the main information we emphasize are the name of the manufacturer of the implant used, name of the implant used, code of the implant, lot number and register number of the product at ANVISA. This information is written on the tracking tags that goes with the product and also on the external label. Other information also must be considered essential, as the date of the surgery, name of the patient who received the implant, name of the surgeon, patient's weight, patient's age and other information requested on the patient's handbook are also to be filled.

Marking: the Special Plates for Bone Synthesis contains the following information laser recorded to allow the tracking and patient's attendance at the post-operative:
-Biomecânica's Logo;
-Manufacturing Lot Number;
-Material Acronym (Stainless Steel) used to manufacture the product (S1);
-Code;
-Dimension;
-Ce Brand

Post Sell (Client's Complain): in case there is a necessity of doing any complaints about the Special Plates for Bone Synthesis or any adverse effect that affects the user's safety, such as non working product, damage of the implantable metallic component, serious problem or death related to these components the responsible surgeon must communicate this adverse event to the competent sanitary agency and to Biomecânica through the e-mail sac@biomecânica.com.br or through the phone 00x14 2104 7900. In case of doubt the responsible surgeon or health professional can communicate the adverse event through the Sanitary Vigilance Notification System at the ANVISA web-site: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

Send of Material for the Manufacturer to Analyze: In case of sending the implant for the manufacturer to analyze, this must be cleaned at the hospital using a bactericide solution and large specter antifungal. Then it must be disinfected or steam sterilized in autoclave or in ethylene oxide. Must be send to Biomecânica inside full package, identified with the clean method, sterilization and data of the product.

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