# PLACAS RETAS

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# STRAIGHT PLATE

Conventional Cortical Screw, non-self-tapping, fully threaded made of austenitic stainless steel according to ASTM F138 // Ø6.5mm Conventional Cancellous Screw, fully threaded, made of austenitic stainless steel according to ASTM F138

## 3.5mm Straight Plates – Small Fragments

Ø3.5mm Conventional Cortical Screw, non-self-tapping, fully threaded, made of austenitic stainless steel according to ASTM F138 // Ø3.5mm Conventional Cortical Screw, self-tapping, fully threaded made of austenitic stainless steel according to ASTM F138

Accessories: The Straight Plates for Osteosynthesis do not have any accessories with the purpose of integrating the medical product.

### Allowable combinations with other materials

The table below "Allowable Contact Alloys" was based on the standard "NBRIISO21534 – Non-active Surgical Implants – Joint Replacement Implants - Specific requirements," but despite the mentioned alloys in the item "In Contact" the combinations must be made with implants of the same brand because they were designed for such combinations observing surface finishing and surface treatment and other factors required in their projects that can interfere with the combinations. Therefore, metal implants from different manufacturers are not recommended because of chemical, physical, biological and functional incompatibility. The table, Allowable Contact Alloys\* describes these alloys.

### Allowable contact alloys

Implant:	* In contact (when there is no articulation)
Straight Plates for Osteosynthesis	Stainless Steel  Wrought Stainless Steel (ASTM F138, NBRIISO5832-1)  Wrought High nitrogen stainless steel (NBRIISO5832-9)  Wrought Titanium 6-aluminum 4-anadium Alloy (NBRIISO 5832-3 and ASTM F136)  Wrought Titanium 6-aluminum 7-niobium alloy (NBRIISO5832-11)

\* For applications where a metal or alloy is in contact with others and articulation is not intended, since adequate attention is given to the design, surface finishing, surface treatment and metallurgical conditions.

### Precautions, restrictions and warnings

Precautions: SINGLE-USE PRODUCT. DESTROY AFTER REMOVAL. DO NOT REUSE THE PRODUCT. PROHIBITED TO REPROCESS IT.

An orthopedic implant may be used in a single patient only once. Although it may appear undamaged, previous stress may create imperfections that may reduce the success of the implant. Improper selection of implant can cause unusual stress on it and may result in subsequent fracture of the implant.

Due to its form and characteristics, this product does not have side effect, however, it is suggested to monitor the patient's recovery through X-ray according to the International Security Standards. Do not use magnetic resonance imaging, because the product manufactured in stainless steel causes interference in the examination harming the diagnosis, and may also cause displacement of the implant causing trouble for the patient.

The use of Straight Plates for Osteosynthesis is designed to be an aid for the fixation and stabilization in the treatment of bone fractures in the lower and upper limbs, not to replace the normal structures of the human skeleton.

The surgeon should be aware of the physical development of the patient, a revision of the implant may be needed. The surgeon should also be familiar with and hold enough knowledge of osteosynthesis and its limitations, including pre and post operative surgical technique used, precautions and potential risks.

The surgeon should be attentive during bone preparation for the introduction of the FB Special Plates, because the stability and success of fixation may be compromised;

During handling of the implant scratches, grooves, or anything that may damage or mark the implant must be avoided, because these defects are stress concentrators and may be sites of crack nucleation and decrease the corrosion resistance and may result in fracture of the implant. In certain situations, it is necessary to mold the implant to perfectly fit to the anatomy of the bone to be treated. If implant modeling is needed, the bone anatomy is the modeling limit.

To ensure proper implantation, Biomecânica instruments must be used, because they are designed and manufactured specifically for the placement of these implants. Using instruments from different manufacturers may compromise the surgery as well as not respect product registration in the National Health Surveillance.

Do not use the product if it is damaged.

NON STERILE PRODUCT - STERILIZE PRIOR TO USE.  
PRODUCT OF SINGLE USE - IT IS STRICTLY PROHIBITED TO RE-USE THIS PRODUCT.

Special storage conditions: Store it in a cool, dry place, protected from light and weather. Do not use the product if the packaging is damaged.

**Note:** do not use implant components from different manufacturers. We, therefore, recommend that the products have the same provenance. Manufacture Date, expiration date and product batch number: SEE LABEL.

### Contraindications:

The use of these implants is contraindicated in the following situations: -sensitivity to metal;-fever or signs of local inflammation;-osteoporosis, obesity, poor-tissue coverage in the region to be operated;-pregnancy;- bone endargement because of diseases and infections, unexplained high sedimentation rate through other diseases, increase of the WBC rate or suspected striking changes in the count; visible infection;

### Adverse Effects

Mechanic loosening may be the result of a defective fixation, unstable reconstruction of the fractured bone or hidden infection. Metal sensitivity reactions in patients were rarely informed. Implantation of a strange material in tissues results in histologic reactions involving various sizes of macrophages and fibroblasts. The importance of this chemical effect is uncertain, but similar changes can occur as a

precursor, or during the healing process.

### Information to be provided to the patient

The patient must be informed about the importance of the postoperative follow-up. The lack of monitoring hampers the detection of problems after surgery. The failure to perform revision surgery when there is a component release or osteolysis may result in progressive loss of periprosthetic bone stock. The patient should be informed about the type of implant used, brand, code, name, batch number, surgeon's name and surgery date.

**Implanted Product Evaluations:** After implantation, in the intraoperative, the professional in charge must perform a radiological control to verify the correct positioning of the product. It is this professional's responsibility to make the clinical and radiological evaluation after the surgical procedure in the frequency stipulated by him/her to check the condition of the implant and the bone development. If the product is in the wrong position or presents any non-compliance, it is the surgeon's responsibility to take the appropriate corrective action.

Useful information to avoid risks derived from the implantation: To lower the risks the information contained in the contraindications, instructions for use and any information contained in "Product Instructions for Use" must be followed strictly.

**Implantation Intrinsic Risk:** Straight Plates for Osteosynthesis are made of materials of biomedical use, following the requirements of ASTM F138 -00 - Standard Specification For Wrought 18 Chromium-14 Nickel-2,5 Molybdenum Stainless Steel Bar And Wire For Surgical Implants (Uns S31673).

**Contamination Risks:** There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the metallic components of Biomecânica Straight Plates for Osteosynthesis come in contact with body tissues and fluids. Explanted products should be treated as highly polluting.

**Sterilization:** This product is supplied non-sterile. Before its use it should be sterilized. We recommend steam sterilization in autoclave in the hospital (ISO 17665-1: 2006 Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices).

**Product discard:** No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the Straight Plates for Osteosynthesis explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecânica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

**TRACEABILITY:** On the packaging of the implantable component there are 2 traceability labels containing the data of the implant used. These labels are the same label contained in the external product labeling. The information contained on this traceability label is related in the "Product Presentation".

This label must be affixed to the patient's records and another is for the hospital control. The hospital must have this information in its files and give them to the patient. We inform that it is mandatory to have the information in the patient's records that permits the traceability of the implanted product.

The attachment of this label in the patient's records enables the traceability of the product used. Among the main information we highlight as essential the name of the manufacturer of the implant used, name of the implant used, the implant code, batch number and registration number of the product in ANVISA. This information is described on the traceability labels that come with the product and is available on the external labeling. Other information must also be considered essential, such as, surgery date, name of patient who received the implant, surgeon's name, patient's weight, patient's age and other information requested in the patient's record must also be filled.

**MARKING:** Biomecânica Straight Plates for Osteosynthesis contain the following information recorded with laser to trace and monitor the patient after surgery:

Biomecânica Logo / Manufacturing Batch Number/ Material Symbol used in the manufacture of the product (S1). Code / Size / CE Mark

**Post-Sales (Customer Complaint):** If there is a need to make a complaint about Biomecânica Straight Plates for Osteosynthesis related to any adverse effect that affects the security of the user, such as problems with the performance of the product, damage to the implantable metal component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the national health authority and Biomecânica by e-mail or by phone sac@biomecânica.com.br 055-14 2104 7926. In cases of doubt, the surgeon in charge or health professional can make the communication of the adverse event through the Public Health Notification System at ANVISA website: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

**Sending the Product for the manufacturer's Analysis:** If it is necessary to send the implant for the manufacturer's analysis, it should be cleaned at the hospital using a bactericidal and broad-spectrum antifungal solution. Then, it should be disinfected or sterilized in a steam autoclave or with ethylene oxide. It should be sent to Biomecânica in an intact packaging, identified with the cleaning and sterilization methods, and product data.

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