FIOS E PINOS ESTÉREIS

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ENGLISH

Detailed description of the medical product as well as its functioning and action basis, content or composition, when applicable, and the list of the accessories intended to compose the product.

The non-absorbable rigid metallic wires and pins – BM are non-absorbable and non-active implants according to NBR ISO 14630 standard and they are dynamically used for external and internal osteosynthesis where destabilizing forces (which moves one fragment away from another) are changed into compression forces among the fragments. The non-absorbable rigid metallic wires and pins – BM are also used to guide the surgical instruments and implants in order to provide successful fixation.

They are available in several models in different dimensions and lengths, in cylindrical shape, smooth or threaded, either calibrated or not and with one of their bullet tips very sharp and resistant.

All models of non-absorbable rigid metallic wires and pins – BM are manufactured from a cold worked round bar in stainless steel and polished according to ASTM F138 standard. Its chemical, mechanical and metallographic composition is in accordance with the requirements of this standard as verified by tests. According to ASTM F138 this type of stainless steel which is used for the manufacturing of orthopedic implants has been successfully used in humans and in contact with tissues and bone. According to the standard it is unknown of an implantable material that hasn’t caused any reaction in the human body; however, this type of stainless steel has presented excellent results. The stainless steel ASTM F138 (18Chrome-14Nickle-2,5Molybdenum) is similar to the stainless steel described in the national standard - NBR ISO 5832-1. The standard - NBR ISO 21534 also highlights the stainless steel as being acceptable for the manufacturing of implants.

The models of non-absorbable rigid metallic wires and pins – BM are available sterile and sterilized in ETO gas or Gamma radiation.

Composition

The non-absorbable rigid metallic wires and pins – BM are manufactured in stainless steel according to ASTM F138 standard.

Ancillary Components

Diaphyseal external Fixator – Radio and Ulna; Articulated External Fixator for the fist; External Fixator – FBM 820 – Radio and Ulna; External Fixator – FBM 800 – Femur and Tibia; External Fixator – FBM 810 – Humerus; External Fixator for the proximal femur; Iliacus External Fixator - Anterior fixation, External Fixator for the metacarpus/ Proximal Phalanx/Distal; JS-Fix External Fixator.

Accessories

The non-absorbable rigid metallic wires and pins – BM do not have any accessory.

Admissible Combinations with other materials

The table bellow “Admissible contact alloys” was based on the standard - “NBRISO21534 - Non-active surgical implants – Joint replacement implants – Particular requirements” however, independently of the alloys mentioned in the item “Contact” it is important that the implants are combined with same brand implants because they have been designed for such combinations based on the surface finishing and treatment and other factors that may interfere with the combinations. Therefore metallic implants from different manufacturers are not recommended due to chemical, physical, biological and functional incompatibility.

<

nucleation and decrease the resistance to corrosion besides resulting in fracture or breakage of the implant.

A correct positioning and selection of the wires and pins is essential to optimize the bone fixation. The surgeon must be familiar with and have extensive knowledge of osteosynthesis and its limitations as well as pre and postoperative techniques, precautions and potential risks.

The surgical and pre-operative procedures as well as the familiarity with the surgical technique and the proper implantation and selection of the implant are important considerations for a successful use of a provisional device for internal fixation.

A periodic medical follow up is necessary in order to note potential modifications in the conditions of the implant and adjacent bone. Only a medical follow up is able to detect a potential looseness of the component or osteolisis.

Do not use the product in case it is damaged.

SINGLE USE PRODUCT – THE RE-USE OF THIS PRODUCT IS STRICTLY PROHIBITED.

Special storage conditions: Keep it in ventilated and dry places, far from the light and weather conditions.

Do not use in case the packaging is damaged.

Observation: Component from different manufacturers should not be used together. Therefore, we recommend that the products have the same provenance. Manufacturing date, expiration date and lot of the product: VIDE LABEL.

Do not use in case the packaging is damaged.

Adverse Effects
The mechanical looseness may be the result of a poor fixation or occult infection. Sensibility reactions to the metal in patients have been rarely reported

The implantation of a foreign material in soft tissues results in histology reactions that involve several sizes of macrophages and fibroblasts. The chemical importance of this effect is uncertain while similar changes may also occur as a precursor during the healing process.

Other potential adverse effects that may occur when a rigid wire or pin is implanted are: postoperative deep or superficial infection; vascular disorders such as thrombosis and pulmonary embolism; bleeding of the blood vessels and/or hematomas; pain, discomfort or abnornal sensations due to the presence or migration of the medical device or surgical procedure; injury to the sensory branch after percutaneous fixation with wires. Other adverse effect that may occur in osteotomy fixation with wires in the treatment of cubitus varus (elbow) is a deformity in the varus of the elbow, pain and limitation of the range of motion amplitude.

Information to be supplied to the Patient

In the items of this instruction: Indication, contraindication, information of use, potential adverse effects, precautions and warnings.

The patient must be informed as the importance of a post-surgical follow up. The failure to monitor the patient prevents the detection of post-surgical problems besides potential looseness of the components and osteolysis occurrence.

The patient must be informed as the type of the used implant, brand, code, name, lot, name of the surgeon and date of the surgery.

Every patient must be informed that before the healing of the fracture they should not walk without the help of crutches or walking stick in order not to impair the patient's recovery.

The patient must be aware that weight-bearing or muscle activity before the period recommended for the bone healing may bring complications such as the rupture of the device.

In case the patient is submitted to Magnetic resonance imaging he must inform the health professional that he has a metallic implant because the implant in stainless steel interferes with the test and impairs the diagnosis besides inducing to the detachment of the implant and bringing problems to the patient. The will and efficiency of the patient to fulfill the instruction of use is extremely important in orthopedic postoperative procedures. Children, elders, patients with mental problems and chemically dependents may generate a high risk of failure of the procedure once they are likely to ignore the instructions and restrictions.

The patient must also be informed that the component does not substitute and neither performs as a healthy bone. Therefore it may deform, break or loose due to stress, excessive physical activity or load on the bone permanently.

Important

The surgeon must orient the patient towards the cares to be taken and the proper physical support during the period of the implant use. Before a complete stabilized fracture, the load may be released. If the fracture is not stabilized, no load can be applied to the implanted system.

Responsibility exemption

Biomechanica as the manufacturer of these devices does not practice medicine nor recommend this or any other surgical technique to be applied to a certain patient. The surgeon that follows any implant procedure is responsible for defining and using the appropriate techniques for the implantation of the medical device in each patient. Biomechanica is not responsible for the choice of the appropriate technique applied to each patient. The surgical technique depends on the technique adopted by the surgeon. The surgeon is also responsible for the choice of the method, type and dimension of the products to be implanted.

Evaluations of the Implanted Product

The professional in charge must do the radiological control after the implantation in the intraoperative in order to verify the correct positioning of the product. The professional in charge must do the clinical and radiological evaluations after the surgical procedures in a frequency stipulated by him in order to verify the implant condition and evaluate the bone healing. The surgeon is responsible for taking the most appropriate corrective actions in case the product is out of the correct positioning or presents any non-conformity.

Useful information for avoiding the risks from implantation

For minimizing the risks resulted from implantation it must be strictly followed: contraindications, instrction of use and all information in the "instruction of use" of the product.

Intrinsic Risks of Implantation

The non-absorbable metallic wires and pins-BM are manufactured with materials of well known biomedical use which comply with the requirements of the standard - ASTM F138 -08- Standard Specification For Wrought 18 Chromium-14 Nickel-2.5 Molybdenum Stainless Steel Bar And Wire For Surgical Implants (Uns S31673).

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There are risks of biological contamination and transmission of viral diseases such as HIV and Hepatitis because the non-absorbable metallic wires and pins-BM get in contact with the tissues and body fluids. The explanted products must be treated as highly contaminants.

Sterilization

This product is supplied sterile. Depending on the client's request it can be sterilized as follows:

Process	Expiration
Sterilization by Ethylene Oxide Gas (ETO)	05 years
Sterilization by Gamma radiation (Gamma Ray)	05 years

Obs.: The difference in the type of sterilization will be defined by the code of the product. The products sterilized in ETO gas have different codes from the products sterilized by Gamma Ray. If the expiration date of the sterilization of the product is expired or the package is damaged, the product must be returned to the responsible distributor or directly to Biomechanica.

Handling of the sterilized material:

When handling the sterilized material by aseptic technique some Standards must be obeyed in order to keep it sterile. It is fundamental to wash the hands with water and soap before handling the sterilized material; the packaging of the material must be intact, dry, unmarked and identified (type of the material and sterilization date); Work facing the material; handle the material at waist level and up, avoid coughing, sneezing or speaking over the exposed material; do not make movements on the sterilized area; check the expiration and packaging adequacy; work in a clean, calm and dry environment where there is no air current; keep a reasonable distance between the body and the material to be handled; follow the other principles of asepsis.

The nursing technique recommended for the handling of the sterilized material is: open it from the opposite edges of the handler, protect the exposed material with the sterile field involving it; touch only the outside of the packaging; do not consider as sterile a package that was previously open.

Guarantee

The guarantee of the non-absorbable rigid metallic wires and pins - BM is related to the observance of this instruction of use. They are: indications and information of use, contraindications, cares and precautions, warnings, potential adverse effects, packaging, sterilization, cleaning and decontamination recommended in this technical report and instruction of use..

Product discard

No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the non-absorbable rigid metallic wires and pins – BM explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Improper for use". According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

Traceability

In the packaging of the implantable component there are 6 traceability tags with the implant data. This tag is the same as the one of the external labeling of the product. The information available in this traceability tag is described in item "Presentation" in this technical report. The tags must be distributed as follows: one must be pasted on the patient's clinical report; one on the report that is handed to the patient; one on the tax documentation that generates the collection; on the AHI in case the patient is supported by SUS or on the invoice in case the patient is supported by a complementary health system; one is available for the control of the supplier (Historical distribution register – RHD); one for the control of the surgeon in charge (main) and the last one for the health insurance in case there is one. It is mandatory the patient's report to have registered all information necessary for the traceability of the implanted product. The pasting of the tag on the patient's report will allow the traceability of the product. Among the main information we consider indispensable the name of the manufacture of the implant, the code and name of the implants besides the lot and register number with ANVISA. Such information is described in the traceability tag that comes with the product and on the outer labeling. Other information that may also be considered indispensable is: date of the surgery, name of the patient that received the implant, name of the surgeon, weight and age of the patient and other information required in the patient's report must be provided.

HILOS Y CLAVIJAS ESTEREIS

MARKING

The non-absorbable rigid metallic wires and pins –BM present the following information laser engraved according to NBR ISO15165 standard to allow the traceability and patient's follow up in the postoperative.

Biomechanica logo / Number of the manufacturing lot / Acronym of the raw-material used in the manufacturing of the product (S1). / Code / Dimension / Ce Mark

Post-sale (Client Complaint)

In case it is necessary to make a complaint about the non-absorbable rigid metallic wires and pins- BM related to some adverse effect that may affect the safety of the user, such as the product is not working or it is damaged, serious problems or death related to those components, the surgeon in charge must inform the authorized health organization and Biomechanica about this adverse event through the email sac@biomechanica.com.br or by the phone 0xx14 2104 7900. In case of doubts, the surgeon in charge or the health professional may inform NOTIVISA about the adverse event on the website of ANVISA: http://www.anvisa.gov.br/hotsite/notivisa/index.htm

Sending of the material for the manufacturer to analyze

In case the implant is sent to the manufacturer for analysis, it must be sanitized at hospital by a bactericide and antifungal solution of broad spectrum. Then it must be disinfected and steam sterilized in autoclave or ethylene oxide. It must be sent to Biomechanica in intact packaging and identified with the cleaning and sterilization methods and the data of the product.

Description of the efficiency and safety of the medical product in accordance with the regulations of ANVISA that determines the Essential Requirements of Efficiency and Safety of the Medical Devices.

The risks associated with the non-absorbable rigid metallic wires and pins – BM are in accordance with the Essential requirements of Efficiency and Safety of the Product according to RDC 56/01 for the products classified in rule 8, class risk 3. The risk management of the product and the process are in accordance with NBRISO14971 as foreseen in PRCE02 procedure – Risk Analysis of Medical Devices.

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