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ENGLISH

PRODUCT DESCRIPTION

The Implant System for Spine Fixation- Pedicular - BM was projected to use in the Spine by via posterior. It is a Group of metallic implantable and modulates components, compatible amongst themselves and anatomical to the bone, used for: Correction, Stabilization, and Fixation permanently or temporarily the segments thoracic, lumbar and sacral of the spine. Having Pedicular Poliaxial screw (articulated), Pedicular Screw, Pedicular Spondylolisthesis Screw, Distraction Stem, Transverse Stem, Connector, Screw Locks Hook, Hook Normal Distractor, Hook Left Distractor Straight, Hook Lateral Left Distractor, Hook Distractor Sloping Right, Hook Distractor Lateral Right, Hook Straight Long Distractor, Hook Distractor Pedicular and Contra -Screw. All the implants are supplied separately. The finish of the implants is by mechanical polishing afterwards this product goes by a residual removing process, that consists of chemical attack of the polished surface of the metal that removes the residual sludge´s of the production process, such as oils, greases and other production materials afterwards this product goes by an ionization process, in other words, superficial treatment to strengthen the superficial layer of the products manufactured in titanium.

PRODUCTION MATERIAL

The products are made through the titanium alloy 6Al 4V ELI NBR ISO 5832-3 or ASTM F-136. It is a biocompatible alloy and it possesses resistance to the corrosion. The resistance to the fatigue and high resistance to the corrosion is intrinsic mechanical properties of 6Al 4V ELI

PURPOSE

It is indicated in the treatment of Lumbar and Thoracic Fractures, Via Posterior , disk hernia, pseudarthrosis, stenosis, Spondylolisthesis, fracture sequels , resection of Tumors, cifosis, idiopathic and congenital Scoliosis and neuro-muscular deformities, Lordosis, Degeneration of Intervertebral Disc and surgical revisions.

PEDICULAR FIXATION SYSTEM

ANCILLARY COMPONENTS

The Implant System for Spine Fixation- Pedicular BM doesn´t possess ancillary components.

ACCESSORY

The Implant System for Spine Fixation- Pedicular BM doesn't possess any accessory with the purpose of integrating the medical product, granting to this product a function or complementarily technical characteristic.

MATERIAL OF SUPPORT

The support materials that accompany the implants of Implant System for Spine Fixation- Pedicular BM are:

- Bull of the product. Use instructions. SIFPC (revision 0);
- Surgical technique (optional). (Revision 0);
- Box of surgical instrumental (optional). COD.:7000

FORM OF PRESENTATION

The Implant System for Spine Fixation- Pedicular BM, is wrapped in plastic envelope, identified with no sterile label, external labeling, containing 1 implant unit, 1 instruction of use and 1 internal label to be labeled in the patient's medical file. The hospital institution is responsible in assuring the identification and the traceability of the medical product, through the data and the necessary information, contained in the labeling (label) of the medical product. The labeling contains the following information:

Manufacturer Data | Description of packing content | Measure | Raw material | Unit | Reference | No Sterile | Sterilization Method | Lot | Amount | Production Date | Validity | Medical Article of unique use | Instructions, Warnings and Precautions: See instructions of use | Technical responsibility | Health Ministry Registration

The patient should be informed as the importance of the post surgical attendance. The lack of attendance impedes the detection of surgical problems.

The patient should be informed as for the type of the implant utilized, marks, code, name, lot, the surgeon’s name and surgery dates.

WARNINGS

Product of Unique Use. Destroy after explanted. Not to reuse the product. An orthopedic implant can only be used in a single patient, a single time. Although it can seem not to be damaged, previous tensions can create imperfections that can reduce the success of the implant. Inappropriate selection of implant can cause uncommon tensions in the implant and they can result in subsequent fracture of the implant.

Due to his presentation form and their characteristics, this product doesn't present side effect, however, we suggests monitoring the patient's re-establishment through the X-Ray, according to International Norms of Safety. Attention when using magnetic resonance in patients that already have column implants manufactured in stainless steel, because it cause interference in the exam harming the diagnoses.

The use of the Implant System for Spine Fixation Pedicular BM is conceived for correction, stabilization and fixation of the segments thoracic, lumbar and sacral of the spine and not to substitute the normal structures of the human skeleton.

Obesity has effect in the load in the spinal system could commit the fixation of the implant for column. Excessive physical activity and traumas called directly in the segments thoracic, lumbar and sacral can commit the success of the implant prematurely. We recommended that the activities are dosed avoiding the excessive “stress” in the operated area. The Implant System for Spine Fixation Pedicular BM cannot suffer loads or levels of similar activities supported in normal and healthy bones.

The mixture of different biomedical materials no compatible and/or similar it can elevate the corrosion risk, taking into account also that in spite of the implant to be manufactured by material biocompatible when implanted in the patient the same is subject to constant environmental changes promoting the concentration of salts, acids and alkalisols that it can promote the corrosion. The clinical results and durability of the implant depend on the existence of a surgical technique needs.

He/She dates from production, period of validity and lot of the product: VIDE LABEL.

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INSTRUCTION OF USE
The results depend on the technique adopted by the surgeon, fitting to this the choice of the method, type and dimension of the products to be inserted.

The mixture of different biomedical materials no compatible and/or similar it can elevate the corrosion risk. The clinical results and durability of the implant depend on the existence of a surgical technique needs.

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ASSEMBLY AND DISMOUNT

The placement of the Implant System for Spine Fixation- BM is accomplished with aid of the Instrumental Kit of Column. Biomecanica

Preparation for screws planting
1. Use the Perforator 8509 for rip off the point of entrance of the pedicle, exposing the spongy structure. The point of entrance of the pedicle is intersected by the vertical line which connects the lateral borders of the bone crest extension of intracurular segments, and the horizontal line which divides the transverse process.
2. Controlling the applied pressure in a series of rotating and firm movements, the meter 8510 is inserted in the pedicle with angulation from 10 to 20°. This meter will follow the itinrary of smaller resistance (spongy structures).
3. The feeder guide 8514 is used to verify the presence of bone plan in all of the directions. The tip of the instrument should slide in contrary sense to the cortical wall in all the plans.
4. In case of doubts, X ray should be accomplished with help of the meter pin 8513.
Insert of the screw
5. Use the square reamer tip 8511 to determine the appropriate length of the screw. The extremity of pin depth meter is put in the channel of the pedicle and it is advanced through the channel until the previous wall of the vertebral body.
6. The screws are inserted using the Initial Key Guide 8524. The tip of the key is inserted in the compartment of the pedicular screw and tight locking the system.
OBS: The pedicular screws are auto-reaming, usually doesn't request thrust preparation. However, the males for this preparation are available in the Instrumental Kit of Column- Biomecanica. It is important to remind that the pedicular screws should not be removed for the depth adjustment. In consequence, should introduce each aligned screw in relation to the adjacent ones. This way checks that the distraction stem will follow a regular line in the correct form.

Insert of the distraction stem
1. After placement of pedicular screws, the distraction stem is selected. The same can be modeled with the aid of Stem Modeler 8518, when necessary. The stems are supplied in varied lengths and they should be extended 5 mm approximately besides the external borders of the screws bodies of superior pedicle, and more inferior.
8. The Fillster Stem Guide 8502 can be used to stabilize the position of the Distraction Stem as the blockade screw of pedicular screw was not inserted.

Insert of transverse stem
9. In certain situations the transverse stem should be applied to provide rigidity to the system and to distribute the loads better among the vertebrae pedicles. The size of the transverse stem should be defined so his compression is appropriated, as well as it can promote the appropriate fixation of the whole group.
10. Use the transverse hook coupled to the transverse stem to establish the fixation of the group.

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RISKS OF CONTAMINATION

Risks of biological contamination and transmission of viral diseases exist as HIV and hepatitis because the metallic components of the Implant System for Spine Fixation Pedicular -BM gets in touch with tissue and corporal fluids. Products explanted should be treated as highly pollutants.

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PRODUCT DISCARD
No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the metallic components of Implant System for Spine Fixation Pedicular - BM explanted from the patients should be properly discarded by the hospital. The complete mischaracterization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomecanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings “Improper for use”. According to Resolution nº 2605 from 08/11/2006, it is forbidden to reprocess implantable devices of any nature classified as single use.

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STERILIZATION

This product is supplied no sterile. Before use should be sterilized. We recommended the steam sterilization in hospital Autoclave (ISO 11134: 1994 *Sterilization of health care products Requirements for validation and routine control Industrial moist heat sterilization.*)

Autoclave
It is equipment for sterilization through saturated steam under pressure.

Picture: Sterilization of the Material for the steam process saturated under pressure (Autoclave), In agreement with the temperature exposition according to (Resolution SS-374, de 15-12-95)

Temperature	Time of Exposition			Observation
Material	120° C	130°C	132° C (High Values)	
Material of stainless steel or other metal thype	30min	15min	4min	Wrapped in paper kraft or field of fabric of double raw cotton.
Metallic instruments inserted on trays or metallic box perforated	30min	15min	4min	Wrapped in paper kraft or field of fabric of double raw cotton. Contra indicated the use of boxes semi open
Cutting blades scissors and saw	30min	15min	4min	Wrapped up individually or conditioned in tray or perforated and wrapped up metallic box in paper kraft or field of fabric of double raw cotton.

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Operational instructions
For departure of the equipment, verify: / if the interrupter is linked; / if the water register is open; / if the discharge register is closed.

1. Open the door of the equipment;
2. Accommodate the material to be sterilized appropriately;
3. Close the door of the equipment;
4. Select the cycle wanted in agreement with the material to be sterilized;
5. Tie the general key;
6. Cycle will elapse automatically, in the sequence;
7. When lighting the “final lamp of cycle”, open the door partially for ten (10) minutes, approximately, for cooling the material.

FOR BETTER EXPLANATION CONSULT THE MANUAL OF INSTRUCTIONS WHICH ACCOMPANIES EACH AUTOCLAVE.

Another sterilization method: besides AUTOCLAVE can be used the following method of sterilization in accordance defined by hospital institution:

Sterilization by oxide of ethylene (E.T.O.) - parameters and procedures established in the validation protocol and in ENS50:1994 *Sterilization of Medical Devices Validation and Routine Control of Ethylene Oxide Sterilization.*

NOTE

It is of entire care the hospitaliar institution responsibility, for the sterilization method, the equipments, the controls, and the instructions about the sterilization method.

CLEANING AND DECONTAMINATION RECOMMENDED

General recommendations

To guarantee the reduction of the infection risks, the implants and surgical instruments must be cleaned after the surgery and before the sterilization. The products should have decontamination using a bactericidal solution and fungicide of wide spectrum, avoiding aggressive agents of cleaning, metal brushes, so the quality of the metallic implants is not harmed. The manual cleaning, always is more recommended by attacking less the implants and surgical instruments. In the manual cleaning, water should be used at room temperature, t use nylon brushes, never of steel or abrasive, because it can harm the protecting layer of the material.

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In case of cleaning machines the implants and surgical instruments should be positioned in drawers, to avoid hit one in another, avoiding damages in the material. It is not recommended to put together materials of different composition during the wash in cleaning machines f.

The responsible for cleaning the implants and surgical instruments should be attentive as for the type of cleaning product used respecting the suitable concentrations by the manufacturer. The implants and surgical instruments should be rinsed carefully until every residue comes out dry immediately after the cleaning

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RE-STERILIZATION

In case of re-sterilization, all the cares above described should be assisted.

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TRACEABILITY

Near the packing of the implantable component, follows one more label containing the data of the used implant. This label should be glued in the patient’s file. The hospital should register that information in their files and to review that patient information. In the file the following information are indispensable:

- Name of the implant / - Surgery dates / - code of the product / - lot number / - patient’s name / - surgeon’s name

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DEMARCACTION

All the implants mentioned in this instruction of use, count the following information recorded by laser to allow the traceability of the product:

- Manufacturer’s logo / - size / - lot / - reference

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SPECIAL CARES IN THE STORAGE AND TRANSPORT OF THE MEDICAL PRODUCT

- Do not use the product if he looks to be damaged. / - The implants should be manipulated carefully. Avoiding damages that can harm the quality of the material and the patient’s safety / -Conserve at air place, dry, to the shelter of the light and far away from the action of bad weather. / - When transporting t, shocks, vibrations and high temperatures should be avoided (above 45°) and inadequate piling up.
* The cares to be taken during the transport are described in the collective packing.

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OBSERVATION

This product should only be manipulated by qualified people.

Components of different manufacturers’ implants should not be used. Therefore we recommended the products have the same origin.

Every implant is of unique use. Therefore, when withdraw from patient, should be clean and discarded.

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COMPLAINT

In case of some complaint of Implant System for Spine Fixation BM these should be reported to Biomecanica authorized distributor. In the case the complaint is related to some adverse effect that affects the user’s safety, as product not working, harm of the metallic component, serious problems or death related with those components Biomecanica should be communicated.

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SENDING THE MATERIAL FOR MANUFACTURER ANALYSIS

In case of sending implants for manufacturer accomplish analysis, this should be higienized in the hospital using a bactericidal solution and fungicide of wide spectrum. Afterwards it should be disinfected or sterilized by oxide of ethylene. They should be dispatch to Biomecanica in integrate packing, identified with the method of cleaning, sterilization and data of the product.

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