Due to its form and characteristics, this product does not have side effect, however, it is suggested to monitor the patient's recovery through X-ray according to the International Security Standards. Do not use magnetic resonance imaging, because the product manufactured in stainless steel causes interference in the examination harming the diagnosis, and may also cause displacement of the implant causing trouble for the patient.

The use of FB Special Plates for Osteosynthesis are designed to be an aid to the fixation and stabilization in the treatment of bone fractures in the lower and upper limbs, not to replace the normal structures of the human skeleton. The surgeon should be aware of the physical development of the patient, a revision of the patient's medical history. The surgeon should also be familiar with and hold enough knowledge of osteosynthesis and its limitations, including pre and post operative surgical technique used, precautions and potential risks. The surgeon should be attentive during the bone preparation for the introduction of the FB Special Plates, because the stability and success of fixation may be compromised;

To ensure proper implantation, Biomechanics instruments must be used, because they are designed and manufactured specifically for the placement of these implants. Using instruments from different manufacturers may compromise the surgery as well as not respect product registration in the National Health Surveillance.

Do not use the product if it is damaged. NON STERILE PRODUCT - STERILIZE PRIOR TO USE. PRODUCT OF SINGLE USE - IT IS STRICTLY PROHIBITED TO RE-USE THIS PRODUCT.

Special storage conditions: Store it in a cool, dry place, protected from light and weather. Do not use the product if the packaging is damaged. **Note:** do not use implant components from different manufacturers. We, therefore, recommend that the products have the same provenance. **Manufacture Date, expiration date and product batch number:** SEE LABEL.

Product Conforming Limit: During handling of the implant scratches, grooves, or anything that may damage or mark the implant must be avoided, because these defects are stress concentration and may sites of crack nucleation and decrease the corrosion resistance and may result in fracture of the implant. In certain situations, it is necessary to mold the implant to perfectly fit to the anatomy of the bone to be treated. If implant modeling is needed, the bone anatomy is the modeling limit.

Loading Restrictions / Mechanical Strength: For the surgery success it is essential to correctly select of implants. It is extremely important the correct handling of the implants.

The FB Special Plates for osteosynthesis are manufactured to help the orthopedic surgeon in internal fixation procedures. Do not expect these implants the same performance of a healthy bone or that it supports intense and continuous mechanical stress. The performance of the implant is related to factors such as body weight, physical activity, obedience to the instructions given by the doctor after surgery, especially in consolidation absence or delayed bone consolidation

Weight, degree of activity and the patient's conditions may impact the performance of the implant. The fact that no implant is as strong as natural bone and, therefore, has limitations in biomechanical applications should be considered.

The implant mechanical contribution should be examined by observing the exact adjustment and compression, because a stable reconstruction of the fractured bone reduces the load to be supported by the implant. Another important factor is related to the implant placement.

Adverse Effects: Mechanic loosening may be the result of a defective fixation, unstable reconstruction of the fractured bone or hidden infection. Metal sensitivity reactions in patients were rarely informed.

Implantation of a strange material in tissues results in histologic reactions involving various sizes of macrophages and fibroblasts. The importance of this chemical effect is uncertain, but similar changes can occur as a precursor, or during the healing process.

Information to be provided to the patient: The following information must be provided to the patient: indications, contraindications, information for use, Possible Adverse Effects, precautions and warnings.

The patient must be informed about the importance of the postoperative follow-up. The lack of monitoring hampers the detection of problems after surgery. The failure to perform revision surgery when there is a component release or osteolysis may result in progressive loss of periprosthetic bone stock.

The patient should be informed about the type of implant used, brand, code, name, batch number, surgeon's name and surgery date.

Contamination Risks: There are risks of biological contamination and transmission of viral diseases such as HIV and hepatitis because the metallic components of the FB Special plates for Osteosynthesis come in contact with body tissues and fluids. Explanted products should be treated as highly polluting.

Decontamination, Cleaning and Sterilization: The product is supplied in its original pre-cleaned packaging, ready to be sterilized. It must be handled with care to avoid contamination. If the product is submitted to conditions that compromise its cleanliness this implant must be cleaned before the sterilization process. Below, we suggested some cleaning procedures: Cleaning can be developed through manual or mechanical methods.

Cleaning Equipment:

•Ultrasonic Washer / Sterilizing Washer / Washer – disinfectant

Manual Cleaning: is the procedure which the dirt is removed by means of physical action with the aid of detergent, water and artifacts such as sponge and brush. Manual cleaning is recommended for causing less harm to implants and surgical instruments. In manual cleaning water should be used at room temperature, using nylon brushes, never steel or abrasive brushes, because they can damage the protective layer of the material.

In the decontamination process, antibacterial and broad-spectrum antifungal solution should be used; aggressive cleaning agents, brushes, metal must be avoided not to damage the metallic implants. In case of cleaning machines, implants and surgical instruments should be organized in drawers not to hit each other and avoid material damage. People responsible for cleaning the implants and surgical instruments should be aware of the type of cleaning product used and respect the concentration recommended by the manufacturer. Implants and surgical instruments must be thoroughly rinsed until all residues are removed and dried immediately after cleaning.

Types of cleaners

Enzyme Cleaners: consist mainly of enzymes, surfactants and solubilizers. A balanced combination of these elements means that the product can remove organic material in a short period of time.

Enzyme solutions: provide excellent cleaning action, but do not have bactericidal and bacteriostatic activity.

Enzymes: are substances produced by living cells and which cause chemical reactions of the process. Once produced by the cells, an enzyme can be isolated and will maintain its catalytic properties, if certain conditions are maintained in its manufacture. Enzymes are classified into three major functional groups depending on the type of substrate that they will affect: proteases, amylases and lipases that act on proteins, substrates, fats and carbohydrates, which tend to solubilize and come off of the products.

Disinfection: is a process where pathogenic or not microorganisms are killed, except for bacterial spores, by physical or chemical methods.

Levels of disinfection:

High level: destroys all microorganisms, with the exception of high numbers of bacterial spores => Glutaraldehyde 2% -20 -30 minutes.

Indication: preferably in hospitals

Intermediate level: eliminates vegetative bacteria, most viruses, fungi and mold/spores => Sodium hypochlorite 1% -30 minutes.

Indication: day care, rest homes, nursing homes.

Low-level: eliminates most bacteria, some viruses and fungi, but does not eliminate microbacteria => Sodium hypochlorite 0.025%

Indication: nutrition.

Sterilization: This product is supplied non-sterile. Before its use it should be sterilized. We recommend steam sterilization in autoclave in the hospital (ISO 11134: 1994 Sterilization of health care products - Requirements for validation and routine control - Industrial moist heat sterilization).

AUTOClave: is sterilizing equipment that uses saturated steam at high pressures.

OPERATING INSTRUCTIONS: To start the equipment, check: if the circuit breaker is switched on; if the drainpipe is open; if the flush valve is closed.

PLACAS ESPECIALES - SPS-FREE-BLOCK

1 - open the door of the equipment; 2 - Load the autoclave properly; 3 - close the equipment door; 4 - Select the desired cycle according to the material to be sterilized; 5 - turn on the key; 6 - the cycle will run automatically, in sequence; 7 - when the light turns on "final cycle", partially open the door for ten (10) minutes or so, for cooling the material down.

It is recommended to apply the following physical parameters of sterilization in autoclave (saturated steam):

Note: Register the time when the heat of the sterilization chamber reaches the desired temperature.

FOR CLARIFICATION CONSULT THE AUTOCLAVE MANUAL.

Another method of sterilization: can be used in addition to the following autoclave sterilization methods as defined by the hospital: ethylene oxide Sterilization (ETO) - parameters and procedures established in the validation protocol and EN550: 1994 Sterilization of Medical Devices - Validation and Routine Control of Ethylene Oxide Sterilization. **NOTE:** The method of sterilization, equipment, controls, and instructions for sterilization used is under the hospital responsibility.

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Conditions of storage of Sterilized Products: Storage Conditions of sterilized products must be clearly dry. It should be a restricted area.

The product: After the sterilization process, do not place it on a cold surface (stone or stainless steel), use baskets or containers containing holes until they cool down; mantle (raw cotton fabric, non-woven fabric, surgical grade paper, crepe paper, paper film, TYVEK or perforated metal boxes) must remain intact and be handled slightly to prevent the packages from tearing or the seal to come out; be stored in closed cabinets with shelves identified shelves in order to facilitate the removal of the material; material should be stored according to the maturity date of sterilization to facilitate the distribution and not keep expired material in stock; store them separately from non-sterile products to reduce the level of external contaminants.

Product discard: No surgical implant should be re-used. Any metallic implant should be discarded after use. Even though the implant doesn't look damaged it may present little defects that may lead to contamination or its breakage due to fatigue. This way the FB Special Plates for osteosynthesis explanted from the patients should be properly discarded by the hospital. The complete sterilization of the implants in order to prevent their re-use is under responsibility of the hospital. The method for the mischaracterization of the implant is also under responsibility of the hospital. Biomechanica recommends the implants to be mechanically destroyed with the help of a hammer or press and immediately identified with the sayings "Impropor for use". According to Resolution n° 2605 from 08/11/2006, it is forbidden to re-use implantable devices of any nature classified as single use.

TRACEABILITY: On the packaging of the implantable component there are 2 traceability labels containing the data of the implant used. This label must be affixed to the patient's records and another is for the hospital control. The hospital must have this information in its files and give them to the patient. It is essential to have the following information in the patient's records: name of the implant used * surgery date product code * batch number name of patient who received the implant surgeon's name.

Post-Sales (Customer Complaint): If there is a need to make a complaint about Biomechanica Special Plates for Bone Synthesis related to any adverse effect that affects the security of the user, such as problems with the performance of the product, damage to the implantable metal component, serious problems or deaths related to these components, the surgeon in charge must report this adverse event to the national health authority and Biomechanica by e-mail or by phone at biomecanica.com.br or 55 14 2104 7900. In cases of doubt, the surgeon in charge or health professional can make the communication of the adverse event through the Public Health Notification System at ANVISA website: Sending the Product for the manufacturer's analysis: If sending the implant for the manufacturer's analysis, it should be cleaned at the hospital using a bactericidal and broad-spectrum antifungal solution. Then, it should be disinfected or sterilized in a steam autoclave or with ethylene oxide. It should be sent to Biomechanica in an intact packaging, identified with the cleaning and sterilization methods, and product data.

Expected performance of the product: Every implant is subject to breakage, wear or loosening and its dimension is determined for a number of requests for a certain time.

The time predicted in traumatology surgeries is until the fracture is consolidated. These requests are related to stress caused by normal activities, such as, walking, getting up or climbing stairs. FB Special Plates for Bone Synthesis were tested to check their load limits, thus presenting fixation values. They are manufactured in stainless steel, according to ASTM F139 and the locking ring manufactured in accordance with ASTM F136 approved through laboratory test reports and quality-product certification. These materials are acceptable for the manufacture of implants as defined in Annex A of NBRISO21534.

According to NBRISO21534, in Annex C.1 the metal combination of non-articulating contacting surfaces between the stainless steel NBRISO5832-1 similar to the ASTM F138 titanium alloy ISO 5832-3 similar to the ASTM F136 are considered suitable combinations.

The manufacture and packaging of FB Special Plates for Bone Synthesis are based on specific procedures. They are distributed non-sterile and packaged in plastic envelope of transparent non-toxic PVC, sealed in a sealing machine in order to ensure its integrity. FB Special Plates for Bone Synthesis are distributed with sufficient information for the correct handling of the product including the safety of people involved in handling it.

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