MINI TWISTER – ÂNCORA DE SUTURA

PORTUGUÊS

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MINI TWISTER – SUTURA ANCHOR

ENGLISH

Quantity and Description of the Traceability Labels that Comes with the Product:
- 6 adhesive labels of identification and traceability of the product;

Description of the Components, Quantity and Raw Material:

- 1 unit of titanium alloy Anchor Screw Ti-6AL-4V; 1 unit of White / Blue Suture Thread of polyethylene fiber UHMWPE; - 1 unit of Stainless Steel Anchor Screw Applicator AISI 304, Polypropylene and Silicone.

Product Presentation:


The product is packed in a blister-type double package, which is manufactured from non-toxic PET (Polietileno Tereftalato), allowing a contact with hospital products. The two blister packages are sealed with a surgical grade paper, without dye and free of holes, through the heat sealing process. The double sealed blister package is also placed in an outer cardboard package, covered with security seal. The implants are sterilized in ETO with 5 years validity.

Traceability:

Along with the product, there are 6 labels containing the necessary data for traceability of the product (Manufacturer's name, trade name of the product and code, technical product name and code, batch number, registration number at ANVISA, quantity, manufacturing date, expiration date, manufacturer's address, technical responsible name and the words: “sterile product,” “do not reprocess,” “single-use product, do not reuse”, “before using read the label instructions” and “warnings / precautions / special care / symbology: see instructions for use”).

A label must be mandatorily fixed to the patient's clinical record, another one in the report delivered to the patient, other in the tax documentation generating the collection, at AIH, in case of a patient treated by SUS (Brazilian Unified Health System), or in the sales invoice, in case of a patient cared by the supplementary health system, another one made available to the supplier's control (historical record of distribution - RHD, another one made available to the control of the responsible surgeon (principal) and the last one to the health insurance, if any, in order to keep the traceability cycle of the product.

Only the Applicator for Anchor Screw 2.0 / 2.7, present in *Mini Twister - Suture Anchor*, receives the laser marking, which is performed in the stainless steel manufactured part, according to the sequence below:

The other components of *Mini Twister - Suture Anchor* do not receive marking due to their anatomy. Even though, the traceability of all components of *Mini Twister - Suture Anchor* is assured through the batch number present in the traceability labels of the product, which can be identified by tracing the batch number of each component.

Product Description

Mini Twister - Suture Anchor consists on a suture medical device, not active, composed of three components, two non-absorbable implants and an instrumental (applicator), assembled and ready for use.


Figure 1: Mini Twister - Suture Anchor

The single component of *Mini Twister - Suture Anchor* showing the variation is the Anchor Screw, and may have a diameter of 2.0 mm or 2.7 mm. In the presence of the Anchor Screw 2.0 mm, the code of *Mini Twister - Suture Anchor* is 3128-20-000; in the presence of the Anchor Screw 2.7 mm, the code of *Mini Twister - Suture Anchor* becomes 3128-27-000.

Both the Anchor Screw Ø 2.0, as the Anchor Screw Ø 2.7, are implants, not absorbable and not active, manufactured from titanium alloy 6-aluminum 4-vanadium (Ti6Al4V), according to the ASTM F136 standard.

Both the Anchor Screw Ø 2.0 and the Anchor Screw Ø 2.7 have a pointed end, which assist the introduction of the screw in the bone wall. The threaded part causes the screw to be introduced in a safe and practical way, providing strength and stability for screw fixation in the bone.

Both the Anchor Screw Ø 2.0 and the Anchor Screw Ø 2.7 have a polished finishing.


Figure 2: Anchor Screw

Table 1: Dimensional Comparison between the Anchor Screw Ø 2.0 and Anchor Screw Ø 2.7

	Thread Diameter	Thread Length	Pointed End Length	Fitting Area Length of the Thread and the Applicator	Total length
Anchor Screw Ø 2.0	2.0 mm	4.5 mm	1.5 mm	3.0 mm	9 mm
Anchor Screw Ø 2.7	2.7 mm	6.0 mm	2.8 mm	3.2 mm	12 mm

The White / Blue Suture Thread is a non-absorbable implant of synthetic origin for surgical suture, composed of a fiber called Dyneema Purity with a high mechanical resistance based on ultra-high molecular weight polyethylene (UHMWPE); multifilament braided; and not active.


Figure 3: Suture Thread

The White / Blue Suture Thread has a length of 900 mm and an average diameter of 0.585 mm, and can vary from 0.470 mm to 0.705 mm.

The Applicator for Anchor Screw 2.0 / 2.7 is an instrumental, non-active applicator, and acts as an applicator of the Anchor Screw Ø 2.0 or Ø 2.7 and the White / Blue Suture Thread.

The Applicator for Anchor Screw 2.0 / 2.7 is manufactured from Stainless Steel Type 304, according to the ASTM A269 and ASTM F899 standard. It also presents a cable made by a plastic injection of Polypropylene. The total length of the injected part is 120 mm (±1) and the largest diameter of the injected part reaches 35 mm (±1). The cable also has a silicone ring with the purpose to attach the suture threads.


Figure 4: Applicator for Anchor Screw 2.0 / 2.7

The stainless steel manufactured part has a polished finish.

The stainless steel manufactured part is a tube with a total length of 185 mm (±1) and an outer diameter of 3.26 mm. At one end, the tube has a hexagonal socket for Screw Anchors Ø 2.0 or Ø 2.7. The White / Blue Suture Thread, attached to the Anchor Screw Ø 2.0 or Ø 2.7, go through the tube.

Ancillary Components

Mini Twister - Suture Anchor does not have ancillary components.

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Indications, Contraindications, Adverse Effects and Warnings

Indications

Mini Twister - Suture Anchor is indicated for the tendon-bone repair of the flexor tendons of the hands.

Average Strength of the Anchor Screw Ø 2.0 to its Loosening in Relation to the Bone	Average Strength of the Anchor Screw Ø 2.7 to its Loosening in Relation to the Bone	Average Strength to Traction Supported by the Suture Thread to Break	Strength Required by the Flexor Tendons of the Hands, According to the Scientific Literature
88,083N	191,884N	142,472N (in the knot) 235,762N (with no knot)	2.95N - 88N

Contraindications

Mini Twister - Suture Anchor is contraindicated in the following cases:

- A.Active infections;
- B.Skeletally immature patients;
- C.Different sites of the body addressed in the use indication;
- D.In the presence of cystic changes, severe osteopenia, or degenerative disease, that can compromise the safety of the suture attachment;
- E.Pathological changes of soft tissues that will be sutured, in order to harm the safety of the suture attachment;
- F.Patients who do not want or are unable to comply with the postoperative instructions due to physical or mental conditions they have;
- G.Sensitivity to a foreign body. Sensitivity tests must be performed to the material in cases of suspected rejection;
- H.Infections and/or restrictions in the blood flow, which can slow the healing and increase the chance of infection and/or rejection of the implant.

Adverse Effects

- A.May occur breakage of the suture;
- B.Occasionally may occur mild erythema and pain at the suture site, due to the trauma of the surgical wound in the tissue;
- C.Anchors can come loose or lose fixation;
- D.Reduction in the bone density and/or bone necrosis due to protection against voltage;
- E.Allergic reaction or sensitivity to the metal of the implant;
- F.Nerve damage due to a surgical trauma;
- G.Disruption of the anchors due to delayed union or disunion of the bone tissue;
- H.Pain, discomfort and/or abnormal sensation due to the presence of the Anchor Screw and Suture Thread;
- I.Shallow and/or deep infection.

Warnings


- A.*Mini Twister - Suture Anchor* is marketed sterile, for a single use. The instrumental present in *Mini Twister - Suture Anchor* should be made useless and discarded after use;
- B.Incomplete insertion of the Anchor Screw will not bring fixation effectiveness;
- C.The postoperative degree of mobility will be determined by the physician and the patient should be informed;
- D.In the case of an implantable device, there is a risk of introducing foreign materials and particles, including gloves talc, materials lint and cleaning agents, in addition to other surface contaminants. All efforts should be made to limit the handling of the attachment devices;
- E.The physician is responsible for the choice and placement of the implant, as well as for the decision to leave or remove implants after the completion of their clinical duties;
- F.Surgical implants must never be reused, and explanted implants must never be implanted again, because it may cause small defects and internal pressure patterns, even if the implant seems to be intact, leading to failure and breakage of the implant.
- G.If the patient is involved in an activity or occupation that may cause stress on the implant (substantial gait, running, weight lifting or muscle strain), there may be failures in the device.
- H.The patient must understand that the Anchor Screw and Suture Thread are not as strong as the normal and healthy bone, and that it may occur Break / Rupture if there is any excess. A patient who can not assimilate this information is in serious danger during the rehabilitation process.

The Patient Must Be Informed Regarding To:

- A.The need for periodic medical follow-up, to observe possible changes of the implant state and adjacent bone. Only the follow-up can detect a possible release of the components or occurrence of osteolysis;
- B.Failure to perform the review surgery when loosening the components or osteolysis can result in a progressive loss of periprosthetic bone stock;
- C.The patient should be instructed to inform that he is an implant holder in the event to undergo the examination of Magnetic Resonance, because the metal-manufactured implant causes interference in the examination damaging the diagnosis, and the exam can cause displacement of the implant causing a problem to the patient;
- D.The patient should understand completely and be warned that the product does not replace and does not have the same performance as the normal bone and, therefore, can break, deform or become loose due to stress or excessive activities, above the established, early load, etc;
- E.The patient should be advised as the need to attend the regular examinations of postoperative follow-up during the whole period when the product remains implanted.

Explanations on Product Use

A. Remove the *Mini Twister - Suture Anchor* of its package keeping its sterile conditions;

- B. Introduce the Screw point in the bone wall;
 - C. Keeping fixed the Screw at the implanta-tion site, turn it clockwise;
 - D. Remove the Applicator from the site;
- E.Make the tying with the Suture Thread.

Special Conditions of Storage, Conservation and/or Handling of the Product

Open Aseptically the Package			
It is Recommended a Temperature of 21 °C (+-4 °C), Environmental Humidity in a Closed, Vented Located and Protected against Weatherproof.			
Do not store directly on the floor (minimum height = 20 cm) and nor in very high places, close to lamps, which could cause dryness of the package or damage to the label.			
Do not store in locations containing contaminant substances such as, for example, cleaning materials, insecticides, pesticides, etc.			
	Do not use if the package is damaged		Keep away from sunlight
	Fragile, handle with care		keep dry
	ESTÉRIL Esterilizado en Gas ETO		"PRODUCTO DE UN SÓLO USO - PROHIBIDO UTILIZAR NUEVAMENTE"

Care Before the Product Usage

Before using the product, the responsible, as designated by the hospital unit, should check the physical conditions of the package and then the product. The responsible shall verify if the cardboard outer package is wrapped by a security seal, film and identified with the information regarding the material present in the package. After this check, the responsible must remove the security seal, and open the cardboard package. The product packed in double blister should be removed from outer package carefully, in a surgical environment. The first blister should be opened by pulling the surgical paper in its spare in the blister edge. The second blister should also be opened by pulling the surgical paper in its spare in the blister edge. With the product in hand, it should be checked if it does not have any cuts, scratches or stains. If there is any adverse effect that affects the safety of the user, as a product not working, damage of the implantable components, serious problems or death related to these components, the surgeon responsible must report this adverse event to the competent health agency and Biomecânica through: e-mail: sac@biomecânica.com.br - phone: (14) 2104-7900. In case of doubts, the surgeon responsible or the health professional may make the communication of the adverse event through the Notifications System on Health Surveillance at the site of ANVISA: <http://www.anvisa.gov.br/hotsite/notvisa/index.htm>

Destruction and Disposal of the Product

All explanted implants should never be implanted again. Stress can cause development of microscopic imperfections and, even if the implant seems to be intact, it can cause its failure. We recommend the pieces to be rasped, bent or cut to their destruction. Explanted devices are considered hospital waste (potentially contaminant product) and should be treated as such, according to the local health authority standards. According to the Resolution RE no. 2605, of 08/11/06, implantable devices of any kind classified as single-use are forbidden from being repro-cessed.

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